

## EXTRACTION

- IN EXTRACTION ROOM -

- Schedule a pre-treatment for mucous samples.
- Homogenize **IC2** and samples to be analysed:

Respiratory specimens

- Add 10 µL **IC2** to the appropriate volume of each sample to be extracted in 1 extraction tube.
- Add 10 µL **IC2** to the appropriate volume of **W0** in 1 extraction tube.
- Perform extraction as below:

INSTRUMENTS	KIT	SAMPLE + IC2 VOLUMES	PROTOCOL	ELUTION VOLUME
	QIAamp® DNA Blood Mini Kit	200 µL of sample + 10 µL IC2		100 µL
	High Pure PCR Template Preparation Kit		Bacteria and Yeast	200 µL
QIAcube	QIAamp® DNA Blood Mini Kit		Blood and body fluid spin protocol V3	100 µL
MagNA Pure Compact®	MagNA Pure Compact Nucleic Acid Isolation Kit I		Total_NA_Plasma_100_400	
MagNA Pure LC System®	MagNA Pure DNA Isolation Kit I		DNA I Blood_Cell High Performance	
NucliSENS® easyMAG™	NucliSENS® Magnetic extraction reagents	Specific A 1.0.2		

Store extracted samples and reagents at -18°C/-22°C

### ANCILLARY REAGENTS:

- **CELL Control r-gene™** - ref.: 71-106
- **Colour Compensation r-gene™** - ref.: 71-103 for LightCycler® 2.0 Instruments use only.

### GLOSSARY:

<b>IC2</b>	Internal Control 2
<b>R11</b>	Bordetella Amplification Premix
<b>CT</b>	Crossing Threshold
<b>W0</b>	Water for extraction
<b>IC2sample</b>	Sample Extraction + Inhibition Control
<b>IC2W0</b>	Reference Extraction + Inhibition Control

## AMPLIFICATION

- IN AMPLIFICATION ROOM -

Validated on:



### AMPLIFICATION PROGRAM

- Enter the following program:

STEPS	TIME	TEMPERATURE	CYCLES	FLUORESCENCE ACQUISITION
Hot Start Taq Polymerase Activation	15 min.	95°C	1	-
Amplification	Denaturation	10 sec.	45	-
	Annealing Elongation	40 sec.		60°C

Note 1: for LightCycler® instruments, set **SEEK TEMPERATURE** to the temperature corresponding to the fluorescence acquisition.

Note 2: for LightCycler® instruments, add a cooling step: 30 sec. / 40°C / 1 cycle at the end of the PCR.

Note 3: for Applied Biosystems instruments, select **NONE** in **PASSIVE REFERENCE**.

Note 4: for Rotor-Gene™, calibrate the signal by clicking on **GAIN OPTIMISATION**.

- Targeted sequence is located in IS481 region. This reagent detects specifically *Bordetella pertussis* and other strains of *Bordetella* containing IS481 region.

### AMPLIFICATION PREPARATION

- Plan number of amplification tubes taking into account:
  - ⇒ number of samples to be amplified;
  - ⇒ positive control (**PC11**);
  - ⇒ negative / inhibition control (**W0**);
- Distribute 15 µL **R11** in all amplification tubes.
- Add 10 µL extracted sample, extracted **W0** or **PC11**.
- Centrifuge for 15 sec.

## RESULTS

### INTERPRETATION

Conditions for test validation:

- Negative Control: NO SIGNAL.
- IC2W0 (560 nm): CT ≤ 32 CYCLES.
- PC11 (530 nm): 25 ≤ CT ≤ 33 CYCLES.

### Sample result interpretation

EXTRACTION and INHIBITION CONTROL	CT [IC2sample] ≤ CT [IC2W0] + 3 cycles		CT [IC2sample] > CT [IC2W0] + 3 cycles	
	NON INHIBITED and correctly extracted sample		INHIBITED or badly extracted sample	
SAMPLE	Calculated CT	Non calculated CT	Calculated CT	Non calculated CT
<b>BORDETELLA detection</b>	Sample validated as positive	Sample validated as negative	Sample validated as positive	NOT VALID

This outlined procedure does not replace the original protocol.

ARGENE tel.: +33 (0)5 61 69 61 00 • infotech@argene.com

## PRINCIPLE for *Bordetella* DETECTION

