

ARGENE

Respiratory Multi Well System r-gene®

Influenza A/B r-gene® - ref.: 71-040 • RSV/hMPV r-gene® - ref.: 71-041
 Rhino&EV/Cc r-gene® - ref.: 71-042 • AdV/hBoV r-gene® - ref.: 71-043
 Chla/Myco pneumo r-gene® - ref.: 71-044 • HCoV/HPIV r-gene® - ref.: 71-045
 Legio pneumo/Cc r-gene® - ref.: 71-046



EXTRACTION

- IN EXTRACTION ROOM -

Respiratory samples* : nasal washes, nasal swabs, nasopharyngeal aspirates, bronchoalveolar liquid (BAL), sputum, tracheobronchial aspirate fluid

* Refer to the technical datasheet of each product to check which respiratory samples are validated

Perform extraction as indicated below :

INSTRUMENTS	KIT	SAMPLE VOLUME	SAMPLE TYPE	PROTOCOL	ELUTION VOLUME	
NucliSENS® easyMAG® (1) (bioMérieux)	NucliSENS® easyMAG® Reagents	200 µL	Respiratory samples	Specific B + 50 µL de diluce	50 µL	
		400 µL			100 µL	
MagNAPure Compact (2) (Roche Diagnostics)	MagNAPure Compact Nucleic acid Isolation Kit I	200 µL		Total_NA_plasma _100_400	Pathogen Complex 200	50 µL
		400 µL				100 µL
QIASymphony SP (3) (Qiagen)	QIASymphony Virus/Bacteria Mini Kit	300 µL			85 µL	

(1) A proteinase K pre-treatment prior to extraction on the NucliSENS® easyMAG® instrument may be required if a sample is considered to be too mucous. In this case, add 10 µL for 200 µL and 20 µL for 400 µL of proteinase K at 20 mg/mL and leave to incubate for 15 minutes at 56°C.

(2) For Legio pneumo/Cc r-gene®, a pre-treatment prior to extraction on the MagNAPure Compact instrument may be required if a sample is considered to be too mucous. In this case, only use 200 µL of sample, add 20 µL of proteinase K at 20 mg/mL (proteinase K, recombinant PCR grade Roche Diagnostics, ref.: 03 115 887 001) and 180 µL of MagNAPure Bacteria Lysis Buffer (Roche Diagnostics ref.: 04 659 180 001). Incubate for 30 minutes at 63°C. Extract the totality, i.e. 400 µL (sample volume to be specified on the instrument), the elution is performed in 50 µL.

(3) Not CE-marked with Legio pneumo/Cc r-gene®

ANCILLARY REAGENTS:

- **CELL Control r-gene®** - ref.: 71-106 kit to check the presence of cells in the sample.
- **DICO Extra r-gene®** - ref.: 71-101
- **RICO Extra r-gene®** - ref.: 71-105

GLOSSARY:

PC4X	Positive Control
W0	Negative Extraction + Amplification Control
R4X	Amplification Premix
CT	Crossing Threshold
RT	Reverse Transcriptase
S	Sample

AMPLIFICATION

- IN AMPLIFICATION ROOM -

LightCycler 480 530 & 560	Applied Biosystems 7500 Fast, 7500 Fast Dx, StepOne, ViiA7 FAM & VIC	Stratagene Mx3005P Agilent Mx3005P Versant kPCR AD FAM & HEX	Dx Real-Time System (Bio-Rad) FAM & HEX	Rotor-Gene Q GREEN & YELLOW
------------------------------	---	---	--	--------------------------------

Refer to the technical datasheet of each product to check which amplification platforms are validated

AMPLIFICATION PROGRAM

STEPS	TIME	TEMPERATURE	CYCLES	FLUORESCENCE ACQUISITION
Reverse Transcription	5 min.	50°C	1	-
Taq Polymerase Activation	15 min.	95°C	1	-
Amplification	Denaturation	10 sec.	45	-
	Annealing	40 sec.		530 + 560 nm
	Elongation	25 sec.		-

Note 1: the temperature variations are parametered by default, which means at 100% or at their maxima.
 Note 2: on LightCycler 480, only "System II" is compatible with the use of the kits. The "System II" features an automatically activated colour compensation inits software.

Note 3: on Applied Biosystems 7500 Fast, 7500 Fast Dx, ViiA7, select **ROX** as **PASSIVE REFERENCE** when programming.

Note 4: on Applied Biosystems StepOne **only**, select **NONE** as **PASSIVE REFERENCE** when programming.

Note 5: on Rotor-Gene Q, calibrate the signal by clicking on **GAIN OPTIMISATION**.

Note 6: on Stratagene Mx3005P or Agilent Mx3005P or Versant kPCR Molecular System AD, select **NONE** in **REFERENCE DYE** when programming.

AMPLIFICATION PREPARATION

Preparation of the premix

- Plan "n" wells for amplification taking into account:
 - number of samples to be amplified;
 - number of Positive Control (**PC4X**);
 - number of Negative Extraction + Amplification Control (**W0**);

Note: for Dx Real-Time System, use transparent plates (ref.: HSP9601) with optical caps (ref.: TCS0803).

Step to be added for RNA premix only

- Pipet 15 µL x (n+1) Amplification Premix (**R4X**);
- Dilute 1 volume **RT** (1 µL) into 9 volumes **H₂O** (9 µL);
- Add 0.15 µL x (n+1) diluted Reverse Transcriptase.

Preparation of PCR assay

- For RNA detection: Distribute 15 µL of the previous premix (**R4X** + diluted **RT**) in the amplification wells;

or

- For DNA detection: Distribute 15 µL of amplification premix (**R4X**) in the amplification wells.
- Add 10 µL of extracted sample, positive control (**PC4X**), extracted **W0**.
- Centrifuge for 15 sec, if applicable.
- Place the plate in the instrument and run the program described above according to the protocol of the instrument used.

RESULTS

TEST VALIDATION CONDITIONS

- **W0** read at 530 nm and 560 nm: **NO SIGNAL.**
- **PC4X** read at 530 nm and 560 nm: **CT < 32 CYCLES.**

INTERPRETATION OF THE RESULTS

RESULTS		Target A - 530 nm	
		CT calculated < 40 cycles	No CT calculated or CT ≥ 40 cycles
Target B - 560 nm	CT calculated < 40 cycles	Target A & B* +	Target A - Target B + if cellular control validated
	No CT calculated or CT ≥ 40 cycles	Target A + Target B - if cellular control validated	Target A & B - if cellular control validated

*For Applied Biosystems StepOne and Rotor-Gene Q apparatus, a crosstalk from 530 nm into 560 nm may be observed.

Note: The evaluation of the efficacy of the extraction and the detection of potential inhibitors is possible using the DICO Extra r-gene® ref.: 71-101 / RICO Extra r-gene® ref.: 71-105 reagents.

Note: Certain batches of transport media used for routing and storage of respiratory samples may contain traces of genomic DNA likely to give a false positive result beyond 35 cycles*. For the same reasons, the use of inhibition controls IC1 (RICO Extra r-gene®) or IC2 (DICO Extra r-gene®) in combination with this cellular control may also engender weakly positive signals beyond 35 cycles. In both cases, **any sample with a CT value for the cellular control target greater than 35 cycles is considered to be a sample with an insufficient number of cells to be interpreted as valid.** This sample must be retested or taken again.

*at 560 nm (cell control channel) for kits containing cellular control (Rhino&EV/Cc r-gene® ref.: 71-042 and Legio pneumo/Cc r-gene® ref.: 71-046)

This outlined procedure does not replace the original protocol
 See Package Insert

Reference 71-040 – 71-040CE_v2 Reference 71-044 – 71-044CE_v2
 Reference 71-041 – 71-041CE_v2 Reference 71-045 – 71-045CE_v2
 Reference 71-042 – 71-042CE_v3 Reference 71-046 – 20682 – A.
 Reference 71-043 – 71-043CE_v2


Principle for Respiratory Multi Well System r-gene® detection

DISTRIBUTION

Example of microplate layout

DISTRIBUTE SAMPLES AND CONTROLS (10 µL PER WELL)

Distribute mix (15 µL per well):
Amplification premix (+ RT if required)



VALIDATION AND INTERPRETATION OF THE RESULTS

e.g.: Detection with the *Influenza A/B r-gene®* - 71-040 kit.

530 nm reading

560 nm reading

TEST VALIDATION CONDITIONS

1st condition: Negative Extraction + Amplification Control (W0)

Well A12



Well A12




NO CT

2nd condition: Positive Control (PC40)

Well A11



Well A11



CT calculated

INTERPRETATION OF THE RESULTS

Well A1



Well A1



If CT calculated:

-> Sample1 Positive for *Influenza A*

If No CT calculated or CT > 40 cycles:

-> Sample1 Negative for *Influenza A*
if cellular control validated

If CT calculated:

-> Sample1 Positive for *Influenza B*

If No CT calculated or CT > 40 cycles:

-> Sample1 Negative for *Influenza B*
if cellular control validated