

CHEMUNEX SCAN CYTOMETRY SOLUTIONS

PIONEERING DIAGNOSTICS



**Advantage of
the Real Time
Microbiology
Testing**

**Chemunex
product range**

**ScanRDI®
analysis**

Applications

Conclusion

Advantage of
the Real Time
Microbiology
Testing

Chemunex
product range

ScanRDI®
analysis

Applications

Conclusion



**Advantage of
the Real Time
Microbiology
Testing**

The challenge

The traditional microbiology

Advantage of the Real Time Microbiology Testing

The challenge




- Microbial contamination is a prime concern → Crucial for industries such as pharmaceutical, biotechnology, personal care, drinking water distribution, breweries...
- Better control on finished product if there is a comprehensive in-process testing at every crucial stage of production
- Delays in microbial testing can directly impact effective consumer protection

The traditional microbiology (1/2)

■ The traditional microbiology : Growth based, agar plate method

Introduced by Robert Koch (1843-1910) and Julius Richard Petri (1852-1921)




Copyright © 2006 Pearson Education, Inc., publishing as Benjamin Cummings.

⇒ 75 % of all microbial testing use the ~130 year old method


The traditional microbiology (2/2)

■ The traditional microbiology : Growth based, agar plate method



➔ **The traditional method is time consuming**

Advantage of the Real Time Microbiology Testing (1/2)



Advantage of the Real Time Microbiology Testing (2/2)



- Avoid the requirement for any cell multiplication
- Potential contamination can be detected immediately
- Storage can be dramatically reduced
- Remedial action can be taken before processes drift out of specifications
- Cleaning actions can be validated immediately
- Process improvements can be evaluated on-line



Chemunex
product range

Flow cytometry

Solid phase cytometry

Why the ScanRDI® ?

Chemunex product range: Flow cytometry



D-Count

- Fully automated flow cytometer
- Up to 350 samples per working shift
- Throughput 50 samples/h (application dependent)
- Capacity up to 64 samples/batch (configuration dependent)



BactiFlow ALS

- Fully automated flow cytometer
- Up to 150 samples per working shift
- Throughput 20 samples/h (application dependent)
- Capacity 1-25 samples/batch



BactiFlow

- Table top flow cytometer for manual use
- Up to 50 samples per day



ChemScan[®] RDI / ScanRDI[®]

The only system allowing
real time
detection & enumeration
of micro-organisms in filterable
samples with a sensitivity
down to one cell

Chemunex product range: Solid phase cytometry (2/2)



ChemScan[®]RDI
ScanRDI[®]

- First generation of ChemScan RDI
- Argon Ion Laser not integrated
- Software 3.4.11



ChemScan[®]RDI SSL
ScanRDI[®] SSL

- Solid state laser (SSL) integrated into the analytical module
- Software 3.4.11



CHEMUNEX[™]
ScanRDI

- Second generation of ScanRDI[®]
- Solid state laser (SSL) integrated into the analytical module
- One holder for CB0.4 and FIFU
- Software 3.4.18

■ *ScanRDI*®: The ultimate combination of speed and sensitivity

■ 1 protocol for all compatible matrices

■ Detection

- Direct detection of bacteria, yeast, molds and spores
- Linear response from 1 to 10^5 cells for bacteria and 1 to 10^4 for yeast and molds
- Spores, stressed cells and fastidious microorganisms detected within minutes
- No multiplication required

■ Sensitivity

- down to one microbial cell in a sample
- independent from the volume filtered (large volume can be tested)

■ Non-destructive test protocol permits microscopic confirmation

■ Robust and easy to use

■ 21 CFR 11 compliant

■ Audit trail



**ScanRDI®
analysis**

A simple three steps procedure

Sample preparation



The ScanRDI® analysis

Data treatment

Results display

Microscope validation

ScanRDI[®] analysis : A simple three steps procedure




Sample preparation : The filtration (1/2)

Sample filtration


Large volumes can be tested using :




CB 0.4



FIFU




0.4µm polyester track-etched membranes (= ChemFilter)




Sample preparation : The filtration (2/2)

The membrane




- **ChemFilter CB 0.4 membrane = A true membrane filter and not a depth filter should be used**
 - 25 mm diameter black track etched polyester membrane
 - 23 μm thick and has a pore specification of 0.36 μm to retain MO on the surface of the filter
 - Pores of accurately controlled dimension are created by ion bombardment of a plastic film followed by controlled chemical etching
 - A tab to facilitate handling to ensure no damage to the membrane during manipulation
- **FIFU : Fluorassure Integral Filtration Unit**
 - Translucent funnel + Blue membrane carrier + Membrane CB04 (0.4 μm black membrane) inserted in square white support + white filtration support + translucent cap
 - 2 white pad supports with one labelling pad each.

Sample preparation : The cell labelling (1/3)



The principle of the cell labelling : Fluorassure viability markers


- ① Viability substrate enters into the cell
- ② Enzymes cleave the viability substrates and release fluorochromes
Enzyme activity
- ③ Fluorochromes accumulate into the cell cytoplasm
Membrane integrity



The Fluorassure viability marker will label all the viable micro-organisms and it is a non destructive method


Sample preparation : The cell labelling (2/3)

Cell labelling – ChemFilter CB0.4





Sample preparation : The cell labelling (3/3)

Cell labelling – ChemFilter FIFU



The ScanRDI® analysis (1/4)




After labeling, the filter is simply placed into the ScanRDI® analyzer and the scan is automatically initiated. All viable microorganisms present are individually detected and counted.






The laser beam (488nm emission wavelength) is guided via a fiber optic cable to a two-axis scanning module which moves the focused laser beam across the entire surface of the membrane.

➔ A laser scans the entire surface of 25 mm diameter membrane in 3 to 5 min

The ScanRDI® analysis (2/4)




- The laser scans the membrane at a speed of 2 meter/sec
- Laser makes 9 545 lines on the membrane




Scanning is fully overlapping → Covering event


The ScanRDI[®] analysis (4/4)



The emission spectrum of the fluorochrome




Datas treatment : Datas processing



The aim of the discriminated setting:
 This discrimination allows to separate the microorganisms from the background


- ***optimise the recovery of organisms***
- ***minimise the possibility of false negative result***

Data Map = Total count



Results Map

= labelled microorganisms




Rejected Background

- Autofluorescent Particles
- Membrane Fluorescence
- Electronic Noise


Results Display

- Results are immediately displayed as direct viable cell counts requiring no operator interpretation.
- In addition, a scan map display shows the precise location of each detected microorganism at the surface of the membrane.
- This enables fast visual result confirmation using an optional microscope.





Microscope Validation (1/2)

- Data on microorganism counts obtained using the ScanRDI® have been through the multi-level discrimination process described above.
- It is possible to confirm that the 'spot' detected is a true organism
- To enable users to perform this task the Cytometer has an epifluorescence microscope with a motorized stage that is directly linked to the Cytometer database.
- The Cytometer can drive the motorized stage to any selected point on the membrane to allow the visual confirmation that the detected event is in fact a microorganism.




1. The membrane holder is placed on the automated microscope stage




2. Validation of the Scan map

Motorized stage


Microscope Validation (2/2)



Candida



Bacillus



Mould

Direct Viability Labelling Demonstrated With Wide Range of Microorganisms

Bacteria Gram -

Achromobacter xylosoxydans
Aeromonas hydrophila
Agrobacterium radiobacter
Alcaligenes eutrophus
Alcaligenes faecalis
Burkholderia cepacia
Burkholderia diminuta
Burkholderia pickettii
Caulobacter sp.
Cedecea lapagei
Citrobacter diversus
Citrobacter freundii
Comamonas terrigena
Edwardsiella hoshinae
Enterobacter aerogenes
Enterobacter agglomerans
Enterobacter cloacae
Enterobacter gergoviae
Enterobacter sakasakii
Enterobacter intermedium
Erwinia Sp.
Escherichia coli
Escherichia coli HB 101
Escherichia coli 0126 :B16
Flavobacterium Sp.
Klebsiella oxytoca
Klebsiella planticola
Klebsiella planticola
Klebsiella pneumoniae
Klebsiella terrigena

Kluyvera Sp.
Moraxella sp.
Pasteurella aerogenes
Proteus mirabilis
Pseudomonas diminuta
Pseudomonas aeruginosa
Pseudomonas alkagenèse
Pseudomonas mesophilica
Pseudomonas putida
Pseudomonas fluorescens
Pseudomonas stutzeri
Salmonella holeraesuis
Salmonella indiana
Salmonella typhimurium
Salmonella eboni
Salmonella sp.
Salmonella virchow
Serratia marcescens
Shigella sonnei
Xanthomonas maltophilia
Yersinia enterocolitica

Bacteria Gram +

Aerococcus viridans
Bacillus anthracis
Bacillus amyloliquefaciens
Bacillus cereus
Bacillus circulans
Bacillus coagulans
Bacillus globigii

Bacillus lentus
Bacillus licheniformis
Bacillus megaterium
Bacillus mycoides
Bacillus pumilus
Bacillus sphaericus
Bacillus stearothermophilus
Bacillus subtilis
Bacillus thuringiensis
Bacteroides fragilis
Bacteroides thetaiotamicron
Bacteroides vulgatus
Clostridium acetobutylicum
Clostridium bifermentans
Clostridium butyricum
Clostridium perfringens
Clostridium sporogenes
Clostridium tyrobutyricum
Corynebacterium aquaticum
Corynebacterium pseudodiphtheriticum
Enterococcus faecium
Enterococcus faecalis
Fusobacterium nucleatum
Lactobacillus acidophilus
Lactobacillus brevis
Lactobacillus buchneri
Lactobacillus bulgaricus
Lactobacillus casei casei
Lactobacillus casei
Lactobacillus cellobiosus
Lactobacillus curvatus

Lactobacillus delbrueckii
Lactobacillus fermentum
Lactobacillus leichmannii
Lactobacillus plantarum
Lactobacillus lactis
Lactobacillus sake
Lactobacillus sp.
Leuconostoc oenos
Leuconostoc Sp.
Listeria innocua
Listeria monocytogenes
Micrococcus luteus
Mycobacterium bovis
Mycobacterium parafortuitum
Mycobacterium mageritense
Mycobacterium tuberculosis oerskovia sp.
Pediococcus damnosus
Pediococcus pentosaceus
Porphyromonas canoris
Porphyromonas gingivalis
Propionibacterium acnes
Staphylococcus aureus
Staphylococcus epidermidis
Staphylococcus hominis
Staphylococcus warneri
Staphylococcus xylosus
Streptococcus faecalis
Streptococcus salivarius
Streptococcus thermophilus
Streptococcus viridans
Thiobacillus ferrooxidans

Yeast

Acremonium kiliense
Candida albicans
Candida ciferii
Candida colliculosa
Candida famata
Candida famata
Candida fumentans
Candida humicola
Candida humicola
Candida krusei
Candida luxitaniae
Candida magnolia
Candida parapsilosis
Candida pelliculosa
Candida tropicalis
Cryptococcus albidus
Debaryomyces hansenii
Galactomyces
geotrichum Geotrichum
candidum
Hansenulasporea uvarum
Hansenula anomala
Kloechera japonica
Kloechera Apis apiculata
Pichia anomala
Pichia guillermondii
Pichia menbrena faciens
Rhodotorula rubra
Saccharomyces baillii
Saccharomyces
bisparus
Saccharomyces
cerevisiae
Saccharomyces rosei

Torulopsis candida
Torulopsis inconspicua
Torulopsis maris
Torulosporea delbrueckii
Zygosaccharomyces baillii
Zygosaccharomyces rouxii

Mould

Acremonium Sp.
Aspergillus versicolor
Aspergillus versicolor
Aspergillus fumigatus
Aspergillus niger
Basydiomycetes Sp.
Bassochlamis fulva
Byssochlamys Sp.
Cladosporium
cladosporioides Epicocum
nigrum ou altenaria Fusarium
oxysporum Fusarium
oxysporum
Fusarium gramineatum
roseum Humicola fuscoatra
Mucor circinelloides
Mucor plumbeus
Mucor racemosus
Mucor Sp.
Neosartoea Sp.
Penicillium decumbens
Penicillium expansum
Penicillium frequentans
Penicillium roquefortii
Rhizopus Sp.
Rhodoturola rubra
Rhizopus oligosporus
Scopulariopsis candida
Trichoderma Sp.

Pharmaceutical applications

Exemple of applications in the field:

Time to results

In-process analysis

TVC Bioburden for in-process

90 min

TVC Bioburden for raw material

90 min



Environmental controls

TVC Bioburden for pharmaceutical water

90 min

Air monitoring using Coriolis

< 3 hours



Surface monitoring using ChemSwab

< 3 hours



Biotechnology

Contaminations of cell cultures

< 2 hours

Control of fermentations

90 min

Finish product testing

Scan Bio II for sterility test

< 3 hours



- One protocol (Scan Bio II) for all compatible matrices

- **Improve productivity and maximize yields**
 - Real-time detection and correction of contamination problems
 - Minimize plant down time for decontamination/cleaning
 - Immediate cleaning validation
 - Reduce stock → Save on warehouse costs

- **Guarantee quality**
 - Reduce the risk for contamination → Increase consumer protection
 - Decrease the likelihood for recall
 - Rapidly test new product developments