

## Plasmid purification products from MACHERY-NAGEL

*MN guide to plasmid purification*

*Find the optimal solution*

**NucleoBond®**

**NucleoSpin®**

**Superior yields**

**Outstanding purities**

**Time-saving procedures**

**... for reliable downstream applications**

**MACHERY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



*Since 1911*





# MN technologies for plasmid purification


## NucleoBond®

### Anion-exchange technology

	Page
NucleoBond® Xtra Midi / Maxi _____	3
NucleoBond® Xtra Midi / Maxi EF _____	4
NucleoBond® Xtra BAC _____	5



Principle of binding: Ionic bond





Principle of elution: pH shift

- Features / Results**
- Ultra-pure, transfection-grade plasmid DNA
  - The new generation of anion exchangers
  - Xtra fast, Xtra high yield, Xtra convenient


## NucleoSpin®

### Silica-membrane technology

NucleoSpin® Plasmid / Plasmid (NoLid) _____	6
NucleoSpin® Plasmid QuickPure _____	7
NucleoSpin® 8 / 96 Plasmid _____	8



Principle of binding: Removing of hydrate shell  
Hydrogen bonds / Salt bridges



Principle of elution: Reconstitution of hydrate shell

- Features / Results**
- Sequencing and PCR-grade plasmid DNA
  - No alcohol precipitation necessary
  - Fast and easy procedure

# NucleoBond® Xtra Midi • NucleoBond® Xtra Maxi

Superior plasmid Midi and Maxi kits – transfection-grade

## ▶ Highest speed – 30 min (Midi)/35 min (Maxi)

Optimal column design and inserted filters – avoid inconvenient syringes

## ▶ Highest yield – typically 250 µg for Midi and 1000 µg for Maxi

Improved silica material – patented design!

## ▶ High purity – transfection-grade plasmid DNA

Established anion-exchange technology

## ▶ NucleoBond® Xtra Plus kits for super high-speed version


NucleoBond® Finalizer for omitting time-consuming centrifugation after DNA precipitation



## Product at a glance

Sample material	<b>Midi</b> < 200 mL (high copy) < 400 mL (low copy)		<b>Maxi</b> < 600 mL (high copy) < 1200 mL (low copy)	
Vector size	< 300 kbp		< 300 kbp	
Typical yield	250 µg		1000 µg	
Precipitation	<b>Midi</b> Centrifugation	<b>Midi Plus</b> NucleoBond® Finalizer	<b>Maxi</b> Centrifugation	<b>Maxi Plus</b> NucleoBond® Finalizer
Preparation time	70 min / prep	30 min / prep	75 min / prep	35 min / prep
Applications	Transfection, cloning, plasmid banks, etc.		Transfection, cloning, plasmid banks, etc.	

## Application data



\* Yield of plasmid DNA is slightly lower due to residual DNA remaining on the desalting tool (compared to kits without desalting tool).

## Highest yield in less time – comparison to anion-exchange competitor kits

Plasmid DNA was isolated following each manufacturer's protocol using maximum culture volume. Plasmid DNA yield (■) was determined after precipitation. The comparison shows the results of the kits with NucleoBond® Finalizer as desalting tool (Fig. A) and without desalting tool (Fig. B).

## Ordering information

Product	Preps	REF
NucleoBond® Xtra Midi	10/50/100	740410.10/.50/.100
NucleoBond® Xtra Midi Plus (incl. Finalizer)	10/50	740412.10/.50
NucleoBond® Xtra Maxi	10/50/100	740414.10/.50/.100
NucleoBond® Xtra Maxi Plus (incl. Finalizer)	10/50	740416.10/.50
NucleoBond® Xtra Combi Rack	-	740415

# NucleoBond® Xtra Midi EF • NucleoBond® Xtra Maxi EF

Superior plasmid kit for transfection of sensitive cell lines – endotoxin-free



▶ **Highest purity – endotoxin-free transfection of sensitive cell lines**

Patented endotoxin removal

▶ **Highest speed – 45 min (Midi)/50 min (Maxi)**

No extra incubation for endotoxin removal

Optimal column design and inserted filters – avoid inconvenient syringes

▶ **Highest yield – typically 250 µg for Midi and 1000 µg for Maxi**

Improved silica material – patented design!


▶ **NucleoBond® Xtra Plus EF kits for super high speed version**

NucleoBond® Finalizer for omitting time-consuming centrifugation after DNA precipitation

## Product at a glance

Sample material	Midi EF		Maxi EF	
		< 200 mL (high copy)	< 400 mL (low copy)	< 600 mL (high copy)
Vector size	< 300 kbp		< 300 kbp	
Typical yield	250 µg		1000 µg	
Precipitation	<b>Midi EF</b> Centrifugation	<b>Midi Plus EF</b> NucleoBond® Finalizer	<b>Maxi EF</b> Centrifugation	<b>Maxi Plus EF</b> NucleoBond® Finalizer
Preparation time	85 min / prep	45 min / prep	90 min / prep	50 min / prep
Applications	Transfection, cloning, plasmid banks, etc.		Transfection, cloning, plasmid banks, etc.	

## Application data



### Transfection of PC12 cells with plasmid DNA purified using NucleoBond® Xtra EF

*Plasmid purification:* 9.2 kbp plasmid, bacterial strain: *E. coli* DH5α, purification with NucleoBond® Xtra Midi EF

*Transfection experiment:* highly sensitive PC12 cells, *in vitro* cultivation, transfection with 0.5 µg plasmid DNA, transfection method: Lipofectamine™ LTX Reagent (1 µL) + 100 µL DMEM, cells were imaged by fluorescence microscopy 24 h post-transfection

*Results:* Cells show transient expression of Phogrin (PTP NP, Protein Tyrosine Phosphatase), localized by red fluorescent antibody, nuclei are localized by DAPI staining (violet)

Data kindly provided by AG Knauer, Duisburg-Essen University, Institute for Molecular Biology II, Essen, Germany

## Ordering information

Product	Preps	REF
<b>NucleoBond® Xtra Midi EF</b>	10 / 50	740420.10 / .50
<b>NucleoBond® Xtra Midi Plus EF (incl. Finalizer)</b>	10 / 50	740422.10 / .50
<b>NucleoBond® Xtra Maxi EF</b>	10 / 50	740424.10 / .50
<b>NucleoBond® Xtra Maxi Plus EF (incl. Finalizer)</b>	10 / 50	740426.10 / .50
<b>NucleoBond® Xtra Combi Rack</b>	-	740415

# NucleoBond® Xtra BAC

## Kit for large construct DNA



### ▶ High yield – up to 150 µg BAC DNA

Optimal silica material for large constructs

### ▶ High speed – 75 minutes per prep

Optimal column design and inserted filters for high flow rates and parallel lysate clearing and loading

### ▶ High purity – transfection-grade BAC DNA

Established anion-exchange technology

## Product at a glance

Sample material	250–750 mL <i>E. coli</i> culture
Vector size	< 300 kbp
Typical yield	10–150 µg
Preparation time	75 min/prep
Applications	Transfection, cloning, restriction analyses, sequencing

## Application data



### Superior yield in less time – comparison to competitor kits

BAC DNA (300 kbp) was isolated in triplicate from 500 mL *E. coli* DH5α using NucleoBond® Xtra BAC and competitor products (Q and I). After precipitation, BAC DNA was reconstituted in 1000 µL TE buffer and 5 µL of each sample was used for analysis on a 1% TAE-agarose gel.

Obtained yields: MN: 150 µg, Q: 44 µg, I: 75 µg



### Suitable for restriction digestion

Restriction digestion of BAC DNA isolated with NucleoBond® Xtra BAC. BAC DNA samples were purified from overnight cultures of *E. coli* DH5α transformed with a pBAC10 clone. Approximately 3 µg of DNA from each sample was digested with 3 units of *MspI*, *HindIII*, or *EcoRI* at 37 °C for 2 hours. Digestions were analysed on a 1% TAE agarose gel. Lane 1: undigested BAC DNA, lane 2: *MspI* digested, lane 3: *HindIII* digested, lane 4: *EcoRI* digested, lane M: marker.

## Ordering information

Product	Preps	REF
NucleoBond® Xtra BAC	10/25	740436.10/.25
NucleoBond® Xtra Combi Rack	-	740415

# NucleoSpin® Plasmid · NucleoSpin® Plasmid (NoLid)

## High-yield plasmid Mini prep – molecular biology-/sequencing-grade

### ▶ Highest yield

Up to 40 µg plasmid DNA

### ▶ Highest quality

Phred score of 20 (up to 1000 bp)

### ▶ Fast procedure

18 Mini preps in 25 min


### ▶ Choose between columns with or without a lid



## Product at a glance

Sample material	1–5 mL <i>E. coli</i> cultures (standard protocol) 5–10 mL <i>E. coli</i> cultures*
Vector size	< 15 kbp
Typical yield	< 25 µg (from 1–5 mL) < 40 µg (from 5–10 mL)*
Elution volume	50 µL
Binding capacity	60 µg
Applications	Sequencing, cloning, PCR, restriction analysis, transformation, enzymatic modifications

## Application data



### High-yield plasmid Mini prep

Plasmid DNA isolation (pUC18) from *E. coli* DH5α using NucleoSpin® Plasmid. For the isolation 2 mL (lane 1–2), 5 mL (lane 3–4) and 8 mL (lane 5–6) LB cultures were used and analyzed on an agarose gel (2 µL of each eluate). Lane M: marker.

## Ordering information

Product	Preps	REF
NucleoSpin® Plasmid	10 / 50 / 250	740588.10 / .50 / .250
NucleoSpin® Plasmid (NoLid)	50 / 250	740499.50 / .250
NucleoSpin® Plasmid Buffer Set*	300	740953

\* If larger culture volumes are frequently used, increased buffer volumes are required. Please see „ordering information“

# NucleoSpin® Plasmid QuickPure

High-speed plasmid Mini prep – molecular biology-/sequencing-grade

▶ **Highest speed**

18 Mini preps in 11 min

▶ **Highest quality**

Phred score of 20 (up to 1000 bp)

▶ **Reliable yield**


Up to 15 µg plasmid DNA



## Product at a glance


Sample material	1–3 mL <i>E. coli</i> cultures
Vector size	< 15 kbp
Typical yield	< 15 µg
Elution volume	50 µL
Binding capacity	15 µg
Applications	Sequencing, cloning, PCR, restriction analysis, transformation, enzymatic modifications

## Application data



### Comparison of plasmid yield

Plasmid DNA isolation (pUC18) from a 3 mL LB *E. coli* DH5α culture. Two isolations were each performed according to the manufacturer's protocol. For the comparison 4 µL of each 50 µL eluate were loaded on a 1 % TAE agarose gel. The yield of NucleoSpin® Plasmid QuickPure is higher or equal to competitors Q, Z, S, I, P. Lane M: marker.



**Just three steps = Less than 10 min per prep**

## Ordering information

Product	Preps	REF
NucleoSpin® Plasmid QuickPure	10/50/250	740615.10/.50/.250

# NucleoSpin® 8/96 Plasmid · NucleoSpin® 8/96 Plasmid Core Kit

## High-throughput plasmid Mini prep – molecular biology-/sequencing-grade

- ▶ **Suitable for up to 5 mL bacterial culture**
- ▶ **Convenient filtration of bacterial lysates**  
NucleoSpin® 8 Filter Strips / 96 Filter Plate  
Optimized filter material minimizes the risk of filter clogging
- ▶ **Minimizes risk of cross-contamination**  
MN Wash Plate technology
- ▶ **Core Kit versions available**  
Cost savings – kits with basic contents\*



### Product at a glance

Format	8-well strips / 96-well plates
Processing	Manual and automated, vacuum or centrifugation**
Lysate clarification	8-well filter strips / 96-well filter plates
Sample material	1–5 mL <i>E. coli</i> culture
Vector size	<15 kbp
Typical yield	4–6 µg/mL <i>E. coli</i> culture
Elution volume	75–150 µL
Preparation time	45 min / 6 strips or 1 plate
Binding capacity	20 µg
Applications	Sequencing, cloning, PCR, restriction analysis, transformation, enzymatic modifications

### Ordering information

Product	Preps	REF
<b>NucleoSpin® 8 Plasmid</b>	12 x 8 / 60 x 8	740621 / .5
<b>NucleoSpin® 8 Plasmid Core Kit*</b>	48 x 8	740461.4
<b>NucleoSpin® 96 Plasmid</b>	1 x 96 / 4 x 96 / 24 x 96	740625.1 / .4 / .24
<b>NucleoSpin® 96 Plasmid Core Kit*</b>	4 x 96	740616.4

\* Core kit contents are focused on automation platforms. Additional accessories can be combined as needed. Please check [www.mn-net.com](http://www.mn-net.com) or contact our Technical Support.

\*\* Support protocol available at Technical Support (tech-bio@mn-net.com)

### Trademarks

**MACHEREY-NAGEL:** NucleoSpin®, NucleoBond®  
**Other companies:** Lipofectamin™ [Invitrogen]

Your local distributor

[www.mn-net.com](http://www.mn-net.com)

**MACHEREY-NAGEL**



**MACHEREY-NAGEL GmbH & Co. KG** · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

**Germany**

**and international:**

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: [info@mn-net.com](mailto:info@mn-net.com)

**Switzerland:**

**MACHEREY-NAGEL AG**

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

**France:**

**MACHEREY-NAGEL EUROL**

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

**USA:**

**MACHEREY-NAGEL Inc.**

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)



Since 1911