## Proteinase K

### 1 Contents

| Proteinase K               |                       |           |           |
|----------------------------|-----------------------|-----------|-----------|
| REF                        | 740506                | 740506.75 | 740506.30 |
| Proteinase K (lyophilized) | 100 mg<br>(2 x 50 mg) | 75 mg     | 30 mg     |
| Proteinase Buffer PB       | 8 mL                  | 8 mL      | 1.8 mL    |

### 2 Product description

Proteinase K is delivered with a stabilizing solution (Proteinase Buffer PB) enabling storage for up to six months at -20 °C without decrease of proteinase activity. Stability is considerably increased if Proteinase K is dissolved in Proteinase Buffer compared to Proteinase K dissolved in water. Proteinase Buffer is compatible with NucleoSpin® and NucleoBond® kits and also suitable for most other applications.

*Note: NucleoSpin® or NucleoBond® kits may require different Proteinase K stock solution concentrations (e.g., NucleoSpin® Tissue: 20 mg/mL; NucleoSpin® 8/96 Tissue: 29 mg/mL). See individual protocols for detailed information.*

### 3 Proteinase K specifications

|   |  |
|---|--|
| <b>Source:</b>                          | <i>Tritirachium album</i>  |
| <b>Form:</b>                            | 2 x 50 mg / 75 mg / 30 mg,<br>white freeze-dried material  |
| <b>Activity:</b>                        | > 2.5 U/mg material (Casein, 37 °C)  |
| <b>Solubility:</b>                      | Readily soluble in a concentration of<br>20 mg/mL in analytical grade water or PB  |
| <b>DNA content:</b>                     | content: < 2 ppm   |
| <b>Contaminating deoxyribonuclease:</b> | No degradation of genomic DNA observed after incubation for at least 6 h at 37 °C in a typical restriction digest buffer: 33 mM Tris-acetate, pH 7.9, 10 mM magnesium acetate, 66 mM potassium acetate, 0.1 mg/mL BSA. |

### 4 Storage conditions and preparation of working solution

Lyophilized Proteinase K and Proteinase Buffer PB can be stored at room temperature (18–25 °C) and are stable up to one year.

Add the indicated volume of Proteinase Buffer PB to dissolve lyophilized Proteinase K in order to obtain a 20 mg/mL solution. Proteinase K solution is stable at -20 °C for up to 6 months.

See the individual protocols for detailed information, if other concentrations are required.

| Proteinase K               |   |  |                                       |
|----------------------------|---|--|---------------------------------------|
| REF                        | 740506  | 740506.75                              | 740506.30                             |
| Proteinase K (lyophilized) | 2 x 50 mg   | 75 mg                                  | 30 mg                                 |
|                            | Add 2.5 mL<br>Proteinase<br>Buffer PB<br>to each vial | Add 3.75 mL<br>Proteinase<br>Buffer PB | Add 1.5 mL<br>Proteinase<br>Buffer PB |

## 5 Safety instructions

Proteinase K contains hazardous contents. Wear gloves and goggles!

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### 5.2 GHS classification

Only harmful features must not be labeled with H and P phrases until 125 mL or 125 g.

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For further information please see Material Safety Data Sheets ([www.mn-net.com](http://www.mn-net.com)).

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