



RNA clean-up

User manual

NucleoSpin® RNA Clean-up

March 2012/Rev. 03

Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Reagents, consumables, and equipment to be supplied by user	5
1.3	About this user manual	5
2	Product description	6
2.1	The basic principle	6
2.2	Kit specifications	6
2.3	Handling, preparation, and storage of starting materials	7
2.4	Elution procedures	8
3	Storage conditions and preparation of working solutions	9
4	Safety instructions – risk and safety phrases	10
4.1	Risk and safety phrases	10
4.2	GHS classification	11
5	Protocols	12
5.1	RNA Clean-up	12
5.2	RNA isolation from up to 10 ⁵ cells	14
6	Appendix	16
6.1	Troubleshooting	16
6.2	Ordering information	18
6.3	Product use restriction/warranty	20

1 Components

1.1 Kit contents

NucleoSpin® RNA Clean-up			
REF	10 preps 740948.10	50 preps 740948.50	250 preps 740948.250
Lysis Buffer RA1	10 mL	25 mL	125 mL
Wash Buffer RA2	15 mL	15 mL	80 mL
Wash Buffer RA3 (Concentrate)*	5 mL	12.5 mL	3 x 25 mL
RNase-free H ₂ O	5 mL	15 mL	65 mL
NucleoSpin® RNA Binding Columns (light blue rings – plus Collection Tubes)	10	50	250
Collection Tubes (2 mL)	10	50	250
Collection Tubes (1.5 mL)	10	50	250
User manual	1	1	1

* For preparation of working solutions and storage conditions see section 3.

1.2 Reagents, consumables, and equipment to be supplied by user

Reagents

- 96–100 % ethanol (to prepare Wash Buffer RA3 and to adjust RNA binding conditions)

Consumables

- 1.5 mL microcentrifuge tubes
- Sterile RNase-free tips

Equipment

- Manual pipettors
- Centrifuge for microcentrifuge tubes
- Personal protection equipment (e.g., lab coat, gloves, goggles)

1.3 About this user manual

It is strongly recommended reading the detailed protocol sections of this user manual if the **NucleoSpin® RNA Clean-up** kit is used for the first time. Experienced users, however, may refer to the Protocol-at-a-glance instead. The Protocol-at-a-glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure.

All technical literature is available on the internet at **www.mn-net.com**.

Please contact Technical Service regarding information about changes of the current user manual compared to previous revisions.

2 Product description

2.1 The basic principle

One of the most important aspects in the isolation and handling of RNA is to prevent degradation of the RNA during the isolation procedure. With the **NucleoSpin® RNA Clean-up** kit, RNA containing samples are mixed with a solution containing large amounts of chaotropic ions. This solution immediately inactivates RNases – which are present in virtually all biological materials – and creates appropriate binding conditions which favor adsorption of RNA to the silica membrane. Simple washing steps remove salts, metabolites, organics like phenol, and macromolecular cellular components. Pure RNA is finally eluted under low ionic strength conditions with RNase-free water (supplied).

The RNA clean-up preparation using **NucleoSpin® RNA Clean-up** kits can be performed at room temperature. The eluate, however, should be treated with care because RNA is very sensitive to trace contaminations of RNases, often found on general lab ware, fingerprints and dust. To ensure RNA stability keep RNA frozen at -20 °C for short-term or -70 °C for long-term storage.

2.2 Kit specifications

- **NucleoSpin® RNA Clean-up** kits are ideal for the clean-up of total RNA from RNA preparations which contain unacceptable amounts of RT-PCR inhibitors (e.g., RNA prepared with phenol-chloroform based methods).
- The kit is further recommended for the isolation of RNA from small amounts of cultured cells whenever copurification of some genomic DNA is acceptable. The kits allow purification of pure RNA with an A_{260}/A_{280} ratio generally exceeding 1.9 (measured in TE buffer (pH 7.5)).
- **NucleoSpin® RNA Clean-up** kits are recommended for the clean-up of RNA from enzymatic reactions like *in vitro* transcribed RNA, amplification reactions, biotinylated RNA, or fluorescent (Cy dye) labeled RNA.
- The purified RNA is ready to use for applications like enzymatic labelling reactions (e.g., dye incorporation), reverse transcriptase-PCR (RT-PCR), and for DNA/RNA based chip hybridisations (e.g., MWG rat microarray, MWG, Ebersberg, Germany or Human Genome U133A Array, Affymetrix, USA).
- Integrity of purified RNA, originally isolated from for example eukaryotic cells, is examined by denaturing agarose gel electrophoresis: rRNA bands are sharp, with the 28S band being about twice as intense as the 18S band.
- The standard protocol (section 5.1) allows the clean-up of up to 200 µg of RNA per NucleoSpin® RNA Binding Column or the isolation of total RNA from up to 1×10^5 cultured cells (section 5.2).

Table 1: Kit specifications at a glance

Parameter	NucleoSpin® RNA Clean-up
Technology	Silica-membrane technology
Format	Mini spin columns
Sample material	< 100 µL RNA sample with single column loading containing up to 200 µg RNA < 200 µL RNA sample with double column loading containing up to 200 µg RNA Up to 10 ⁵ cells
Fragment size	> 200 nt
Typical recovery (0.1–200 µg RNA input)	85–95 %
A ₂₆₀ /A ₂₈₀	1.9–2.1
Elution volume	40–120 µL
Preparation time	Approx. 20 min/6 preps
Binding capacity	200 µg

2.3 Handling, preparation, and storage of starting materials

RNA intended to be used as sample for the **NucleoSpin® RNA Clean-up** procedure should be handled with the same care as any RNA sample. The stability of prepurified RNA samples (e.g., RNA isolated with phenol based protocols) depends very much on the performed procedure. RNA in biological samples is not protected against digestion until the sample material is flash frozen or disrupted in the presence of RNase inhibiting or denaturing agents. Therefore it is important that biological samples are flash frozen in liquid N₂ immediately and stored at -70 °C or processed as soon as possible. Samples can be stored in lysis buffer after disruption at -70 °C for up to one year, at +4 °C for up to 24 hours or up to several hours at room temperature. Frozen samples are stable up to 6 months. Frozen samples in lysis buffer should be thawed slowly before starting with the isolation of total RNA.

Wear gloves at all times during the preparation. Change gloves frequently.

2.4 Elution procedures

It is possible to adjust the elution method and the volume of RNase-free water used for the subsequent application of interest. In addition to the standard method described in the individual protocols (recovery rate about 70–90 %) there are several modifications possible:

- **High yield:** Perform two elution steps with the volume indicated in the individual protocol. About 90–100 % of bound nucleic acid will be eluted.
- **High yield and high concentration:** Elute with the standard elution volume and apply the eluate once more onto the column for reelution.

Eluted RNA should immediately be placed and always kept on ice for optimal stability because almost omnipresent RNases (general lab ware, fingerprints, dust) will degrade RNA. For short-term storage freeze at -20 °C, for long-term storage freeze at -70 °C.

3 Storage conditions and preparation of working solutions

Attention: Buffers RA1 and RA2 contain chaotropic salt. Wear gloves and goggles!

CAUTION: Buffers RA1 and RA2 contain guanidinium thiocyanate which can form highly reactive compounds when combined with bleach (sodium hypochlorite). **DO NOT** add bleach or acidic solutions directly to the sample-preparation waste.

- All kit components should be stored at room temperature (18–25 °C) and are stable up to one year. Storage at lower temperatures may cause precipitation of salts.
- Check that 96–100 % ethanol is available as additional solution in the lab.

Before starting any **NucleoSpin® RNA Clean-up** protocol, prepare the following:

- **Wash Buffer RA3:** Add the indicated volume of 96–100 % ethanol (see table below) to Wash Buffer RA3 Concentrate. Mark the label of the bottle to indicate that ethanol was added. Store Wash Buffer RA3 at room temperature (18–25 °C) for up to one year.

NucleoSpin® RNA Clean-up			
REF	10 preps 740948.10	50 preps 740948.50	250 preps 740948.250
Wash Buffer RA3 (Concentrate)	5 mL Add 20 mL ethanol	12.5 mL Add 50 mL ethanol	3 x 25 mL Add 100 mL ethanol to each bottle

4 Safety instructions – risk and safety phrases

The following components of the **NucleoSpin® RNA Clean-up** kits contain hazardous contents.

Wear gloves and goggles and follow the safety instructions given in this section.

4.1 Risk and safety phrases

Component	Hazard contents	Hazard symbol	Risk phrases	Safety phrases
<i>Inhalt</i>	<i>Gefahrstoff</i>	<i>Gefahrstoff-symbol</i>	<i>R-Sätze</i>	<i>S-Sätze</i>
RA1	Guanidinium thiocyanate 30–60 % <i>Guanidiniumthiocyanat 30–60 %</i>	✘ Xn*	R 20/21/22-32-52/53	S 13-61
RA2	Guanidinium thiocyanate 30–60 % + ethanol 20-35 % <i>Guanidiniumthiocyanat 30–60 % + Ethanol 20–35 %</i>	✘ Xn*	R 10-20/21/22-32-52/53	S 13-16-61

Risk phrases

- R 10 Flammable.
Entzündlich.
- R 20/21/22 Harmful by inhalation, in contact with skin and if swallowed.
Gesundheitsschädlich beim Einatmen, Verschlucken und Berührung mit der Haut.
- R 32 Contact with acids liberates very toxic gas.
Entwickelt bei Berührung mit Säure sehr giftige Gase.
- R 52/53 Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
Schädlich für Wasserorganismen, kann in Gewässern längerfristig schädliche Wirkungen haben.

Safety phrases



- S 13 Keep away from food, drink and animal foodstuffs.
Von Nahrungsmitteln, Getränken und Futtermitteln fernhalten.
- S 16 Keep away from sources of ignition – No smoking.
Von Zündquellen fernhalten – Nicht rauchen!
- S 61 Avoid release to the environment. Refer to special instructions / safety data sheet.
Freisetzung in die Umwelt vermeiden. Besondere Anweisungen einholen / Sicherheitsdatenblätter zu Rate ziehen.

* Hazard labeling not necessary if quantity per bottle below 125 g or ml (certificate of exemption according to 67/548/EEC Art. 25, 1999/45/EC Art. 12 and German GefStoffV § 20 (3) and TRGS 200 7.1). For further information see Material Safety Data Sheet.

4.2 GHS classification

Only harmful features do not need to be labeled with H and P phrases until 125 mL or 125 g.

Mindergefährliche Eigenschaften müssen bis 125 mL oder 125 g nicht mit H- und P-Sätzen gekennzeichnet werden.

Component	Hazard contents	GHS symbol	Hazard phrases	Precaution phrases
Inhalt	Gefahrstoff	GHS Symbol	H-Sätze	P-Sätze
RA1	Guanidinium thiocyanate 30–60 % <i>Guanidiniumthiocyanat</i> 30–60 %	 Warning Achtung	302, 412, EUH031	260, 273, 301+312, 330
RA2	Guanidinium thiocyanate 30–60 % + ethanol 20- 35 % <i>Guanidiniumthiocyanat</i> 30–60 % + <i>Ethanol</i> 20–35 %	 Warning Achtung	226, 302, 412, EUH031	210, 233, 260, 273, 301+312, 330, 403+235

Hazard phrases

H 226	Flammable liquid and vapour. <i>Flüssigkeit und Dampf entzündbar.</i>
H 302	Harmful if swallowed. <i>Gesundheitsschädlich bei Verschlucken.</i>
H 412	Harmful to aquatic life with long lasting effects. <i>Schädlich für Wasserorganismen, mit langfristiger Wirkung.</i>
EUH031	Contact with acids liberates toxic gas. <i>Entwickelt bei Berührung mit Säure giftige Gase.</i>

Precaution phrases

P 210	Keep away from heat / sparks / open flames / hot surfaces – No smoking. <i>Von Hitze / Funken / offener Flamme / heißen Oberflächen fernhalten.</i>
P 233	Keep container tightly closed <i>Behälter dicht verschlossen halten.</i>
P 260	Do not breathe vapours. <i>Dampf nicht einatmen.</i>
P 273	Avoid release to the environment. <i>Freisetzung in die Umwelt vermeiden.</i>
P 301+312	IF SWALLOWED: Call a POISON CENTER or doctor /physician if you feel unwell. <i>BEI VERSCHLUCKEN: Bei Unwohlsein GIFTINFORMATIONSZENTRUM oder Arzt anrufen.</i>
P 330	Rinse mouth. <i>Mund ausspülen.</i> <i>Bei Symptomen der Atemwege: GIFTINFORMATIONSZENTRUM oder Arzt anrufen.</i>
P 403+235	Store in a well ventilated place. Keep cool. <i>Kühl an einem gut belüfteten Ort aufbewahren.</i>

For further information please see Material Safety Data Sheets (www.mn-net.com).
Weiterführende Informationen finden Sie in den Sicherheitsdatenblättern (www.mn-net.com).

5 Protocols

5.1 RNA Clean-up

Before starting the preparation:

- Check if Wash Buffer RA3 was prepared according to section 3.

1 Sample preparation

Fill up RNA samples smaller than 100 μL with RNase-free water to **100 μL** .

RNA samples from 100–200 μL should be filled up with RNase-free water to 200 μL .



Fill up RNA sample to 100 μL with water

2 Preparation of lysis-binding buffer premix

Prepare a Buffer RA1-ethanol premix with a ratio of 1:1.

For each **100 μL RNA sample** mix **300 μL Buffer RA1** and **300 μL of ethanol (96–100 %)**.

If multiple samples are processed, the preparation of a master-premix is recommended (e.g., 2 mL Buffer RA1 + 2 mL 98% ethanol for approximately 6 preparations).

Prepare premix:

Mix 300 μL RA1 with 300 μL ethanol (96–100 %)

3 Adjust RNA binding conditions

To **100 μL RNA sample** add **600 μL (6 volumes) of Buffer RA1-ethanol-premix**. Mix sample with premix by vortexing.

If a 200 μL RNA sample is processed, add 1200 μL Buffer RA1-ethanol premix.

After addition of ethanol a stringy precipitate may become visible which will not affect the RNA isolation. Be sure to mix thoroughly and apply sample as homogeneous solution onto the column.



+ 6 vol. premix

Mix

4 Bind RNA

For each preparation, take one NucleoSpin® RNA Binding Column (light blue ring) placed in a Collection Tube and load the lysate (700 µL).

Centrifuge for **30 s** at **8,000 x g**. Discard Collection Tube with flow-through and place the column in a new Collection Tube.

Maximal loading capacity of NucleoSpin® RNA Binding Columns is 750 µL. Repeat the procedure if larger volumes are to be processed.



Load 700 µL lysate



**8,000 x g,
30 s**

5 Wash and dry silica membrane

1st wash

Add **700 µL Buffer RA3** to the NucleoSpin® RNA Binding Column. Centrifuge for **30 s** at **8,000 x g**. Discard flow-through and reuse Collection Tube.

+ 700 µL RA3

**8,000 x g,
30 s**

2nd wash

Add **350 µL Buffer RA3** to the NucleoSpin® RNA Binding Column. Centrifuge for **2 min** at **8,000 x g**.

Transfer the NucleoSpin® RNA Binding Column to a nuclease-free Collection Tube (1.5 mL, supplied). Open the lid of the column and let the membrane dry for 3 min.

If for any reason, the liquid level in the Collection Tube has reached the NucleoSpin® RNA Binding Column after centrifugation, discard flow-through and centrifuge again.

The procedure ensures complete removal of ethanol from the column.



+ 350 µL RA3



**8,000 x g,
2 min**

6 Elute RNA

Elute the RNA in **60 µL RNase-free H₂O**, (supplied) and centrifuge at **8,000 x g** for **1 min**.

If higher RNA concentrations are desired, elution can be done with 40 µL. Overall yield, however, will decrease when using smaller volumes.

For further alternative elution procedures see section 2.4.



**+ 60 µL
RNase-free
H₂O**



**8,000 x g,
1 min**

5.2 RNA isolation from up to 10⁵ cells

Before starting the preparation:

- Check if Wash Buffer RA3 was prepared according to section 3.

1 Sample preparation

As sample material use **up to 10⁵ cells** in a volume of up to **100 µL**.



Fill up sample to 100 µL (e.g. with PBS)

2 Cell lysis

Add **300 µL Buffer RA1** and vortex vigorously in order to lyse the cells.

+ 300 µL RA1

Vortex

3 Adjust RNA binding conditions

Add **300 µL ethanol (96 – 100 %)** to the lysate and mix by vortexing or pipetting up and down.



+ 300 µL ethanol (96–100%)

Mix

After addition of ethanol a stringy precipitate may become visible which will not affect the RNA isolation. Be sure to mix thoroughly and apply sample as homogeneous solution onto the column.

4 Bind RNA

For each preparation, take one NucleoSpin® RNA Binding Column (light blue) placed in a Collection Tube and load the lysate (700 µL).



Load lysate

Centrifuge for **30 s** at **8,000 x g**. Discard Collection Tube with flow-through and place the column in a new Collection Tube.

8,000 x g, 30 s



Maximal loading capacity of NucleoSpin® RNA Binding Columns is 750 µL. Repeat the procedure if larger volumes are to be processed.

5 Wash and dry silica membrane

1st wash

Add **250 µL Buffer RA2** to the NucleoSpin® RNA Binding Column. Centrifuge for **30 s** at **8,000 x g**. Discard flow-through and reuse Collection Tube

+ 250 µL RA2


**8,000 x g,
30 s**

2nd wash

Add **700 µL Buffer RA3** to the NucleoSpin® RNA Binding Column. Centrifuge for **30 s** at **8,000 x g**. Discard flow-through and reuse Collection Tube.

+ 700 µL RA3

**8,000 x g,
30 s**



3rd wash

Add **350 µL Buffer RA3** to the NucleoSpin® RNA Binding Column. Centrifuge for **2 min** at **8,000 x g**.

+ 350 µL RA3

Transfer the NucleoSpin® RNA Binding Column to a nuclease-free Collection Tube (1.5 mL, supplied). Open the lid of the column and let the membrane dry for 3 min.

**8,000 x g,
2 min**

If for any reason, the liquid level in the Collection Tube has reached the NucleoSpin® RNA Binding Column after centrifugation, discard flow-through and centrifuge again.

The procedure ensures complete removal of ethanol from the column.

6 Elute RNA

Elute the RNA in **60 µL RNase-free H₂O**, (supplied) and immediately centrifuge at **8,000 x g** for **1 min**.

**+ 60 µL
RNase-free
H₂O**

If higher RNA concentrations are desired, elution can be done with 40 µL. Overall yield, however, will decrease when using smaller volumes.



**8,000 x g,
1 min**

For further alternative elution procedures see section 2.4.

6 Appendix

6.1 Troubleshooting

Problem	Possible cause and suggestions
RNA is degraded/no RNA obtained	<p data-bbox="331 343 552 365"><i>RNase contamination</i></p> <ul data-bbox="331 379 983 533" style="list-style-type: none"> • Create an RNase-free working environment. Wear gloves during all steps of the procedure. Change gloves frequently. Use of sterile, disposable polypropylene tubes is recommended. Keep tubes closed whenever possible during the preparation. Glassware should be oven-baked for at least 2 hours at 250 °C before use.
Poor RNA quality or yield	<p data-bbox="331 572 754 595"><i>Reagents not applied or restored properly</i></p> <ul data-bbox="331 609 983 727" style="list-style-type: none"> • Sample and reagents have not been mixed completely. Always vortex vigorously after each reagent has been added. • No ethanol has been added. Binding of RNA to the silica membrane is only effective in the presence of ethanol. <p data-bbox="331 770 443 793"><i>Kit storage</i></p> <ul data-bbox="331 807 983 922" style="list-style-type: none"> • Store kit components at room temperature. Storage at low temperatures may cause salt precipitation. • Keep bottles tightly closed in order to prevent evaporation or contamination. <p data-bbox="331 965 499 987"><i>Sample material</i></p> <ul data-bbox="331 1002 983 1129" style="list-style-type: none"> • Sample material not stored properly. Whenever possible, use fresh material. If this is not possible, flash freeze the samples in liquid N₂. Samples should always be kept at -70 °C. Never allow tissues to thaw before addition of lysis buffer. Perform disruption of samples in liquid N₂.
Contamination of RNA with genomic DNA	<ul data-bbox="331 1169 983 1326" style="list-style-type: none"> • The NucleoSpin® RNA Clean-up procedure does not comprise a DNA digestion step. Therefore the extent of DNA contamination mainly depends on the sample material. If lowest level of DNA contamination is desired, use one of the rDNase containing NucleoSpin® RNA kits (see ordering information).

Problem**Possible cause and suggestions**

Suboptimal performance of RNA in downstream experiments	<p><i>Carry-over of ethanol or salt</i></p> <ul style="list-style-type: none">• Do not let the flow-through touch the column outlet after the second wash using Wash Buffer RA3. Be sure to centrifuge at the corresponding speed for the respective time in order to remove ethanolic Wash Buffer RA3 completely.• Check if Wash Buffer RA3 has been equilibrated to room temperature before use. Washing at lower temperatures lowers efficiency of salt removal by Wash Buffer RA3.• A 2 min centrifugation with a subsequent 3 min drying with open lid is sufficient for an extensive removal of ethanol from the column. Residual ethanol will typically be around 1%. Increasing the drying step with open lid from 3 min to 20 min will decrease the residual ethanol content commonly to below 0.1%, but also RNA recovery will be reduced 5–20%. <p><i>Store isolated RNA properly</i></p> <ul style="list-style-type: none">• Eluted RNA should always be kept on ice for optimal stability since trace contaminations of omnipresent RNases (general lab ware, fingerprints, dust) will degrade the isolated RNA. For short term storage freeze at -20 °C, for long term storage freeze at -70 °C. <p><i>RNA concentration is too low</i></p> <ul style="list-style-type: none">• For highest RNA concentration and most sensitive downstream applications, NucleoSpin® RNA Clean-up XS is recommended. NucleoSpin® RNA Clean-up XS allows elution in only 5–20 µL volume (see ordering information).
Higher RNA yield than theoretically possible	<ul style="list-style-type: none">• If performing clean-up of samples containing less than approximately 300 ng, RNA subsequent quantification by A_{260} measurement may simulate yields larger than the RNA input. This may be due to absorbance of silica abrasion. In order to prevent incorrect A_{260} quantification of small RNA amounts, centrifuge the elution tube for 30 s at 8.000–11.000 x <i>g</i> and withdraw an aliquot for measurement without disturbing any sediment or use a silica abrasion insensitive RNA quantification method (e.g., RiboGreen fluorescent dye).

6.2 Ordering information

Product	REF	Pack of
NucleoSpin® RNA Clean-up	740948.10	10 preps
	740948.50	50 preps
	740948.250	250 preps
NucleoSpin® RNA Clean-up XS	740903.10	10 preps
	740903.50	50 preps
	740903.250	250 preps
NucleoSpin® RNA II	740955.20	20 preps
	740955.50	50 preps
	740955.250	250 preps
NucleoSpin® RNA XS	740902.10	10 preps
	740902.50	50 preps
	740902.250	250 preps
NucleoSpin® RNA L	740962.20	20 preps
NucleoSpin® RNA Blood	740200.10	10 preps
	740200.50	50 preps
NucleoSpin® miRNA	740971.10	10 preps
	740971.50	50 preps
	740971.250	250 preps
NucleoSpin® RNA Plant	740949.20	20 preps
	740949.50	50 preps
	740949.250	250 preps
NucleoSpin® FFPE RNA	740969.10	10 preps
	740969.50	50 preps
	740969.250	250 preps
NucleoSpin® 8 RNA	740698	12 x 8 preps
	740698.5	60 x 8 preps
NucleoSpin® 96 RNA	740709.2	2 x 96 preps
	740709.4	4 x 96 preps
	740709.24	24 x 96 preps
NucleoMag® 96 RNA	744350.1	1 x 96 preps
	744350.4	4 x 96 preps

* DISTRIBUTION AND USE OF NUCLEOSPIN® TRIPREP IN THE USA IS PROHIBITED FOR PATENT REASONS.

Product	REF	Pack of
NucleoSpin® TriPrep*	740966.10	10 preps
	740966.50	50 preps
	740966.250	250 preps
NucleoSpin® RNA/Protein	740933.10	10 preps
	740933.50	50 preps
	740933.250	250 preps
Buffer RA1	740961	50 mL
	740961.500	500 mL
rDNase Set	740963	1 set

Visit www.mn-net.com for more detailed product information.

6.3 Product use restriction/warranty

NucleoSpin® RNA Clean-up kit components are intended, developed, designed, and sold FOR RESEARCH PURPOSES ONLY, except, however, any other function of the product being expressly described in original MACHEREY-NAGEL product leaflets.

MACHEREY-NAGEL products are intended for GENERAL LABORATORY USE ONLY! MACHEREY-NAGEL products are suited for QUALIFIED PERSONNEL ONLY! MACHEREY-NAGEL products shall in any event only be used wearing adequate PROTECTIVE CLOTHING. For detailed information please refer to the respective Material Safety Data Sheet of the product! MACHEREY-NAGEL products shall exclusively be used in an ADEQUATE TEST ENVIRONMENT. MACHEREY-NAGEL does not assume any responsibility for damages due to improper application of our products in other fields of application. Application on the human body is STRICTLY FORBIDDEN. The respective user is liable for any and all damages resulting from such application.

DNA/RNA/PROTEIN purification products of MACHEREY-NAGEL are suitable for *IN-VITRO*-USES ONLY!

ONLY MACHEREY-NAGEL products specially labeled as IVD are also suitable for *IN-VITRO*-diagnostic use. Please pay attention to the package of the product. *IN-VITRO*-diagnostic products are expressly marked as IVD on the packaging.

IF THERE IS NO IVD SIGN, THE PRODUCT SHALL NOT BE SUITABLE FOR *IN-VITRO*-DIAGNOSTIC USE!

ALL OTHER PRODUCTS NOT LABELED AS IVD ARE NOT SUITED FOR ANY CLINICAL USE (INCLUDING, BUT NOT LIMITED TO DIAGNOSTIC, THERAPEUTIC AND/OR PROGNOSTIC USE).

No claim or representations is intended for its use to identify any specific organism or for clinical use (included, but not limited to diagnostic, prognostic, therapeutic, or blood banking). It is rather in the responsibility of the user or - in any case of resale of the products - in the responsibility of the reseller to inspect and assure the use of the DNA/RNA/protein purification products of MACHEREY-NAGEL for a well-defined and specific application.

MACHEREY-NAGEL shall only be responsible for the product specifications and the performance range of MN products according to the specifications of in-house quality control, product documentation and marketing material.

This MACHEREY-NAGEL product is shipped with documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. MACHEREY-NAGEL's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Supplementary reference is made to the general business terms and conditions of MACHEREY-NAGEL, which are printed on the price list. Please contact us if you wish to get an extra copy.

There is no warranty for and MACHEREY-NAGEL is not liable for damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; defects in products or

components not manufactured by MACHEREY-NAGEL, or damages resulting from such non-MACHEREY-NAGEL components or products.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, REPRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS.

In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, indirect, incidental, compensatory, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications. This warranty is exclusive and MACHEREY-NAGEL makes no other warranty expressed or implied.

The warranty provided herein and the data, specifications and descriptions of this MACHEREY-NAGEL product appearing in MACHEREY-NAGEL published catalogues and product literature are MACHEREY-NAGEL's sole representations concerning the product and warranty. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agent or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

Product claims are subject to change. Therefore please contact our Technical Service Team for the most up-to-date information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Applications mentioned in MACHEREY-NAGEL literature are provided for informational purposes only. MACHEREY-NAGEL does not warrant that all applications have been tested in MACHEREY-NAGEL laboratories using MACHEREY-NAGEL products. MACHEREY-NAGEL does not warrant the correctness of any of those applications.

Last updated: 07/2010, Rev. 03

Please contact:
MACHEREY-NAGEL GmbH & Co. KG
Tel.: +49 24 21 969-270
tech-bio@mn-net.com