DILUMAT S[®]

MANUEL D'UTILISATION
USER'S MANUAL

TABLE DES MATIERES

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ENGLISH

CONTENTS

1. PREFACE	43
3. UNPACKING	45
3.1 Inventory.....	45
3.2 Assembly	46
4. DESCRIPTION	47
4.1 Front	47
4.2 Back	47
4.3 Peristaltic Pump	48
5. SETUP	49
5.1 Installation.....	49
5.2 Power	49
5.3 Turning On	49
6. USER INTERFACE	50
6.1 Keyboard	50
6.2 Screen	50
6.3 Sounds	51
7. USE	52
7.1 Installing the tube into the pump	52
7.2 Installing the distribution nozzle and tube	53
7.3 Dilution	54
7.4 Dispensing	56
7.5 Interrupting dispensing or dilution	58
7.6 Traceability	58
8. SETTINGS	61
8.1 Time.....	61
8.2 Date	61
8.3 Selecting the operating mode (dilution/dose dispensing).....	62
8.4 Selecting the Operating Mode	62
8.5 Adjustment	63
8.6 Verifying the Weight	65
9. MAINTENANCE	66
9.1 Cleaning	66
9.2 Accessing the Counters.....	66
10. REFERENCES	67
10.1 Options and Accessories.....	67
10.2 Which tube set for which application?	68
11. TROUBLESHOOTING	69
11.1 Symptoms.....	69
11.2 Error Codes	69
12. TECHNICAL CHARACTERISTICS	70

Pictograms used in the manual

This symbol indicates that the message is informative.

This symbol indicates that the message is important.

● Disposal and Recycling

This device is composed of metal, plastic materials, electric components, and electronic components.

The 2002/96/CE directive on recovering and recycling Waste from Electric and Electronic Equipment (WEEE) requires that Electric and Electronic Equipment (EEE) be recycled at the end of its life.

Therefore, end users must dispose of electric and electronic waste at a waste management facility or a dedicated collection site.

● AES Laboratoire® is a registered trademark of AES CHEMUNEX.

● Information Contained in this Manual

Due to the informative nature of the information and diagrams contained in this manual, AES CHEMUNEX reserves the rights to make changes to the document without prior notice.

● Document History

INDEX	DATE	NOTES
A	06-09-08	Creation
B	08-09-08	Modification
C	07-09-10	Modification of chapter 8.5 + replacement of « calibration” by “adjustment”

DECLARATION OF CONFORMITY **CE**

AES CHEMUNEX certifies that the **DILUMAT S®**
(reference: AESAP1082)
complies with the following European directives:

- **Directive 89/336/CEE** (Electromagnetic Compatibility)
Updated by directives 92/31/CEE and 93/68/CEE

Reference Standard:

- Emission

Generic Standard:

NF EN 50 081 – 1 (Ed. June 1992)

Basic Standard:

NF EN 55011 (Ed. 2006) (conduit and radiated)

NF EN 61 000-3-3 (Ed.1995 + A1 Ed 2001)

- Immunity

Generic Standard:

NF EN 50 082-1 (June 1992)

Basic Standard:

NF EN 61 000-4-2 (Static Discharges) Ed 1995 + A1 Ed. 1998 + A2 Ed. 2001

NF EN 61 000-4-3 (Immunity to Radiated Fields) Ed. 2006

NF EN 61 000-4-4 (Fast Transient/Burst) Ed. 2004

NF EN 61 000-4-5 Ed. 2006 (Immunity to Shock Waves)

NF EN 61 000-4-6 Ed. 2007 (Immunity to Interference Induced by Magnetic Fields)

NF EN 61 000-4-11 Ed. 2004 (Immunity to Voltage Dips and Interruptions)

- **Directive 2006/95/CE** (Electric Device Safety)
Updated by directive 93/68 CEE

Reference Standard:

- NF EN 61010-1 (Ed. 2001)

- **Directive 98/37/CEE** (Machine Safety)

Chrystelle Grudet

Quality Manager

1. PREFACE

This diluter is designed to deliver a very high level of performance combined with true simplicity of use.

The DILUMAT S® was developed to meet your expectations. It has undergone extensive testing to ensure quality and reliability.

Its operation and technical characteristics comply with the CE directives. (See the certificate of compliance provided in this manual.)

Beyond its technical abilities, the DILUMAT S® was developed with careful attention to user comfort.

WARRANTY :

AES CHEMUNEX certifies that the device corresponds to the characteristics defined in this manual at the time of its delivery. Machines are guaranteed against manufacturing and component defects for a period of twelve months from the shipping date.

AES CHEMUNEX can repair parts at its facilities during the warranty period. Any malfunctioning part or device should be returned to the AES factory at the expense of the AES representative or the customer, if necessary.

This warranty does not extend to problems due to transportation, an accident, improper use, or internal interferences caused by a bad electric network or individuals who are not authorized to use the device.

2. INTRODUCTION

The DILUMAT S® is an easy-to-use diluter.

Obtaining a mother dilution of solid and semi-solid samples is the first step in microbiological analysis. This step's accuracy therefore affects the outcome of the entire analysis. The manual method requires a long, tedious process of weighing the sample, using many bottles of diluent (ex. one per sample).

The Dilumat S allows this step to be secure, optimized, and automated.

The user can quickly insert the sample quantity. The Dilumat S accurately weighs it and adds exactly the right amount of diluent for the desired dilution factor.

The Dilumat S® can be connected to a larger volume of diluent (glass bottles; sterile pockets or diluent processor).

Produced by AES Chemunex for over 20 years, the Dilumat S® lets you take advantage of the latest innovations to make your job easier.

The LCD screen is user friendly. Its icons help you to understand the various menus and features.

The simplified keyboard allows the device to be operated quickly.

The automatic **rotating arm**,

- enables you to start the dilution automatically, ensuring that the distribution nozzle is protected against manual contamination.
- completely comes apart from the bag's opening and eliminates cross-contamination between samples.

The quick-load **pump** guarantees a reproducible distribution. Its placement in the front of the device makes it easier to insert tubes.

The DILUMAT S® is **compact** so that it can be set up under laminar flow hoods.

The accuracy of the DILUMAT S® fully meets the recommendations under the ISO 7218 and ISO 6887-1 international standards.

The DILUMAT S® can dilute your 25g sample in about 15 seconds.

The solid cover allows for **easy maintenance**, excellent resistance to cleaning products, and durability for your Dilumat S.

The 2nd pump option can easily be installed, if desired.



The printer option allows for greater traceability.

The bar code reader option can be added to increase the simplicity and speed of use.

The DILUMAT S® can be **connected to a PC** to retrieve data.

3. UNPACKING

3.1 Inventory


It is important to save the packaging in case the device ever needs to be moved.

<input checked="" type="checkbox"/>	#	Contents
<input type="checkbox"/>	1	Dilumat S monopump
<input type="checkbox"/>	1	Bag holder
<input type="checkbox"/>	1	CE22 power cord
<input type="checkbox"/>	1	Distribution mast
<input type="checkbox"/>	1	Mast fixing screw
<input type="checkbox"/>	1	1 Box containing <ul style="list-style-type: none"> - 1 Set of distribution tubes - 50 Standard unfiltered bags - 1 Adhesive sleeve - 1 Inspection report - 1 French-English user manual - 1 CD containing user manuals in French, English, German, Italian, and Spanish


Check that the device has not been damaged during transport. Ensure that all of the accessories listed above are included.

3.2 Assembly


4. DESCRIPTION

4.1 Front


Label	Description
1	Bag holder
2	Holding clips for the distribution nozzle
3	Distribution head
4	Distribution nozzles
5	Distribution mast
6	Right peristaltic pump
7	Screen
8	Keyboard
9	Left peristaltic pump (option)

4.2 Back

Label	Description
1	Bar code reader connector
2	Printer or PC connector
3	Voltage selector and fuse box
4	Power socket
5	ON/OFF switch

4.3 Peristaltic Pump

Label	Description
1	Cover
2	Tube fitting
3	Adjustment knob for the internal diameter of the tube
4	Rollers

This equipment is powered up. If you have to inspect the inside, Disconnect the power from the device before handling it.

5. SETUP

5.1 Installation

i Recommendations:

- ✓ The DILUMAT S® must be set up near a power outlet.
- ✓ The diluter must sit on a flat, stable, and level surface.
- ✓ It should be used on a clean and dry work surface.
- ✓ The placement of the device should allow the on/off switch to be accessible.

5.2 Power



	110-120V	220-240V
Voltage	115V +/-10%	230V +/- 10%
Maximum Current	1A	0.5A
Frequency	50-60Hz	50-60Hz
Power Consumed	60 W	60 W
Fuse Rating	1 AT	1 AT

Before plugging in the device, check that the voltage selector on the back is correctly set for the country's power characteristics (triangular cursor pointing to the selection).

i At delivery, the cursor is set up for 230V +/-10%.


5.3 Turning On

! **Before turning on the device, it is recommended that the device remain at room temperature for 30 minutes to stabilize the sensor.**

	<p>Connect the power cord to the DILUMAT S®.</p> <p>i If using a printer, PC, or bar code reader, please check that the peripheral is connected to the Dilumat S® (connection at the back of the device).</p>
	<p>Press the ON/OFF key to turn it on.</p>
	<p>An initialization sequence displays.</p>


6. USER INTERFACE

6.1 Keyboard

Label	Description
1	Settings key
2	Weigh key / verification
3	Menu navigation arrows
4	OK key
5	Cancel key
6	Change Pump key (if option included)

6.2 Screen

The screen contains 3 lines.

The information displayed varies according to the Dilumat S® use sequences (initializing, using, configuring, etc.)

Label	Description
1	Appears when the left pump is being used (if the 2 nd pump option is installed)
2	Appears when the right pump is being used.
3	Counter indicating the weight of the diluent distributed using the current pump.
4	Dilution speed or amount to be dispensed by the current pump.
5	Pictogram / icon indicating the status of Dilumat S® (see the table below)
6	Selected user code
7	Weight
8	Sample ID (if option installed)

Icons

Displayed Icon	Meaning
	User code
	Dilution / dispensing in progress
	Dilution / dispensing paused
	Dilution / dispensing cancelled
	Do not touch the device: allowing for the tare or adjustment in progress
	Operation completed successfully
	Problem / error
	Adjustment sequence in progress
	Weight checking sequence in progress.
	Date setting
	Time setting
	Sets the number of reports to be printed or exported
	Weighing mode
	Dilution mode

6.3 Sounds

	A short beep indicates that the action was completed successfully.
	A long beep indicates a problem or abnormal action.


7. USE

7.1 Installing the tube into the pump


To ensure that the pump operates properly and that diluent is dispensed accurately, **it is essential that the pre-sterilized distribution tube be placed into the device correctly.**

Follow steps 1 through 3 in the diagrammed procedure below:

4. Lift the cover of the pump head.

5. Turn the knob to adjust for the interior diameter of the tube used (**dilution: 6.4mm**).


6. Insert the tube as shown below.
 Next, close the pump cover and fasten the distribution nozzle to the end of the arm.
 Stretch the tube to each side of the rotor, and avoid closing the cover onto the tube.


Direction of pump rotation:

The right pump on the Dilumat S® 4 turns counterclockwise.
 The left pump (if included), turns clockwise.


7.2 Installing the distribution nozzle and tube

7.3 Dilution

<p>9. An initialization screen appears when the device is powered on.</p>	
<p>10. Select a user code (1-99) using the up/down arrows. Then confirm the user code.</p>	
<p>11. The "Ready to weigh" screen appears.</p>	
<p>NOTE: The device retains the most recently used operating mode. To change the operating mode, please refer to section 8.3.</p>	

Accessing the Pre-Dilution Functionalities:

Second pump option


If a second pump has been installed, it is possible to select the distribution pump.

 Press the "Change Pump" key.


- A left or right (< or >) cursor on the screen indicates which pump is currently selected.
- Pressing the key changes the pump.
- There is a counter and a use parameter (dilution factor or volume to dispense) for each pump, displayed at the top of the screen.


Bar code reader option or numeric keypad

The sample's bar code can be scanned during the dilution sequence or input with the keypad. Only the first 16 digits of the bar code are read.

<p>1. Place the bag on the holder. Next, pinch the edges of the support so that the sides of the bag stick to the adhesive.</p>	
<p>2. Measure the weight. The hourglass indicates that nothing should be done.</p>	
<p>3. The "Ready to weigh" screen appears.</p>	

i The keyboard functionalities are identical to those described above for placing the bag onto the holder.


<p>1. Insert the sample into the bag.</p>	
<p>2. The screen displays the weight of the sample. Confirm to start the dilution.</p>	


Dilution completed successfully (accompanied by a short beep)

Error encountered during dilution (accompanied by a long beep)

7.4 Dispensing

<p>4. An initialization screen appears when the device is powered on.</p>	
<p>5. Select a user code (1-99) using the up/down arrows. Then confirm the user code.</p>	
<p>6. The "Ready to dispense" screen appears.</p>	
<p>NOTE: The device retains the most recently used operating mode. To change the operating mode, please refer to section 8.3.</p>	

Accessing the Functionalities:

Second Pump Option


If a second pump has been installed, it is possible to select the distribution pump.

 Press the "Change Pump" key.


- A left or right (< or >) cursor on the screen indicates which pump is currently selected.
- Pressing the key changes the pump.
- There is a counter and a use parameter (dilution factor or volume to dispense) for each pump, displayed at the top of the screen.


Bar Code Reader Option or numeric keypad

The sample's bar code can be scanned during the dilution sequence or input with the keypad. Only the first 16 digits of the bar code are read.

<p>Place the bag on the holder.</p> <p>Next, pinch the edges of the support so that the sides of the bag stick to the adhesive.</p>	
<p>Measure the weight.</p> <p>The hourglass indicates that nothing should be done.</p>	
<p>The "Ready to dispense" screen appears.</p>	

i The keyboard functionalities are identical to those described above for placing the bag onto the holder.

<p>Insert the sample into the bag.</p>	
<p>The screen displays the weight of the sample.</p> <p>Confirm to start dispensing.</p>	

Dispensing completed successfully (accompanied by a short beep)

Error encountered during dispensing (accompanied by a long beep)

7.5 Interrupting dispensing or dilution

A dispensing or dilution operation may be suspended or cancelled using the "Cancel" key.

Suspending the Operation


Pressing the "Cancel" key suspends the dispensing.
There is a short beep.

The screen displays this symbol.
The operation may be resumed by pressing the OK key.

Cancelling the Operation

Pressing the "Cancel" key a second time interrupts the dispensing.
The long beep indicates that the operation has been interrupted.

The screen displays this symbol to indicate the interruption.
The arm returns to its starting position.

7.6 Traceability

A PC or printer connection allows you to publish traceability reports concerning dispensing, dilution, adjustments, and verification of weights taken.

Please refer to the options listed in section 10.

Reports generated using the Zebra printer:

Dilution Report	Dispensing Report
<pre> DILUTION REPORT ----- 01 - DILUMAT S ID : XXXXXXXX (8 characters) 02 - Last calibration : XXXX/XX/XX 03 - Dilution Result : OK ----- 04 - Dilution ID : XXXXXXXXXX (9 characters) 06 - User ID : XX 07 - Date : (YYYY/MM/DD) : XXXX/XX/XX 08 - Time : (24h) : XX:XX:XX 09 - Dilution factor : XX 10 - Sample weight (g) : XX.XX 11 - Dilution weight (g) : XXX.X 12 - Accuracy (%) : XX.X 13 - Sample ID: XXXXXXXXXXXXXXXXXX (16 characters) </pre>	<pre> DISPENSING REPORT ----- 01 - DILUMAT S ID : XXXXXXXX (8 characters) 02 - Last calibration : XXXX/XX/XX 27 - Dispensing Result : OK ----- 05 - Dispensing ID : XXXXXXXXXX (9 characters) 06 - User ID : XX 07 - Date : (YYYY/MM/DD) : XXXX/XX/XX 08 - Time : (24h) : XX:XX:XX 14 - Set weight (g) : XX.XX 15 - Dispensed weight (g) : XXX.X 12 - Accuracy (%) : XX.X 13 - Sample ID : XXXXXXXXXXXXXXXXXX (16 characters) </pre>

Adjustment Report	Check Report
CALIBRATION REPORT ----- 01 - DILUMAT S ID : XXXXXXXXX (8 characters) 02 - Last calibration : XXXX/XX/XX ----- 06 - User ID : XX 07 - Date : (YYYY/MM/DD) : XXXX/XX/XX 08 - Time : (24h) : XX:XX:XX 16 - Set weight (g) : XXXX.XX 17 - Checked weight (g) : XXXX.XX 18 - Accuracy (%) : XX.X 19 - Calibration : OK	CHECK REPORT ----- 01 - DILUMAT S ID : XXXXXXXXX (8 characters) 02 - Last calibration : XXXX/XX/XX ----- 06 - User ID : XX 07 - Date : (YYYY/MM/DD) : XXXX/XX/XX 08 - Time : (24h) : XX:XX:XX 16 - Set weight (g) : XXXX.XX 17 - Checked weight (g) : XXXX.XX 18 - Accuracy (%) : XX.X 20 - Check : OK

Definition of Report Data

01 Dilumat s ID	Serial number for the device
02- Last calibration	Date of last calibration
03- Dilution Result	Result of the dilution
04- Dilution ID	ID for the dilution (counter)
05- Dispensing Id	ID for the dispensing (counter)
06- User ID	ID for the user
07- Date (YYYY/MM/DD)	Date of operation (yyyy/mm/dd format)
08- Time (24H)	Time of operation
09- Dilution factor	Dilution factor
10- Sample weight (g)	Weight of the sample in grams
11- Dilution weight (g)	Final dilution weight in grams
12- Accuracy (%)	Accuracy of the dilution or dispensing
13- Sample ID	Bar code ID for the sample or sample ID (available only with the bar code reader or keypad option)
14- Set weight (g)	Weight of the dispensing set in grams
15- Dispensed weight (g)	Weight of the dispensed diluent in grams
Adjustment and Check Report Data	
16- Set weight (g)	Weight of the calibration set in grams
17- Checked weight (g)	Actual weight of the calibration in grams
18- Accuracy (%)	Accuracy of the calibration or verification
19- Calibration	Conclusion of the calibration
20- Check	Conclusion of the verification
Simulation Report Data	
21- Set cycles	Number of test cycles programmed
22- Completed cycles	Number of test cycles completed
23- Counter 1	Number of dilutions completed since the last reset
24- Counter 2	Number of dispensing operations completed since the last reset
25- Counter 3	Number of hours and minutes powered on
26- Status	Conclusion of the test phase

Note:

- These reports are published only in English.
- KO indicates a failure.

Reports generated with a PC connection:

Dilution Report

Order of data:

[weft id];[Dilumat S ID];[Dilution #];[User ID];
[Last calibration];[Date];[Time];[Dilution factor];[Initial weight];
[Final weight];[Accuracy];[Status];[Sample #]

Data format:

DL;XXXXXXXX;XXXXXXXXXX;XX;XXXX/XX/XX;XXXX/XX/XX;XX:XX:XX;XXX;XX.XX;XXX.X;XX.X;XX;XXXXXXXXXXXXXXXXXX

Dispensing Report

Order of data:

[weft id];[Dilumat S ID];[Dispensing #];[User ID];
[Last calibration];[Date];[Time];[Set weight];[Dispensed weight];
[Accuracy];[Status];[Sample #]

Data format:

DS;XXXXXXXX;XXXXXXXXXX;XX;XXXX/XX/XX;XXXX/XX/XX;XX:XX:XX;XXX.XX;XXX.XX;XX.X;XX;XXXXXXXXXXXXXXXXXX

Calibration Report

Order of data:

[weft id];[Dilumat S ID];[User ID];[Prev. calibration];
[Date];[Time];[Set weight];[checked weight];[Accuracy];[Status]

Data format:

CL;XXXXXXXX;XX;XXXX/XX/XX;XXXX/XX/XX;XX:XX:XX;XXX.XX;XXX.XX;XX.X;XX

Check Report

Order of data:

[weft id];[Dilumat S ID];[User ID];[Prev. check];[Date];
[Time];[Set weight];[checked weight];[Accuracy];[Status]

Data format:

CK;XXXXXXXX;XX;XXXX/XX/XX;XXXX/XX/XX;XX:XX:XX;XXX.XX;XXX.XX;XX.X;XX

SIMUL Test Report

Order of data:



[weft id];[Dilumat S ID];[Date];[Time];[Set cycles];
[Completed cycles];[Counter 1];[Counter 2];[Counter 3];[Status]

Data format:

CK;XXXXXXXX;XXXX/XX/XX;XX:XX:XXX;XXX;XXXXXXXXXX;XXXXXXXXXX;XXXXXX:XX;XX

8. SETTINGS

8.1 Time

 Press the "Settings" key.	
Using the "Navigation" arrows, select the  icon and then press OK.	  
Use the left/right arrows to select the digit to change.	 
Use the up/down arrows to increase or decrease the selected digit.	 
 Confirm by pressing the "OK" key.	

8.2 Date

 Press the "Settings" key.	
Using the "Navigation" arrows, select the  icon and then press OK.	  
Use the left/right arrows to select the digit to change. (yyyy/mm/dd)	 
Use the up/down arrows to increase or decrease the selected digit.	 
 Confirm by pressing the "OK" key.	

8.3 Selecting the operating mode (dilution/dose dispensing)

 Press the "Settings" key.	
Using the "Navigation" arrows, select the  icon and then press OK.	  
Using the "Change Pump" key, select the pump to change.	 
Use the up/down arrows to select the operating mode to use (dilution or dose dispensing).	 
The left/right arrows can enable a dilution factor or amount to distribute (depending on the selected operating mode).	 
Use the up/down arrows to increase or decrease the dilution factor or amount to distribute (depending on the mode selected). Press the OK key.	  

8.4 Selecting the Operating Mode

 Press the "Settings" key.	
Using the "Navigation" arrows, select the  icon and then press OK.	  
Use the up/down arrows to select the operating mode and number of tickets. PRINT 1 to PRINT 9 allow you to select the number of tickets to print. The LIMS selection allows you to export to a PC that is connected to the Dilumat S.	  

8.5 Adjustment

Note : The adjustment must be done cold. No dilution should be done beforehand. Use the same reference weight for the adjustment and checking sequences.


1. Press the "Settings" key.	
2. Using the right/left arrows, select the "Adjustment" icon  and press OK  .	  
3. Make sure that the reference weight is actually placed on the diluter.	
The  pictogram indicates that the device is in "adjustment" mode.	 
4. Press the "Tare" key to start adjustment.	
5. Place the reference weight on the diluter and check that small characters indicate the reference mass. Big characters indicate the measured mass. Example: with a 500g reference mass: 	
<p> <i>When if you realize an adjustment with a 200g modify the small characters using the arrows and change the reference weight's mass. Remember confirming it with the OK key. By default, the standard weight is 500.00g. It is possible to enter a standard weight from 200.00 to 1500.00g. The most recent value entered is kept.</i></p>	
6. Press OK to start the measurement. Wait for the end of the calculation and the  symbol display.	
7. Remove the reference weight from the holder.	
8. Press the "Tare" key  . Wait until the hourglass disappears and the symbol  reappears.	
9. Place the reference weight back on he holder.	

- | |
|--|
| <p>10. Press OK . And wait until the hourglass disappears.</p> |
| <p>11. Check that the  symbol appears, indicating that adjustment is compliant.</p> |
| <p>12. Press once more OK  to go back to the menu and print the adjustment report when a printer is connected to the diluter.</p> |

	<p>"OK" indicates that the verification is valid; weighing error less than 1%.</p>
	<p>The exclamation point indicates that an error occurred during the verification; unstable weight, error greater than 1%, abnormal gap between the standard weight and the weighed amount.</p>
	<p>Press OK to validate the verification and return to the "Ready to dispense / dilute" screen. A verification report is sent to the printer or PC (if one of these options is installed).</p>
	<p>Press the Cancel key to return to the "Settings" screen without saving the new data. A verification report is sent to the printer or PC (if one of these options is installed).</p>

8.6 Verifying weighing

Note : For this verification, the mass used may be different.

<p>From the "Ready to dispense / dilute" screen, press and hold the "Tare" button to access the weighing verification.</p>	
<p>Press the "Tare" key to start the verification.</p>	
<p>The hourglass appears on the left side of the screen to indicate that you must wait while weighing is in progress.</p>	
<p>If needed, use the navigation arrows to change the value of the verification mass.</p> <p>Place the verification weight in the device.</p> <p>Press OK to start measuring the verification mass.</p>	 <p> <i>By default, the standard weight is 500.00g. It is possible to enter a standard weight from 200.00 to 1500.00g. The most recent value entered is kept.</i></p>
<p>The hourglass indicates that you must wait while the diluter memorizes the value.</p>	
	<p>"OK" indicates that the verification is valid; weighing error less than 1%.</p>
	<p>The exclamation point indicates that an error occurred during the verification; unstable weight, error greater than 1%, abnormal gap between the standard weight and the weighed amount.</p>
	<p>Press OK to validate the verification and return to the "Settings" screen. A verification report is sent to the printer or PC (if one of these options is installed).</p>
	<p>Press the Cancel key to return to the "Ready to dispense / dilute" screen without saving the new data. A verification report is sent to the printer or PC (if one of these options is installed).</p>

9. MAINTENANCE

9.1 Cleaning

The DILUMAT S® does not require strict daily maintenance.

- ✓ To be able to work in sterile conditions, you must autoclave the distribution tubes in moist heat.
- ✓ The bag holder and bowl holder must be removed from the device and cleaned daily with a damp cloth and hot water. Do not use a scouring pad on the bag holder.
- ✓ The cover must also be cleaned daily.
- ✓ Cleaning with 70% alcohol should be done weekly.

In the event of a spillage, be careful not to let liquid run over the sensor, even though the bag holder has been specially designed to limit the risk of leaks.

9.2 Accessing the Counters

<p>Press the "Validation" key when powering on the device.</p>	 <div style="background-color: black; color: white; padding: 10px; text-align: center;"> <p>A E S L a b o r a t o i r e D I L U M A T S V 1 . 1 2 b</p> </div>
<p>The following screen appears.</p>	<div style="background-color: black; color: white; padding: 10px; text-align: center; font-size: 2em;"> <p>0 0 0</p> </div>
<p>Press the "Validation" key again to access the counters for the device.</p>	 <div style="background-color: black; color: white; padding: 10px;"> <p>S / N : 3 1 2 3 4 5 6 7 (1) : 0 0 0 0 0 0 1 2 3 (2) : 0 0 0 0 0 0 1 2 3 (3) : 0 0 0 0 0 0 : 0 0</p> </div> <p>S/N = serial number (4) = Number of dilutions performed (5) = Number of dispensings performed (6) = Number of hours and minutes powered on</p>
<p>To return to operating mode, press the "Cancel" key.</p>	 <div style="background-color: black; color: white; padding: 10px; text-align: center;">  2 6 </div>

10. REFERENCES

10.1 Options and Accessories

Reference	Description
AESDI002	Altuglas bag holder
AESDI0043	Bowl holder plate
AESDI0045	Full Altuglas bag holder
AESDI0048	Puradisc 0.2 (50u) filters
AESDI0063	Adhesive
AESDI0098	H=70mm wedge for bowl plate
AESDI0099	H=125mm wedge for bowl plate
AESDI0103	Distribution nozzle holding clips (2u)
AESDI0314	Dilusafe rack for Dilubag
AESDL0038	Set of precut tubes (5u) L: 16mm
AESDL0110	Distribution tube holding clip (2u)
AESDL0300	Second pump kit
AESDL0317	Fuse (10u)
AESDL0324	Distribution nozzle
AESDL0355	Full packaging
AESDL0381	Dilumat S cover
AESDL0394A	Distribution set + tube
AESDL0401	PC connection cable
AESDL0402	DB9/USB Connector Dilumat S
AESDL0403	USB numeric keypad
AESLM0004A	Bar code reader
AESLM0029	Thermal printer labels
AESLM0037	Labels for thermal printer
AESMA108201	User manual
AESMI0388	Right angle power cord
AESMN3928	Silicone tube 6x10 (25m)
AESMN3928A	Pump rotor Silicone tube 6x10 (5m)
AES3500/25	ST STOMACHER 3500/25 bags
AES400/50G	STOMACHER IRRAD /50 bags
AES400FB/25G	Full page filter bags

10.2 Which tube set for which application?




	Distribution Set	Sample Set	Contents
	AESDL0394A	AESDI0062	Bottles Type GL45 (schott)
	AESDL0394A	AESDI0308 + AESDI0062	Bottles AEB611499L (Salt Tryptone) AEB110308M (Peptoned Water) AEB610419L (Half Fraser)
	AESDL0394A	AESAT00117A	Masterclave 09 Self-Processor AESAP01080
	AESDL0394A	AESMB0016	Masterclave 528 Self-Processor AESAP1081
	AESDL0394A	AESDI0313	DILUBAG

11. TROUBLESHOOTING

11.1 Symptoms

Symptom	Possible Cause	Solution(s)
Nothing appears on the screen when the device powers on.	Power problem	Check that the power cable is properly attached. Check the fuse at the back of the device.
The calibration is out of tolerance. Or A weight on the bag holder is not being recognized.	The calibration is wrong.	Turn the device off and on. Perform a calibration (see procedure in section 8.5). If the problem persists, reset the device's program. (Turn off the device and then turn it back on by pressing the Validation and Cancel buttons.) If the problem persists, please contact Customer Service.

11.2 Error Codes

Code	Cause(s)	Solution(s)
	Weight too high	Use a different weight.
	Weight too high	Change the dilution / dispensing settings or the sample size.
	Sample too large	Lighten the sample.
	Incorrect bar code read	Try entering or scanning the bar code again.

12. TECHNICAL CHARACTERISTICS

Power consumed	60 Watts maximum	
Voltage	115V +/- 10%, 50-60Hz or 230V, +/-10%, 50-60Hz	
Fuse rating	1 AT	
Room temperature	15 to 35°C	
Relative humidity	From 10 to 80%	
Atmospheric pressure	700 hPa to 1100 hPa	
Altitude	Up to 2000 meters	
Pollution grade	2	
Installation category	Type II following directive CEI664	
Dimensions of the device	Device alone: (LxHxD) 442 x 455 x 450 mm	Packaged: (LxHxD) 600 x 300 x 60mm
Weight	Device alone: 11 Kg	Packaged: 14.5 Kg
Alphanumeric display: 4x 20 characters, backlit		
9-key membrane keyboard		
Bridge (weights 0-2.5 kg)		
<p>Device accuracy: The Dilumat guarantees the following level of accuracy (in a stable environment): Weighing: <1% starting at 5g Dilution/Dispensing: < 5% of the final weight starting at 10g</p>		
Accuracy on weight		
Dilution / Dispensing : Accuracy on the final weight		

RECOMMENDATIONS ASSOCIATED WITH BIOLOGICAL CONTAMINATION RISKS

Purpose

To protect people involved with analysis systems, instruments and all of their parts must be decontaminated before use. This goal is based on the competence of biologists and the recommendations of providers.

Methods

- ✓ The biologist will decontaminate the device as specified in the user guide within a timely manner and record the decontamination in the maintenance log.
- ✓ The biologist must identify the level of risk associated with their work (HIV, mycobacteria, prions, etc.) and note it in the maintenance log for those outside the laboratory.
- ✓ It is recommended to protect the outside of the instrument during transport (box, film, cover).

Process Application for Instruments

- For external actions involving the instruments (maintenance technician, biomedical engineer, etc.)
- When repairing instruments at the user site before movement by a carrier (returns, moves)

The declaration of decontamination must include:

- ✓ The name of the material and its serial number
- ✓ The individual performing the decontamination
- ✓ The protocol used

It must be signed and dated.

References

Framework directive 89/391/CEE on implementing measures for improving on-the-job health and safety for workers

Directive 90/679/CEE on protecting workers against risks associated with on-the-job exposure to biological agents.

Law dated December 31, 1991 (Article L 230-1 to 5 of the Labor Code)

Decree dated May 4, 1994 (Article R 231-60 to R 231-65.3 of the Labor Code)

Decree dated July 16, 2007 (J.O. dated August 4, 2007)

Declaration of Decontamination

To be faxed or shipped with the device sent for repairs

In an effort to protect our partners against the dangers of contaminated devices, we ask that you complete, sign, and return this declaration to us.

Without this form properly completed, we cannot repair or calibrate the device.

Device:

Serial Number:

We declare:

- That the shipped devices have been carefully rinsed and decontaminated before shipment.
- That the shipped devices present no risk of bacteriological, virological, chemical or radioactive contamination.
- That we are aware that shipping a contaminated device may be a personal liability or a liability of the represented company.

Last Name:

First Name:

Title:

Company:

Date:

Signature:

Corporate Stamp

Phone:

Fax:

Rue Maryse Bastié
Ker Lann – CS17219
35172 BRUZ cedex
FRANCE

Telephone: +33 (0)2 23 50 12 12
Fax: +33 (0)2 23 50 12 00
aeschemunex.com
contact@aeschemunex.com