



LSC

Author:

J. Thomson and D. A. Burns

PerkinElmer, Inc.

Downers Grove, IL

Environmental Sample Preparation for LSC

Introduction

The advance of the nuclear industry in all its forms coupled with growing concerns for

possible environmental contamination has led to an increased interest in the quantification of radioisotopes in the environment. Radionuclides currently present in the environment originate from a variety of natural sources such as fallout from nuclear weapons testing, and in discharges from both nuclear and non-nuclear industries. Natural sources and fallout from nuclear devices provide the main input to terrestrial ecosystems, except for certain radionuclides emanating from nuclear installations. In the late 1950's and early 1960's, during and immediately after the period of most frequent above ground nuclear weapons testing, numerous studies were performed to determine the distribution and movement of fallout radionuclides in air, precipitation, agricultural produce, animals and soils. Around this time the nuclear power industry was developing, and consequently, both the diversity and amount of radioactive species in the environment increased. Currently, the majority of the high level waste from the nuclear industry is stored for ultimate disposal in sites classified as stable, such as deep geological strata. However, as many nuclear facilities are situated in coastal areas, the bulk of low level radioactive waste is discharged to the sea.



Because of both nuclear fallout and discharges from the nuclear industry (including releases from the Chernobyl accident in 1986), certain radionuclides are studied more than others. This is due to either their radio toxicity, increased presence in the environment or ease of entry into the food chain^{1, 2, 3, 4, 5, 6, 7}; a selection of these is presented in Table 1.

Table 1. Environmental radionuclides of specific interest.

| Radionuclide | Environmental Location |
|--|---|
| ³H | Milk, crops, groundwater, and precipitation |
| ¹⁴ C | Milk, crops, animals and sea water |
| ³⁵ S | Milk, crops, animals and soil/sediments |
| ⁹⁰ Sr | Milk, crops, animals and soil/sediments |
| ^{134,137} Cs | Milk, crops, animals and soil/sediments |
| ²³⁸ Pu, ^{239,240} Pu, ²⁴¹ Pu, ²⁴¹ Am | Milk, crops, animals and soil/sediments |

The separation and isolation of these radionuclides from the complex sample matrices often encountered is presenting researchers with a myriad of problems; however many of these have been eased by the introduction of novel chromatographic separation technology (Eichrom Industries Inc. Darien, Illinois, USA)

By employing this new technology, previously difficult and time consuming radionuclide separations are completed more effectively and efficienctly. In combination with advances in liquid scintillation counting (LSC) technology by PerkinElmer, it is now possible to consider LSC as an alternative screening tool to alpha spectrometry and gas flow proportional counting.

Chromatographic Sample Preparation

A range of products that employ the technique of extraction chromatography to efficiently pre-concentrate and separate radionuclides from a variety of different matrices have recently been introduced by Eichrom. Extraction chromatography combines the power and selectivity of solvent extraction with the ease of use of a chromatographic column. Through careful selection of the extractant, which is bound to the resin support, a material is produced which is highly specific for a particular radionuclide or group of radionuclides. In the majority of cases, either retention or elution of the radionuclide of interest is achieved by modification of the concentration of the acidic eluent used.

Table 2. Compatibility of Ultima Gold cocktails with Eichrom eluents.

| Method | Isotope | Emitter | Strip Volume (ml) | Eluent | Recommended Ultima Gold Cocktail |
|----------|----------------------|---------|----------------------|------------------------------|-------------------------------------|
| ACS06 | 234-238 U | α, γ | 15 | 0.02 M HCl | AB/XR/LLT |
| | ²³⁴ Th | β, γ | 20 | 5 M HCl | AB ¹ |
| ACW01 | 234-238 U | α, γ | 20 | 0.02 M HCl | AB/XR/LLT |
| ACVVUT | ¹³⁴ Th | β, γ | 20 | 6 M HCI | AB¹/LLT¹ |
| | ²³⁴⁻²³⁸ U | α, γ | 15 | 0.01 M HCl | AB/XR/LLT |
| ACW03 | ²⁴¹ Pu | β | 10 | 0.1 M Ammonium Oxalate | AB/XR |
| | ²⁴¹ Am | α, γ | 3 + 20 | 9 M + 4 M HCl | AB¹/LLT¹ |
| ACW06 | 234-238 _U | α, γ | 15 | 0.02 M HCl | AB/XR/LLT |
| ACVVUb | ²³⁴ Th | β, γ | 15 | 5 M HCl | AB ¹ |
| OTW01 | ²¹⁰ Pb | β, γ | 20 | Water | LLT |
| SRW01 | ^{89,90} Sr | β | 10 | 0.05 M HNO ₃ | AB/XR |
| SRS01 | ^{89,90} Sr | β | 10 | 0.05 M HNO₃ | AB/XR |
| TCS01 | ⁹⁹ Tc | β | 2 | TEVA Resin | AB/XR/LLT |
| TCW01 | ⁹⁹ Tc | β | 2 | TEVA Resin | AB/XR/LLT |
| SRU01 | ^{89,90} Sr | β | 10 | 0.05 M HNO₃ | AB/XR |
| ACW04 | ²⁴¹ Am | α, γ | 15 | 2 M HCl | AB¹/LLT¹ |
| | ²³⁴⁻²³⁸ U | α, γ | 15 | 0.02 M HCl | AB/XR/LLT |
| ACU02 | ²⁴¹ Pu | β | 15 | 3 M HCI-0.25 M Ascorbic Acid | AB/LLT |
| | ²⁴¹ Am | α, γ | 3 + 20 | 9 M + 4 M HCl | AB¹/LLT¹ |
| A.C.W.00 | ²⁴¹ Pu | β | 10 | 0.1 M Ammonium Oxalate | AB/XR |
| ACW09 | ²⁴¹ Am | α, γ | 3 + 20 | 9 M + 4 M HCl | AB¹/LLT¹ |
| ACW07 | ²⁴¹ Pu | β | 10 | 0.1 M Ammonium Oxalate | AB/XR |
| H3W1 | ³H | β | 25 | Non acidified water sample | LLT |
| H3VVI | ⁶³ Ni | β | 15 | 3 M HNO ₃ | AB¹/LLT¹ |

¹ Indicates limited sample uptake capacity (see Table 3 for further details)

Employing such technology has resulted in an accurate and reproducible sample preparation method for the determination of radionuclides in environmental samples. Eichrom has published methods for radionuclide separations using their chromatography technology and these can be conveniently explained as follows:

| Method | Refer to: |
|---------------|------------------------------------|
| ACW, ACS, ACU | Actinides in Water, Soil and Urine |
| SRW, SRS, SRU | Strontium in Water, Soil and Urine |
| TCW, TCS | Technetium in Water and Soil |
| H3W | Tritium in Water |
| OTW | Other (e.g. Lead in Water) |

A summary of these methods together with the recommended PerkinElmer Cocktail for isotopic determination by LSC is presented in Table 2. Ultima Gold cocktails⁹ are recommended for these applications due to the use of diisopropyl naphthalene (DIN) as the solvent base which enhances alpha/beta resolution in LSC⁷. The information presented in Table 2 shows both the eluent used in the final stripping of the radionuclide from the

chromatographic column and the appropriate Ultima Gold cocktail⁹ which will accommodate either all, or at least a sizeable aliquot of the eluent. Table 2 should be used in conjunction with Table 3 which demonstrates the maximum capacity of each of the Ultima Gold cocktails with each of the eluents. It should be remembered that the use of alpha/beta LSC for alpha determination will only provide a gross alpha measurement and is capable of limited alpha-alpha resolution. An example of such an alpha/beta LSC is the PerkinElmer TriCarb 3180TR (comparable with TriCarb B3180 TR/SL) including alpha/beta separation option which uses Time-Resolved Pulse Decay Analysis (TR-PDA)⁷ technology to separate the alpha spectrum from the beta spectrum.

An alternative method to acidic stripping of the radionuclide from the column is to elute with Isopropyl alcohol (IPA). IPA effectively strips the resin coating (containing the radionuclide) and this can be counted with 4π geometry in a suitable Ultima Gold cocktail. The information in both these tables is presented to help researchers investigate the use of LSC as an alternative technology for gross alpha determination.

Table 3. Sample capacity of Ultima Gold cocktails for chromatographic eluents.

| Eluent | Strip Volume (ml) | ml Sample/10 ml Ultima Gold AB @ 20 °C | ml Sample/10 ml Ultima Gold LLT @ 20 °C | ml Sample/10 ml Ultima Gold XR @ 20 °C |
|------------------------------------|-------------------|---|--|---|
| 0,01 M HCl | 15 | 10.0 | 8.0 | 10.0 |
| 0,02 M HCl | 15-20 | 9.0 | 7.0 | 10.0 |
| 2,0 M HCl | 15 | 3.5 | 3.5 | 1.0 |
| 5,0 M HCl | 15 | 2.0 | 1.5 | <0.5 |
| 6,0 M HCl | 20 | 1.0 | 1.5 | <0.5 |
| 4,65 M HCl (9 M + 4 M Mixture) | 3 + 20 | 1.5 | 2.0 | <0.5 |
| 9,0 M HCl | 20 | 1.0 | 1.0 | <0.25 |
| 3 M HCl/0,25 M Ascorbic acid | 15 | 2.0 | 2.0 | 0.5 |
| 0,05 M HNO ₃ | 10 | 8.0 | 7.0 | 9.0 |
| 3,0 M HNO ₃ | 15 | 2.0 | 2.25 | 1.0 |
| 0,02 M HNO ₃ /0,02 M HF | 10 | 8.0 | 10.0 | 10.0 |
| 0,1 M Ammonium oxalate | 10 | 8.0 | 6.0 | 9.0 |
| Water | 25 | 10.0 | 10.0 | 10.0 |

Aqueous Sample Preparation

Many of the radioactive species of interest to low level researchers are present in an aqueous medium, usually water. Therefore, any suitable LSC cocktail must not only have a high capacity for water but also be compatible with water from a variety of different sources. These include distilled, deionized, tap, rain, river and even sea water. In addition to high sample capacity, other preferred requirements for liquid scintillation counting include a very low background contribution and high counting efficiency. Ideally the LSC cocktail should also be based on the high flash point, safer solvent DIN. Ultima Gold LLT is

such a cocktail and is primarily designed for the low level Tritium monitoring and research sectors. Additionally, Ultima Gold LLT has other unique performance characteristics that set it apart from currently available cocktails. Ultima Gold LLT meets additional requirements for a low level counting cocktail such as long term stability and sub-ambient temperature stability, and it can accept the important mineral acid species normally encountered in alpha/beta counting applications. A selection of the more important properties of Ultima Gold LLT is illustrated in Tables 4 and 5.

Table 4. Suitability of Ultima Gold LLT for low level tritium measurements.

| Water/Cocktail Ratio (ml) | Cocktail | ³H Efficiency (0,4 – 4,5 keV) | Background (CPM) | E²/B | E ² V ² /B | MDA (Bq/l) |
|------------------------------|-----------------|----------------------------------|------------------|------|----------------------------------|------------|
| 8:12 | Ultima Gold LLT | 24.6 | 1.15 | 526 | 33680 | 1.22 |
| 10:10 | Ultima Gold LLT | 21.2 | 1.11 | 405 | 40490 | 1.11 |
| 11:9 | Ultima Gold LLT | 18.1 | 0.95 | 345 | 41730 | 1.06 |

Table 5. Sample capacity of Ultima Gold LLT cocktail for different acids.

| Temperature (°C) | 1 M HCl (ml) | 2 M HCl (ml) | 1 M HNO ₃ (ml) | 2 M HNO ₃ (ml) | 1 M H ₂ SO ₄ (ml) | 2 M H ₂ SO ₄ (ml) | 1 M H ₃ PO ₄ (ml) | 2 M H ₃ PO ₄ (ml) |
|---------------------|-----------------|-----------------|------------------------------|------------------------------|--|--|--|--|
| 14 | 4.0 | 2.0 | 3.0 | 2.25 | 5.5 | 3.5 | 4.5 | 4.0 |
| 16 | 4.0 | 2.5 | 3.5 | 2.25 | 6.5 | 4.0 | 5.0 | 4.0 |
| 18 | 4.0 | 2.5 | 3.5 | 2.25 | 7.0 | 4.0 | 7.0 | 4.0 |
| 20 | 5.0 | 3.0 | 3.5 | 2.5 | 7.0 | 4.0 | 4.0 | 4.0 |

Table 4 shows some of the properties of Ultima Gold LLT measured in Low Level Mode in a TriCarb 2770 TR/SL (comparable with TriCarb 3180 TR/SL) comparable with TriCarb B3180 TR/SL) at 15 °C. All samples were prepared in duplicates in polyethylene vials. The measurement time for the determination of the background was 500 minutes. The minimum detectable activity (MDA) was calculated using the formula from Currie⁸.

Table 5 indicates that Ultima Gold LLT is capable to take up amounts of acids generally used for environmental samples. Ultima Gold LLT also shows a high sample capacity for urine samples, and therefore is ideally suited for bioassays with ³H samples in urine. As a quick guide Table 6 shows the most

important environmental radionuclides from nuclear power plants and the nuclear industry.

Conclusion

Environmental sample preparation encompasses a wide variety of techniques, including extraction chromatography, adic extraction, ashing and solvent extraction, and chemical separation. Advances in liquid scintillation technology, together with new and emerging sample preparation techniques, now enable researchers to consider LSC as an alternative environmental sample radionuclide counting method, or as a potentially useful screening tool.

Table 6. The main radionuclides originating from nuclear devices and the nuclear industry.

| Element | Radionuclide | Half-life | Emitter | Energy (MeV) | Sources |
|------------|------------------|-----------|--------------|----------------------------|---------------------------|
| Hydrogen | ³H | 12.3 y | β | 0.02 | Fallout, nuclear industry |
| Carbon | ¹⁴ C | 5730 y | β | 0.16 | Fallout, nuclear industry |
| Phosphorus | ³² P | 14.3 d | β | 1.71 | Fallout, nuclear industry |
| Sulphur | ³⁵ S | 88 d | β | 0.17 | Nuclear industry |
| Argon | ⁴¹ Ar | 1.83 h | β | 1.20; 2.49 1.29 | Nuclear industry |
| Calcium | ⁴⁵ Ca | 165 d | β | 0.26 | Nuclear industry |
| Chromium | ⁵¹ Cr | 27.8 d | γ EC | 0.32 0.75 | Nuclear industry |
| Manganese | ⁵⁴ Mn | 303 d | γ | 0.84 | Fallout, nuclear industry |
| | ⁵⁵ Fe | 2.6 y | EC | 0.23 | Fallout, nuclear industry |
| Iron | ⁵⁹ Fe | 45.1 d | β | 0.27; 0.47 0.14-1.3 | Nuclear industry |
| Cobalt | ⁵⁸ Co | 71.3 d | β γ EC | 0.47 0.81; 0.86 2.31 | Nuclear industry |
| | ⁶⁰ Co | 5.3 y | β | 0.32; 1.49 1.17; 1.33 | Fallout, Nuclear industry |
| Nickel | ⁶³ Ni | 92 y | β | 0.07 | Nuclear industry |
| Zinc | ⁶⁵ Zn | 244 d | β γ EC | 0.33 1.11 1.11 | Nuclear industry |
| Arsenic | ⁷⁶ As | 1.1 d | β | 0.35-2.96 0.51-2.66 | Nuclear industry |
| Krypton | ⁸⁵ Kr | 10.8 y | β | 0.67 0.14 | Fallout, Nuclear industry |
| Strontium | ⁸⁹ Sr | 52 d | β | 1.46 | Fallout, Nuclear industry |
| SHOHLIUH | ⁹⁰ Sr | 28.1 y | β | 0.55 | Fallout, Nuclear industry |

Table 6. Continued

| Element | Radionuclide | Half-life | Emitter | Energy (MeV) | Sources |
|------------|--------------------|------------------------|-------------------|--------------------------|---------------------------|
| | 90γ | 2.67 d | β | 2.27 | Nuclear industry |
| Yttrium | 91γ | 58.8 d | β γ | 0.33; 1.55 1.21 | Nuclear industry |
| Zirconium | ⁹⁵ Zr | 65 d | β γ | 0.36-1.13 0.24; 0.72 | Fallout, Nuclear industry |
| Niobium | ⁹⁵ Nb | 35.1 d | β γ | 0.16 0.77 | Fallout, Nuclear industry |
| Technetium | ⁹⁹ Tc | 2.12·10⁵ y | β | 0.29 | Fallout, Nuclear industry |
| Ruthenium | ¹⁰³ Ru | 39.6 d | β γ | 0.226 0.04-0.61 | Fallout, Nuclear industry |
| | ¹⁰⁶ Ru | 367 d | β | 0.04 0.51-2.64 | Nuclear industry |
| Silver | ^{110m} Ag | 253 d | β γ | 0.09; 0.53 0.66; 0.82 | Nuclear industry |
| Antimony | ¹²⁴ Sb | 60.3 d | β γ | 0.06-2.32 0.04-2.3 | Nuclear industry |
| Antimony | ¹²⁵ Sb | 2.7 y | β γ | 0.10-0.62 0.04-0.67 | Fallout, Nuclear industry |
| | ^{125m} Te | 58 d | γ | 0.04; 0.11 | Nuclear industry |
| Tellurium | ¹³² Te | 3.25 d | β γ | 0.22 0.05-0.23 | Nuclear industry |
| Iodine | 129 | 1.7·10 ⁷ y | β | 0.19 0.04 | Fallout, Nuclear industry |
| iodine. | 131 | 8.07 d | β | 0.26-0.81 0.08-0.72 | Fallout, Nuclear industry |
| V | ^{131m} Xe | 11.8 d | γ | 0.16 | Nuclear industry |
| Xenon | ¹³³ Xe | 5.27 d | β γ | 0.27; 0.35 0.08; 0.38 | Nuclear industry |
| | ¹³⁴ Cs | 2.05 y | β | 0.09; 0.41 0.48-1.4 | Nuclear industry |
| Cesium | ¹³⁶ Cs | 13 d | β | 0.34; 0.56 0.07; 1.24 | Fallout |
| | ¹³⁷ Cs | 30.2 y | β γ | 0.51; 1.18 0.66 | Fallout, Nuclear industry |
| Barium | ¹⁴⁰ Ba | 12.8 d | β γ | 0.47-1.02 0.14-0.54 | Fallout |
| Lanthanum | ¹⁴⁰ La | 1.67 d | β γ | 1.25-2.17 0.11-2.55 | Fallout, Nuclear industry |
| Cerium | ¹⁴¹ Ce | 33 d | β γ | 0.44; 0.58 0.15 | Fallout, Nuclear industry |
| Cellulli | ¹⁴⁴ Ce | 285 d | β γ | 0.18-0.31 0.03-0.13 | Fallout, Nuclear industry |
| Neodymium | ¹⁴⁷ Nd | 11.1 d | β γ | 0.38; 0.82 0.09-0.69 | Fallout |
| Promethium | ¹⁴⁷ Pm | 2.5 y | β | 0.23 | Nuclear industry |
| Europium | ¹⁵⁴ Eu | 16 y | β | 0.27-1.86 0.06-1.6 | Nuclear industry |
| | ¹⁵⁵ Eu | 1.81 y | β | 0.10-0.25 0.04-0.11 | Nuclear industry |
| Polonium | ²¹⁰ Po | 138.4 d | α γ | 5.30 0.80 | Nuclear industry |
| Thorium | ²³⁴ Th | 24.1 d | β γ | 0.10; 0.19 0.03-0.09 | Nuclear industry |
| | ²³⁴ U | 2.47·10⁵ y | α γ | 4.60; 4.72 0.05-0.58 | Nuclear industry |
| Uranium | ²³⁵ U | 7.1·10 ⁸ y | α γ | 4.16-4.60 0.07-0.37 | Nuclear industry |
| Oranium | ²³⁶ U | 2.39·10 ⁷ y | α γ | 4.33; 4.44 0.05 | Nuclear industry |
| | ²³⁸ U | 4.51·10 ⁹ y | $lpha$ γ | 4.14; 4.15 0.05 | Nuclear industry |

Table 6. Continued

| Element | Radionuclide | Half-life | Emitter | Energy (MeV) | Sources |
|-----------|-------------------|------------------------|---------|--------------------------------|---------------------------|
| Neptunium | ²³⁷ Np | 2.14·10 ⁶ y | α γ | 4.40-4.87 0.02-0.24 | Nuclear industry |
| | ²³⁸ Pu | 86 y | α γ | 5.36; 5.46; 5.50 0.04 | Fallout, nuclear industry |
| | ²³⁹ Pu | 2.44·10⁴ y | α | 5.01-5.16 0.04-0.77 | Fallout, nuclear industry |
| Plutonium | ²⁴⁰ Pu | 6580 y | α | 5.02; 5.12; 5.17 0.05; 0.10 | Fallout, nuclear industry |
| | ²⁴¹ Pu | 13.2 y | α β | 4.80-5.05 0.02 | Fallout, nuclear industry |
| | ²⁴² Pu | 3.79·10⁵ y | α | 4.86; 4.90 | Nuclear industry |
| Americium | ²⁴¹ Am | 458 y | α | 5.39-5.55 0.03-0.06 | Fallout, nuclear industry |
| Curium | ²⁴² Cm | 163 d | α | 5.97-6.11 0.04 | Fallout, nuclear industry |

References

- 1. Ministry of Agriculture, Fisheries and Food. *Radionuclides in Foods, Food Surveillance,* Paper No. 43 (1994), ISBN 0 11 242975 0.
- 2. Ministry of Agriculture, Fisheries and Food. *Directorate of Fisheries Research, Aquatic Environment Monitoring Report,* No. 42 (1994).
- 3. K. L. Bocock; *Radionuclides in Terrestial Ecosystems*, 1981, ISBN 0904282422.
- 4. National Radiological Protection Board. *Radiation Exposure* of the UK Population, 1993, Review (NRPB-R263), HMSO, ISBN 0 85951 364 5.
- Her Majesty's Inspectorate of Pollution. Monitoring Programme

 Radioactive Substances, Report for 1993, HMIP, Cameron
 House, Lancaster LA1 4XO.

- Ministry of Agriculture, Fisheries and Food, Terrestial Radioactivity Monitoring Programme (TRAMP), Report for 1993, Radioactivity food and agricultural products in England and Wales (TRAMP/8).
- C. J. Passo Jr., G. T. Cook; Handbook of Environmental Liquid Scintillation Spectrometry 1994, Packard Instrument Company, Meriden, CT 06450 USA.
- 8. L. A. Currie; Limits of qualitative detection and quantitative determination. *Analytical Chemistry*, 40, 586-593 (1968).
- PerkinElmer(Germany) GmbH, 63110 Rodgau-Jügesheim, LSC Application Note, Cocktails for Liquid Scintillation Counting, October 2004.

The application note information enclosed is used to illustrate the technique and may not represent the latest instrument, reagents and cocktails. Customers should validate the technique in their laboratory. Contact Customer Care at www.perkinelmer.com/contact to find the current PerkinElmer instruments, reagents and cocktails.

PerkinElmer, Inc. 940 Winter Street Waltham, MA 02451 USA P: (800) 762-4000 or (+1) 203-925-4602 www.perkinelmer.com

