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System for the identification of *Listeria*

SUMMARY AND EXPLANATION

API® Listeria is a standardized system for the identification of *Listeria*, which uses miniaturized tests, as well as a database. The complete list of those organisms that it is possible to identify with this system is given in the Identification Table at the end of this package insert.

PRINCIPLE

The API® Listeria strip consists of 10 microtubes containing dehydrated substrates which enable the performance of enzymatic tests or sugar fermentations. During incubation, metabolism produces color changes that are either spontaneous or revealed by the addition of reagents.

The reactions are read according to the Reading Table and the identification is obtained by consulting the profile list in the package insert or using the identification software.

CONTENT OF THE KIT (Kit for 10 tests) :

- 10 API® Listeria strips (**STR**)
- 10 ampules of API® Suspension Medium, 2 ml (**MED**)
- 1 ampule of ZYM B reagent (**ZYM B**)
- 10 incubation boxes (**INCUB**)
- 10 result sheets (**SHEET**)
- 1 Package Insert provided in the kit or downloadable from www.biomerieux.com/techlib

COMPOSITION

Strip

The composition of the API® Listeria strip is given in the Reading Table of this package insert.

Medium

MED (API® Suspension Medium) 2 ml	Demineralized water
--	---------------------

Reagent

ZYM B * reagent 8 ml	Fast Blue BB (active ingredient)	0.12 g
	Methanol	40 ml
	Dimethylsulfoxide (DMSO)	60 ml

The quantities indicated may be adjusted depending on the titer of the raw materials used.

* Signal word : **DANGER**

Hazard statements

- H226** : Flammable liquid and vapour.
- H302** : Harmful if swallowed.
- H311** : Toxic in contact with skin.
- H315** : Causes skin irritation.
- H319** : Causes serious eye irritation.
- H331** : Toxic if inhaled.
- H335** : May cause respiratory irritation.
- H370** : Causes damage to organs.

Precautionary statements

- P210** : Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P280** : Wear protective gloves/protective clothing/eye protection/face protection.
- P301 + P312** : IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
- P302 + P352** : IF ON SKIN: Wash with plenty of water.
- P304 + P340** : IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- P305 + P351 + P338** : IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For more detailed information consult the material safety data sheet.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents

- McFarland Standard (Ref. 70900), No. 1 on the scale
- **apiweb™** identification software (Ref. 40011) (consult bioMérieux)


Material

- Pipettes or PSipettes
- Ampule rack
- Ampule protector
- General microbiology laboratory equipment

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- Refer to the hazard statements "H" and precautionary statements "P" indicated above.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI® M29-A, *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Current revision*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents after the expiry date.
- Before use, check that the packaging and components are intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, etc.

- Allow reagents to come to room temperature before use.
- Open ampules carefully as follows :

- Place the ampule in the ampule protector.
- Hold the protected ampule in one hand in a vertical position (white plastic cap uppermost).
- Press the cap down as far as possible.
- Position the thumb tip on the striated part of the cap and press forward to snap off the top of the ampule.
- Take the ampule out of the ampule protector and put the protector aside for subsequent use.

* For ampule with no dropper-cap :

- Carefully remove the cap.

* For ampule with dropper-cap :

- Turn the ampule upside down and maintain it in a vertical position.
- Squeeze gently on the cap to release a drop to transfer all the reagent into the dropper-bottle.

NOTE: For subsequent use, it is recommended to squeeze on the cap before turning the ampule upside down in order to suck back any spare drops of reagent and thus avoid sprinkling cap or fingers with the reagent.

- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.

STORAGE CONDITIONS

Strips STR :

The strips should be stored at 2-8°C until the expiry date indicated on the packaging.

Media MED :

The media may be stored at 2-30°C until the expiry date indicated on the packaging.

Reagent :

The reagent should be stored in the dark at 2-8°C until the expiry date indicated on the packaging.

The reagent may be kept for up to 1 month after the ampule has been opened and the reagent transferred into the dropper-bottle, (or until the expiration date if this is earlier) : **record the date opened on the bottle label.**

The ZYM B reagent is very sensitive to light : check the appearance of the reagent before transferring it into the dropper-bottle. The ZYM B reagent is normally yellow to amber in color. Dispose of the reagent if any tint of pink (sign of deterioration) is observed. Exposure to laboratory lighting for a short period of time (approx. 1 hour) will damage the reagent.

Make sure that the reagent is put back in the refrigerator immediately after use.

USE OF THE REAGENT

- Open the ampule of reagent as indicated in the paragraph "Warnings and Precautions" (ampule with dropper-cap).
- Before using each ampule of reagent, it is recommended to perform a quality control test with **the strain ATCC® 33090™** indicated in the Quality Control paragraph, in order to eliminate any defective reagents.
- Dispense one drop of reagent.
- Carefully close the bottle after use and store it as indicated in the paragraph "Storage conditions".

SPECIMENS

API® Listeria is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Selection of the colonies

- Check that the strain belongs to the genus *Listeria* (short Gram-positive rods which are polymorphic, motile at 25°C but not at 37°C, catalase-positive and oxidase-negative).
- As specimens often contain a mixture of several types of *Listeria*, it is preferable to make a subculture on blood agar using a well-isolated colony. Incubate the plate for 24 hours at 36°C ± 2°C.

NOTE : The following media may be used to take colony samples before using the API Listeria strip :

- non-selective blood agar media with a Columbia or TSA base and with or without antibiotics ;
- selective media for *Listeria*, excluding the McBride medium which inhibits the enzymatic manifestation of bacteria on the API® Listeria strip. If this medium is used for isolation purposes, make a subculture on blood agar.

Preparation of the strip

- Prepare the incubation box (tray and lid) and distribute about 3 ml of distilled water or demineralized water [or any water without additives or chemicals which may release gases (e.g. Cl₂, CO₂, etc.)] into the honeycombed wells of the tray to create a humid atmosphere.
- Record the strain reference on the elongated flap of the tray. (Do not record the reference on the lid as it may be misplaced during the procedure).
- Remove the strip from its individual packaging.
- Place the strip in the incubation box.
- Discard the desiccant.

Preparation of the inoculum

- Open an ampule of API® Suspension Medium (2 ml) as indicated in the paragraph "Warnings and Precautions" (ampule with no dropper-cap).
- Using a pipette or PSlipette, pick up a few well-isolated colonies. It is recommended to use young cultures (18-24 hours old).
- Prepare a suspension with a turbidity equivalent to **1 McFarland**. This suspension must be used immediately after preparation.
- Observe the type of hemolysis and record it on the result sheet. This qualifies as an additional test.

Inoculation of the strip

- Distribute the bacterial suspension into each tube, avoiding the formation of bubbles (tilt the strip slightly forward and place the tip of the pipette or PSlpette against the side of the cupule) :
 - Fill tube and cupule of the [DIM] test (about 100 µl), **avoiding the formation of a convex meniscus**.
 - Only fill the tube part of tests ESC to TAG (about 50 µl).

NOTE : The quality of the filling is very important : tubes which are insufficiently or excessively full may cause false positive or false negative results.

- Close the incubation box.
- Incubate for 18-24 hours at 36°C ± 2°C in aerobic conditions.

READING AND INTERPRETATION

Reading the strip

- Add a drop of ZYM B reagent to the [DIM] test.
- Read all the reactions within 3 minutes by referring to the Reading Table in this package insert.
- Record whether the reactions are positive or negative (+/-) on the result sheet.
- Also record the type of hemolysis on the result sheet. The result of this test is not taken into account for interpretation of the strip.

Interpretation

Identification is obtained with the **numerical profile**.

- Determination of the numerical profile :

On the result sheet, the tests are separated into groups of three and a value 1, 2 or 4 is assigned to each. By adding together the values corresponding to positive reactions within each group, a 4-digit numerical profile is obtained which constitutes the numerical profile.

- Identification :


This is performed using the database (V1.2)

* with the numerical profile :

- Look up the profile in the list in this package insert ; this list is not exhaustive, so if the profile is missing, please consult the software program indicated below or the bioMérieux Technical Assistance service.

* with the **apiweb™** identification software:

- Enter the 4-digit numerical profile manually via the keyboard.

QUALITY CONTROL

The media, strips and reagents are systematically controlled at various stages of their manufacture.

Streamlined quality control may be used to confirm acceptable performance of the API® LISTERIA system after shipping/storage. This methodology may be performed by following the instructions above for testing and meeting the criteria stated in CLSI® M50-A Quality Control for Commercial Microbial Identification Systems.

Testing may be conducted using *Listeria innocua* ATCC® 33090™ to evaluate the performance of the DIM and XYL tests. Testing performed by bioMérieux has shown that the DIM and XYL tests are the most labile on the API® LISTERIA strip. When testing the strip, *Listeria innocua* ATCC® 33090™ can be used to detect degradation.

For those users who are required to perform **comprehensive quality control** testing with the strip, the following three strains should be tested to demonstrate positive and negative reactivity for most of the API® LISTERIA tests.

1. *Listeria innocua* ATCC® 33090™
2. *Listeria ivanovii* ATCC® BAA-139™
3. *Listeria monocytogenes* ATCC® 19115™

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

	[DIM]	ESC	αMAN	DARL	XYL	RHA	MDG	RIB	G1P	TAG
1.	+	+	+	+	-	+	+	-	-	-
2.	+	+	-	+	+	-	V	-	+	-
3.	-	+	+	+	-	+	+	-	-	-

* Occasionally (-).

Profiles obtained after culture on Columbia sheep blood agar.

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

- The API® Listeria system is intended uniquely for the identification of bacteria belonging to the genus *Listeria* included in the database (see Identification Table at the end of this package insert). It cannot be used to identify any other microorganisms or to exclude their presence.
- Only pure cultures of a single organism should be used.

RANGE OF EXPECTED RESULTS

Consult the Identification Table at the end of this package insert for the range of expected results for the various biochemical reactions.

PERFORMANCE

643 collection strains and strains of various origins belonging to species included in the database were tested :

- 98.91 % of the strains were correctly identified (with or without supplementary tests).
- 0.62 % of the strains were not identified.
- 0.47 % of the strains were misidentified.

WASTE DISPOSAL

Dispose of unused ZYM B reagent following procedures for hazardous chemical waste.

Dispose of all used or unused reagents (other than the ZYM B reagent) as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

READING TABLE

TESTS	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS	RESULTS	
				NEGATIVE	POSITIVE
DIM	Enzymatic substrate	0.106	Differentiation <i>L. innocua</i> / <i>L. monocytogenes</i>	<u>ZYM B / < 3 min</u>	
				pale orange pink beige grey beige	orange
ESC	Esculin Ferric citrate	0.16 0.024	hydrolysis (ESCulin)	pale yellow	black
αMAN	4-nitrophenyl-αD-mannopyranoside	0.045	α-MANnosidase	colorless	yellow
DARL	D-ArabitOL	0.4	acidification (D-ARabitOL)	red / orange-red	yellow / yellow-orange
XYL	D-Xylose	0.4	acidification (XYLose)		
RHA	L-Rhamnose	0.4	acidification (RHAmnose)		
MDG	Methyl-αD-glucopyranoside	0.4	acidification (Methyl-αD-Glucopyranoside)		
RIB	D-Ribose	0.4	acidification (RIBose)		
G1P	Glucose-1-Phosphate	0.4	acidification (Glucose-1-Phosphate)		
TAG	D-Tagatose	0.4	acidification (TAGatose)		

- The quantities indicated may be adjusted depending on the titer of the raw materials used.
- Certain cupules contain products of animal origin, notably peptones.

PROCEDURE	p. I
LIST OF NUMERICAL PROFILES	p. II
IDENTIFICATION TABLE	p. II
LITERATURE REFERENCES	p. III
INDEX OF SYMBOLS	p. IV

REVISION HISTORY

Change type categories

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
Note:	<i>Minor typographical, grammar, and formatting changes are not included in the revision history</i>

Release date	Part Number	Change Type	Change Summary
2015/03	07887Q	Administrative	REVISION HISTORY CONTENT OF THE KIT, SPECIMENS, INDEX OF SYMBOLS
		Technical	COMPOSITION, WARNINGS AND PRECAUTIONS

BIOMERIEUX, the blue logo, API and **apiweb** are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute Inc.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

Any other name or trademark is the property of its respective owner.

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