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Identification system for coryneform bacteria

SUMMARY AND EXPLANATION

API Coryne is a standardized system for the identification of coryneform bacteria in 24 hours, which uses miniaturized tests and a specially adapted database. The complete list of those organisms that it is possible to identify with this system is given in the Identification Table at the end of this package insert.

PRINCIPLE

The API Coryne strip consists of 20 microtubes containing dehydrated substrates for the demonstration of enzymatic activity or the fermentation of carbohydrates.

The enzymatic tests are inoculated with a dense suspension of organisms, which reconstitutes the enzymatic substrates. During incubation, metabolism produces color changes that are either spontaneous or revealed by the addition of reagents.

The fermentation tests are inoculated with an enriched medium (containing a pH indicator) which reconstitutes the sugar substrates. Fermentation of carbohydrates results in acidification which is detected by a spontaneous color change in the pH indicator.

The reactions are read according to the Reading Table and the identification is obtained by referring to the Analytical Profile Index or using the identification software.

CONTENT OF THE KIT (Kit for 12 tests) :

- 12 API Coryne strips
- 12 ampules of API GP Medium
- 12 ampules of API Suspension Medium, 3 ml
- 1 McFarland ampule, point 6
- 12 result sheets
- 12 incubation boxes
- 1 package insert

COMPOSITION

Strip

The composition of the API Coryne strip is given in the Reading Table of this package insert.

Media

API GP Medium 2 ml	L-cystine	0.5 g
	Tryptone (bovine/porcine origin)	20 g
	Sodium chloride	5 g
	Sodium sulfite	0.5 g
	Phenol red	0.17 g
	Demineralized water	to make 1000 ml
pH : 7.4 - 7.8		
API Suspension Medium 3 ml	Demineralized water	
McFarland Standard 6	BaSO ₄	2.88 · 10 ⁻⁴ mol/l

The quantities indicated may be adjusted depending on the titer of the raw materials used.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents

- Reagents : NIT 1 + NIT 2 (Ref. 70 442)
 - ZYM A (Réf. 70 494)
 - ZYM B (Réf. 70 493)
 - PYZ (Réf. 70 492)
- Mineral oil (Ref. 70 100)
- Hydrogen peroxide (3 %)
- DENSIMAT (Ref. 99 234)
- API Coryne Analytical Profile Index (Ref. 20 990) or **apiweb**™ identification software (Réf. 40 011) (consult bioMérieux)
- Columbia blood agar with CNA (Ref. 43 071) or without CNA (Ref. 43 041) or Trypcase Soy + blood (Ref. 43 001)


Material

- Pipettes or PSipettes
- Ampule rack
- Ampule protector
- Swabs
- General microbiology laboratory equipment

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI M29-A, *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Current revision*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiry date.
- Before use, check that the packaging and components are intact.
- Do not use strips which have been damaged : cupules deformed, etc.
- It is recommended to perform a quality control test when a new ampule of ZYM B reagent is opened.

- Open ampules carefully as follows :
 - Place the ampule in the ampule protector.
 - Hold the protected ampule in one hand in a vertical position (white plastic cap uppermost).
 - Press the cap down as far as possible.
 - Position the thumb tip on the striated part of the cap and press forward to snap off the top of the ampule.
 - Take the ampule out of the ampule protector and put the protector aside for subsequent use.
 - Carefully remove the cap.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.

STORAGE CONDITIONS

The strips and media should be stored at 2-8°C until the expiry date indicated on the packaging.

The ampules of API Suspension Medium may be stored at 2-30°C until the expiry date indicated on the packaging.

SPECIMENS (COLLECTION AND PREPARATION)

API Coryne is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Selection of the colonies

Once the strain to be identified has been isolated and verified to be a member of the group of Gram-positive, non-sporeforming, facultatively aero-anaerobic rods :

- Note the type of hemolysis.
- Pick up a well-isolated colony and suspend it in 0.3 ml of sterile water.
- Flood an agar plate (Trypcase Soy agar + 5% sheep blood or Columbia agar + 5% sheep blood with or without CNA) with this suspension (or aseptically swab the entire surface of the agar).
- Incubate the plate for 24 - 48 hours at 37°C.

Preparation of the strip

- Prepare an incubation box (tray and lid) and distribute about 5 ml of distilled water or demineralized water [or any water without additives or chemicals which may release gases (e.g. Cl₂, CO₂, etc.)] into the honey-combed wells of the tray to create a humid atmosphere.
- Record the strain reference on the elongated flap of the tray. (Do not record the reference on the lid as it may be misplaced during the procedure.)
- Remove the strip from its individual packaging.
- Place the strip in the incubation box.

Preparation of the inoculum

- Open an ampule of API Suspension Medium as indicated in the paragraph "Warnings and Precautions".
- Using a swab, harvest all the culture from the previously prepared subculture plate. It is recommended to use young cultures (24-48 hours old).
- Make a dense suspension with a turbidity greater than **6 McFarland**. This suspension must be used immediately after preparation.

Inoculation of the strip

- In the first eleven tests of the strip (tests NIT to GEL) distribute the suspension obtained, avoiding the formation of bubbles (tilt the strip slightly forwards and place the tip of the pipette or PSIPette against the side of the cupule) :
 - For the tests NIT to ESC : distribute approximately 100-150 µl into each cupule.
 - For the URE test : fill the tube only.
 - For the GEL test : fill both tube and cupule.
- In the last nine tests of the strip (tests 0 to GLYG) :
 - Open an ampule of API GP Medium as indicated in the paragraph "Warnings and Precautions" and transfer the rest of the suspension into it (approx. 0.5 ml). Mix well.
 - Distribute this new suspension into the tubes only.
- Fill the cupule of the underlined tests (URE, and 0 to GLYG) with mineral oil to form a slight convex meniscus.
- Place the lid on the tray.
- Incubate for 24 hours (± 2 hours) at 36°C ± 2°C in **aerobic** conditions.

READING AND INTERPRETATION

Reading the strip

After the incubation period :

- Add the reagents :
 - NIT test : 1 drop each of NIT 1 and NIT 2
 - PYZ test : 1 drop of PYZ
 - PyrA, PAL, βGUR, βGAL, αGLU, βNAG tests : 1 drop each of ZYM A and ZYM B (*).
- (*) **It is recommended to control** each ampule of ZYM B before using for the first time.
To do this, it is recommended to use **the strain ATCC® 27402** indicated in the Quality Control paragraph in order to eliminate any defective reagents.
- Wait 10 minutes, then read all the reactions according to the Reading Table. If necessary, expose the strip to a strong light (1000 W) for 10 seconds to decolorize any excess reagents in tubes PyrA to βNAG.
- Perform the CATalase test (used as the 21st test) : add 1 drop of hydrogen peroxide (3 %) to the ESC or GEL test. Wait 1 minute. The appearance of **bubbles** corresponds to a **positive** reaction.
- Record all the reactions on the result sheet.

Interpretation

Identification is obtained with the **numerical profile**.

- Determination of the numerical profile :
On the result sheet, the tests are separated into groups of 3 and a value of 1, 2 or 4 is indicated for each. By adding together the values corresponding to positive reactions within each group, a 7-digit profile number is obtained which constitutes the numerical profile. The catalase reaction constitutes the 21st test and has a value of 4 if it is positive.

READING TABLE

TESTS	ACTIVE INGREDIENTS	QTY (mg/cup.)	REACTIONS/ENZYMES	RESULTS	
				NEGATIVE	POSITIVE
NIT	potassium nitrate	0.136	reduction of NITrates	NIT 1 + NIT 2 / 10 min colorless very pale pink	
PYZ	pyrazinocarboxamide	0.56	PYraZinamidase	PYZ / 10 min colorless very pale brown very pale orange	
PYRA	pyroglutamic acid- β -naphthylamide	0.0256	PYRrolidonyl Arylamidase	ZYM A + ZYM B (PyrA \rightarrow β NAG) / 10 min colorless pale orange	
PAL	2-naphthyl-phosphate	0.0244	ALkaline Phosphatase	colorless beige - pale purple pale orange	purple
β GUR	Naphthol ASBI-glucuronic acid	0.0548	β -GlucURonidase	colorless pale grey pale beige	blue
β GAL	2-naphthyl- β D-galactopyranoside	0.0312	β -GALactosidase	colorless beige - pale purple	purple
α GLU	2-naphthyl- α D-glucopyranoside	0.0308	α -GLUcosidase	colorless beige - pale purple pale green	purple
β NAG	1-naphthyl-N-acetyl- β D-glucosaminide	0.0348	N-Acetyl- β -Glucosaminidase	colorless beige - pale purple pale brown pale grey	brown
ESC	esculin ferric citrate	0.546 0.078	β -glucosidase (ESCulin)	colorless grey	black
URE	urea	0.76	UREase	yellow orange	red pink
[GEL]	gelatin (bovine origin)	0.6	Hydrolysis (GELatin)	no diffusion of black pigment	diffusion of black pigment
0	Negative control	-	Fermentation	red orange	yellow yellow - orange
GLU	D-glucose	1.56	Fermentation (GLUcose)		
RIB	D-ribose	1.4	Fermentation (RIBose)		
XYL	D-xylose	1.4	Fermentation (XYLose)		
MAN	D-mannitol	1.36	Fermentation (MANnitol)		
MAL	D-maltose	1.4	Fermentation (MALtose)		
LAC	D-lactose (bovine origin)	1.4	Fermentation (LACtose)		
SAC	D-saccharose (sucrose)	1.32	Fermentation (SACcharose)		
GLYG	glycogen	1.28	Fermentation (GLYcoGen)		
CAT	(ESC or [GEL] test)	-	CATalase	H ₂ O ₂ (3 %) / 1 min no bubbles bubbles	

- The quantities indicated may be adjusted depending on the titer of the raw materials used.
- Certain cupules contain products of animal origin, notably peptones.

PROCEDURE	p. I
IDENTIFICATION TABLE	p. II
LITERATURE REFERENCES	p. III
INDEX OF SYMBOLS	p. IV

bioMérieux, the blue logo, API and **apiweb** are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.
ATCC is a trademark belonging to American Type Culture Collection.
Any other name or trademark is the property of its respective owner.

bioMérieux SA
RCS LYON 673 620 399
69280 Marcy-l'Etoile / France
Tél. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
www.biomerieux.com

bioMérieux, Inc
Box 15969,
Durham, NC 27704-0969 / USA
Tél. (1) 919 620 20 00
Fax (1) 919 620 22 11
Printed in France

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