

ID 32 C

IVD

Identification system for yeasts

SUMMARY AND EXPLANATION

ID 32 C is a standardized system for the identification of yeasts, which uses 32 miniaturized assimilation tests and a database. The complete list of those yeasts that it is possible to identify with this system can be found in the Identification Table at the end of this package insert.

Reading and interpretation are carried out automatically or manually.

PRINCIPLE

The ID 32 C strip consists of 32 cupules, each containing a dehydrated carbohydrate substrate.

A semi-solid, minimal medium is inoculated with a suspension of the yeast organism to be tested. After 24-48 hours of incubation, growth in each cupule is read either using the ATB™ Expression™ or *mini API*® instruments, or visually.

Identification is obtained using the identification software.

CONTENT OF THE KIT (Kit for 25 tests)

- 25 ID 32 C strips
- 25 incubation lids
- 25 ampules of API C Medium
- 1 package insert provided in the kit or downloadable from www.biomerieux.com/techlib

COMPOSITION

Strip

The composition of the ID 32 C strip is given below in the list of tests :

CUPULES	TESTS	SUBSTRATES	QTY (mg/cup.)
1.0	GAL	D-GALactose	0.70
1.1	ACT	cycloheximide (ACTidione)	0.014
1.2	SAC	D-SACcharose (sucrose)	0.66
1.3	NAG	N-Acetyl-Glucosamine	0.64
1.4	LAT	LacTic acid	0.64
1.5	ARA	L-ARAbinose	0.70
1.6	CEL	D-CELLobiose	0.66
1.7	RAF	D-RAFFinose	2.34
1.8	MAL	D-MALtose	0.70
1.9	TRE	D-TREhalose	0.66
1.A	2KG	potassium 2-KetoGluconate	1.09
1.B	MDG	Methyl- α -D-Glucopyranoside	1.92
1.C	MAN	D-MANnitol	0.68
1.D	LAC	D-LACtose (bovine origin)	0.70
1.E	INO	INOsitol	0.70
1.F	0	No substrate	-
0.0	SOR	D-SORbitol	2.72
0.1	XYL	D-XYLose	0.70
0.2	RIB	D-RIBose	0.70
0.3	GLY	GLYcerol	0.82
0.4	RHA	L-RHAMnose	0.68
0.5	PLE	PaLatinosE	0.66
0.6	ERY	ERYthritol	1.44
0.7	MEL	D-MELibiose	0.66
0.8	GRT	sodium GlucuRonaTe	0.76
0.9	MLZ	D-MeLeZitose	0.66
0.A	GNT	potassium GlucoNaTe	0.92
0.B	LVT	levulinic acid (LeVulinaTe)	0.48
0.C	GLU	D-GLUcose	0.78
0.D	SBE	L-SorBosE	0.70
0.E	GLN	GLucosamiNe	0.68
0.F	ESC	ESCulin ferric citrate	0.28 0.069

- The numbers indicated correspond to those printed on the strip.
- The quantities indicated may be adjusted depending on the titer of the raw materials used.

Medium

API C Medium		
7 ml	Ammonium sulfate	5 g
	Monopotassium phosphate	0.31 g
	Dipotassium phosphate	0.45 g
	Disodium phosphate	0.92 g
	Sodium chloride	0.1 g
	Calcium chloride	0.05 g
	Magnesium sulfate	0.2 g
	L-Histidine	0.005 g
	L-Tryptophan	0.02 g
	L-Methionine	0.02 g
	Gelling agent	0.5 g
	Vitamin solution	1 ml
	Trace elements	10 ml
	DeminerIALIZED water	to make 1000 ml
	Final pH : 6.4-6.8 at 20-25°C	

- Although API® C Medium contains gelling agent, it **requires no prior heating** and may be as easily pipetted as a liquid medium. It is preferable to warm it at room temperature a few hours before use. **Do not shake.**
- The quantities indicated may be adjusted depending on the titer of the raw materials used.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

Reagents / Instrumentation

- API Suspension Medium, 2 ml (Ref. 70 700)
- DENSIMAT (Ref. 99 234) or ATB Densitometer or McFarland Standard (Ref. 70 900), point 2 on the scale
- ATB Expression or *mini API*, or *apiweb*™ identification software (Ref. 40 011) (consult bioMérieux)
- ATB Electronic Pipette (consult bioMérieux) or ATB Inoculator and Tips (Ref. 15 710)


Material

- Pipettes or PSipettes
- Ampule rack
- Ampule protector
- Air-tight box
- General microbiology laboratory equipment

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).

- All specimens, yeast cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling yeasts should be observed throughout this procedure. Refer to "CLSI® M29-A, *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline - Current revision*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiry date.
- Before use, check that the packaging and components are intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, etc.
- Open ampules carefully as follows :
 - Place the ampule in the ampule protector.
 - Hold the protected ampule in one hand in a vertical position (white plastic cap uppermost).
 - Press the cap down as far as possible.
 - Position the thumb tip on the striated part of the cap and press forward to snap off the top of the ampule.
 - Take the ampule out of the ampule protector and put the protector aside for subsequent use.
 - Carefully remove the cap.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed, particularly the antimicrobial susceptibility patterns.



STORAGE CONDITIONS

The strips and media should be stored at 2-8°C until the expiry date indicated on the packaging.

SPECIMENS (COLLECTION AND PREPARATION)

ID 32 C is not for use directly with clinical or other specimens.

The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Preparation of the strip

- Remove the strip from its packaging.
- Discard the desiccant.
- Place the lid on the strip.
- Record the strain reference on the elongated flap of the strip. (Do not record the reference on the lid as it may be misplaced during the procedure).

Preparation of the inoculum

- Open an ampule of API® Suspension Medium (2 ml) as indicated in the paragraph "Warnings and Precautions" or use any tube containing 2 ml of distilled water without additives.
- Remove 1 or several identical colonies from the culture medium. It is recommended to use young cultures (24-48 hours old).

- Prepare a suspension with a turbidity equivalent to 2 McFarland : measure with the ATB™ Densitometer, or the DENSIMAT, or compare with a turbidity control (McFarland Standard).

NOTE : If the strip is to be read AUTOMATICALLY, the ATB Densitometer or the DENSIMAT must be used to adjust the turbidity of the yeast suspension.

- Open an ampule of API C Medium as indicated in the paragraph "Warnings and Precautions" and transfer approximately 250 µl of the previous suspension into the ampule. This suspension must be used immediately after preparation.

Inoculation of the strip

- AUTOMATIC inoculation :
 - Place the strip, the inoculated ampule of API C Medium and a Tip on the ATB Inoculator tray.
 - The inoculator will automatically homogenize the contents of the ampule and fill the cupules (135 µl / cupule).
- MANUAL inoculation :
 - Homogenize the inoculated ampule of API C Medium and dispense 135 µl of the suspension into each cupule of the strip using the ATB Electronic Pipette.
- Place the lid on the strip.
- Incubate at 29°C ± 2°C for 24-48 hours.

NOTE : Some ventilated incubators may completely dehydrate the medium in the cupules. In this case, place the strip in an air-tight box with a receptacle containing a small volume of water. A humid atmosphere is thus created, which prevents the tests from drying out.

READING AND INTERPRETATION

Reading the strip

- AUTOMATIC reading using ATB Expression™ or *mini API*:
 - check that the middle part of the strip is clean so that the reader can recognize the strip code,
 - check that the name printed on the strip corresponds to the strip name displayed by the software.
 The reader searches for detectable growth in each cupule and transmits the information to the computer.
- VISUAL reading :
 - Compare each cupule to the control (0) and record as positive any cupule that is more turbid.

Interpretation

Identification is obtained using the database (V3.0) :

- AFTER AUTOMATIC READING :
 - the results transmitted to the computer are interpreted by the ATB Expression or *mini API* identification software.
- AFTER VISUAL READING :
 - the reactions obtained are coded into a **numerical profile** :
 - On the result sheet, the tests are separated into groups of 3 and a number 1, 2 or 4 is indicated for each. The values corresponding to positive reactions are then added together within each group.

METHODOLOGIE / PROCEDURE / METHODIK / TECNICA / PROCEDIMENTO / ΔΙΑΔΙΚΑΣΙΑ / METOD / METODE / METODYKA

Tests morphologiques
Morphology tests
Morphologische Tests
Tests morfológicos
Tests morfológici
Testes morfológicos
Εξετάσεις μορφολογίας
Morfologiska tester
Morfologisk test
Ocena morfologii



API® Suspension Medium 2 ml



250 µl



API C Medium



32 x 135 µl




ID 32 C



24:00 - 48:00



29°C ± 2°C



ID 32 C



+ - + - + -



2 McF

**TABLEAU D'IDENTIFICATION / IDENTIFICATION TABLE / PROZENTTABELLE / TABLA DE IDENTIFICACION / TABELLA DI IDENTIFICAZIONE /
QUADRO DE IDENTIFICAÇÃO / ΠΙΝΑΚΑΣ ΤΑΥΤΟΠΟΙΗΣΗΣ / IDENTIFIERINGSTABELL / IDENTIFIKATIONSTABEL / TABELA IDENTYFIKACYJNA**

% de réactions positives après 24/48 H à 29°C ± 2°C / % of positive reactions after 24/48 hrs. at 29°C ± 2°C / % der positiven Reaktionen nach 24/48 h bei 29°C ± 2°C /
% de las reacciones positivas después de 24/48 H a 29°C ± 2°C / % di reazioni positive dopo 24/48 ore a 29°C ± 2°C / % de reacções positivas após 24/48 H a 29°C ± 2°C /
% θετικών αντιδράσεων μετά από 24/48 ώρες στους 29°C ± 2°C / % positiva reaktioner efter 24/48 timmar vid 29°C ± 2°C / % positive reaktioner efter 24/48 timer ved 29°C ± 2°C /
% pozytywnych reakcji po 24/48 godzinach w 29°C ± 2°C









ID 32 C	V3.0	GAL	ACT	SAC	NAG	LAT	ARA	CEL	RAF	MAL	TRE	2KG	MDG	SOR	XYL	RIB	GLY	RHA	PLE	ERY	MEL	GRT	MLZ	GNT	LVT	MAN	LAC	INO	GLU	SBE	GLN	
<i>Candida albicans 1</i>		98	99	100	100	96	0	0	1	100	97	100	98	99	98	1	10	0	100	0	0	0	0	2	1	100	1	1	98	1	99	
<i>Candida albicans 2</i>		100	100	0	50	0	0	0	0	99	67	100	0	67	67	0	0	0	33	0	0	0	0	0	0	99	0	0	100	0	50	
<i>Candida boidinii</i>		5	100	0	86	75	0	0	0	0	0	0	0	99	99	100	99	0	0	100	0	0	0	0	1	100	0	0	100	1	67	
<i>Candida catenulata</i>		100	74	1	67	75	0	0	0	33	33	0	0	74	2	0	91	0	0	0	0	0	0	17	0	100	0	0	100	0	74	
<i>Candida colliculosa</i>		26	19	78	0	82	0	15	78	31	74	80	31	72	4	4	59	4	33	0	4	0	27	31	0	78	1	0	96	1	0	
<i>Candida dattila</i>		43	0	100	0	0	0	14	100	99	100	1	100	99	20	0	86	0	100	0	0	0	83	0	0	71	1	0	100	33	0	
<i>Candida dubliniensis</i>		100	100	100	90	10	0	0	0	100	11	99	0	100	3	0	0	0	1	0	0	0	0	0	0	100	0	0	100	0	60	
<i>Candida famata</i>		100	13	100	97	43	83	91	93	100	99	99	100	100	75	2	99	38	100	66	19	40	100	52	1	100	31	0	100	75	75	
<i>Candida glabrata</i>		0	0	1	0	3	0	0	0	1	99	0	0	0	1	0	30	1	0	0	1	3	0	31	0	0	0	1	100	0	0	
<i>Candida globosa</i>		0	0	60	33	0	0	0	60	50	0	75	60	80	0	0	60	0	60	0	0	0	0	60	0	60	0	0	60	60	25	
<i>Candida guilliermondii</i>		100	53	100	97	30	99	99	100	100	99	97	97	92	99	11	99	3	97	0	72	0	97	9	2	94	0	0	100	85	97	
<i>Candida hellenica</i>		100	100	100	100	33	67	100	67	99	100	0	0	100	100	67	99	100	99	0	0	99	0	67	0	100	0	100	100	100	100	
<i>Candida holmii</i>		88	50	88	0	0	0	0	88	0	75	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	88	0	0
<i>Candida inconspicua/norvegensis</i>		1	0	0	0	94	0	0	0	0	0	0	0	6	0	0	91	0	0	0	0	0	0	0	1	6	0	0	99	0	86	
<i>Candida intermedia</i>		99	0	100	100	50	10	85	85	100	100	100	67	100	99	5	10	50	100	0	1	50	100	85	33	100	75	0	100	85	99	
<i>Candida kefyr</i>		99	100	99	1	99	72	12	99	1	1	5	0	85	83	1	51	0	0	0	1	0	1	0	0	77	96	0	100	1	1	
<i>Candida krusei</i>		1	3	1	90	99	0	0	0	1	0	0	0	1	1	0	95	1	1	1	1	0	0	3	1	1	1	1	100	1	5	
<i>Candida lambica</i>		0	0	0	96	100	0	0	0	0	0	0	0	4	97	0	89	0	0	0	0	0	0	0	0	0	0	0	100	0	89	
<i>Candida lipolytica</i>		0	100	0	97	97	0	0	0	0	0	0	0	54	0	11	100	0	0	100	0	0	0	77	0	84	0	0	100	0	3	
<i>Candida lusitanae</i>		94	1	100	91	14	5	99	7	100	100	91	93	99	94	11	82	100	100	5	0	0	100	84	1	100	0	0	100	95	86	
<i>Candida magnoliae</i>		0	0	98	0	0	0	0	75	0	0	63	0	100	0	20	100	0	0	0	0	0	0	99	0	100	0	0	100	83	0	
<i>Candida melibiosica</i>		100	0	100	50	1	0	100	50	100	100	100	2	100	99	50	50	35	100	0	50	1	100	89	1	100	0	0	100	50	50	
<i>Candida membranifaciens</i>		100	0	100	99	50	100	100	100	100	99	100	80	100	100	99	100	0	100	100	99	20	100	100	1	100	20	0	100	60	60	
<i>Candida norvegica</i>		0	0	0	0	83	0	80	0	0	0	0	0	50	67	0	83	33	0	0	0	0	0	0	1	83	0	0	83	0	0	
<i>Candida parapsilosis</i>		100	25	100	100	1	96	2	1	100	96	93	98	100	96	1	93	1	100	0	0	1	99	92	70	100	1	0	100	72	97	
<i>Candida pelliculosa</i>		43	0	100	0	96	0	70	99	100	97	0	96	97	85	9	100	0	96	90	0	1	96	19	0	100	0	0	100	0	0	
<i>Candida pulcherrima</i>		99	0	100	100	0	0	83	0	100	99	100	99	100	83	80	100	20	100	0	0	0	100	83	0	100	0	0	100	100	99	
<i>Candida rugosa</i>		75	0	1	36	75	1	0	0	0	0	0	0	99	75	0	67	0	0	0	0	0	0	18	0	99	0	0	99	64	30	
<i>Candida sake</i>		80	26	90	90	30	0	28	2	90	85	80	31	85	74	20	83	0	85	2	0	1	70	30	5	85	0	0	85	67	70	

ID 32 C	V3.0	GAL	ACT	SAC	NAG	LAT	ARA	CEL	RAF	MAL	TRE	2KG	MDG	SOR	XYL	RIB	GLY	RHA	PLE	ERY	MEL	GRT	MLZ	GNT	LVT	MAN	LAC	INO	GLU	SBE	GLN	
<i>Candida silvicola</i>		40	100	90	60	50	75	100	0	75	100	1	20	100	80	100	100	99	75	60	0	0	75	40	0	100	0	0	100	1	60	
<i>Candida sphaerica</i>		88	96	100	4	96	0	51	92	67	74	4	61	98	43	0	90	0	65	0	0	0	64	4	0	96	92	0	100	64	4	
<i>Candida tropicalis</i>		99	30	100	100	50	4	91	9	100	99	99	99	100	99	21	9	3	100	0	0	0	99	33	7	100	5	0	100	5	95	
<i>Candida utilis</i>		0	0	100	0	100	0	100	100	100	99	0	57	25	88	0	100	0	100	0	0	0	100	75	0	50	0	0	100	0	0	
<i>Candida valida</i>		0	0	0	79	14	0	0	0	0	0	0	0	0	0	0	62	0	0	0	0	0	0	0	0	0	0	0	0	100	0	86
<i>Candida zeylanoides</i>		0	46	0	99	0	0	0	0	0	53	87	0	97	0	0	100	0	0	0	0	0	0	1	0	100	0	0	100	53	50	
<i>Cryptococcus albidus</i>		40	0	91	0	0	91	91	82	91	82	91	64	91	91	10	0	70	82	0	0	82	91	50	0	55	30	50	82	10	0	
<i>Cryptococcus curvatus</i>		100	17	100	100	75	17	100	75	50	67	100	33	83	100	100	99	33	50	67	0	99	50	100	0	60	100	67	100	0	83	
<i>Cryptococcus humicola</i>		99	71	100	92	100	86	100	43	100	100	100	100	100	100	99	99	100	99	100	100	100	50	100	71	99	100	99	100	71	99	
<i>Cryptococcus laurentii</i>		100	20	100	75	25	100	100	100	99	100	77	100	88	99	40	99	100	60	100	100	80	99	0	99	100	99	100	20	60		
<i>Cryptococcus neoformans</i>		99	0	100	95	0	75	65	81	100	75	100	100	100	81	71	0	91	100	52	0	100	96	76	0	100	0	99	100	50	64	
<i>Cryptococcus terreus</i>		83	0	0	100	0	100	100	0	1	58	100	0	99	100	99	0	67	0	0	0	99	1	100	0	100	85	69	100	91	33	
<i>Cryptococcus uniguttulatus</i>		0	0	100	100	0	99	0	99	100	100	100	100	99	100	0	33	33	100	0	0	100	100	67	0	99	0	100	100	0	99	
<i>Debaryomyces etchellsii/carsonii</i>		90	0	100	91	36	18	50	0	100	64	99	100	100	60	12	60	0	100	12	0	0	99	12	0	100	0	0	100	91	30	
<i>Debaryomyces polymorphus</i>		100	98	100	100	1	50	99	100	100	100	100	100	100	95	25	40	1	100	60	60	0	100	63	0	100	60	0	100	99	99	
<i>Geotrichum capitatum</i>		31	92	0	0	33	0	0	0	0	0	0	0	0	0	0	100	0	0	0	0	0	0	0	0	0	0	0	92	25	0	
<i>Geotrichum spp</i>		90	100	0	0	94	5	0	0	0	0	0	0	99	99	0	99	0	0	0	0	0	0	0	1	95	0	0	100	99	0	
<i>Kloeckera apis/apiculata</i>		0	99	0	0	0	0	100	0	0	0	33	0	0	0	0	0	0	0	0	0	0	0	99	0	0	0	0	100	0	0	
<i>Kloeckera japonica</i>		0	100	0	0	0	0	75	0	0	74	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	85	0	0	
<i>Kodamaea ohmeri</i>		99	0	100	100	0	0	99	99	100	100	100	100	100	0	20	100	0	100	0	0	0	0	0	0	100	0	0	100	80	99	
<i>Pichia farinosa</i>		99	0	1	100	0	0	20	0	0	99	0	0	99	40	40	100	0	0	99	0	0	0	75	1	100	0	0	100	0	50	
<i>Rhodotorula glutinis</i>		80	9	99	0	0	61	28	92	84	97	30	25	55	70	68	66	9	74	0	0	3	84	28	0	55	0	0	99	20	0	
<i>Rhodotorula minuta</i>		8	0	86	57	14	86	71	0	3	86	97	0	43	86	29	100	0	0	0	0	99	86	99	1	43	8	0	99	1	0	
<i>Rhodotorula mucilaginosa</i>		95	0	100	0	0	85	15	100	72	85	0	0	38	90	92	50	5	71	0	0	0	71	26	0	57	0	0	100	8	0	
<i>Saccharomyces cerevisiae</i>		63	1	89	4	65	1	7	89	91	34	0	28	2	0	2	5	2	28	1	5	0	15	4	0	2	2	2	100	1	2	
<i>Saccharomyces kluyverii</i>		80	0	100	40	60	0	5	100	80	75	70	60	87	20	0	0	0	87	0	75	20	25	40	0	99	1	20	100	20	0	
<i>Sporobolomyces salmonicolor</i>		20	0	30	0	0	0	0	30	0	85	0	0	85	0	0	0	0	0	0	0	0	0	85	0	85	0	0	85	0	0	
<i>Stephanoascus ciferrii</i>		99	100	100	100	99	100	75	99	100	100	99	0	100	100	99	100	100	75	100	75	50	0	100	0	99	0	100	100	100	100	
<i>Trichosporon inkin</i>		95	50	100	99	100	0	100	0	100	97	100	100	1	100	100	50	0	100	70	1	75	95	100	1	97	100	97	100	25	70	
<i>Trichosporon asahii</i>		99	99	100	100	100	98	100	0	100	95	100	99	30	100	100	75	75	100	98	1	98	70	100	1	17	100	25	100	17	98	
<i>Trichosporon mucoides</i>		100	99	100	100	100	100	99	99	100	100	100	99	99	100	100	99	99	100	99	100	100	94	100	1	99	100	100	100	94	99	
<i>Williopsis saturnus</i>		0	0	100	0	90	0	100	100	20	20	0	20	60	100	0	99	20	20	0	0	0	15	70	0	20	0	0	100	0	0	
<i>Zygosaccharomyces spp</i>		17	1	1	0	6	0	1	6	6	1	1	1	65	1	0	59	1	6	0	1	1	6	1	0	63	0	0	99	0	6	

**BIBLIOGRAPHIE / LITERATURE REFERENCES / LITERATUR /
BIBLIOGRAFIA / ΑΝΑΦΟΡΕΣ ΑΡΘΡΟΓΡΑΦΙΩΝ / REFERENSLITTERATUR /
LITTERATURHENVISNINGER / PIŚMIENICTWO**

1. BARNETT J.H., PAYNE R.W., YARROW D.
Yeasts : Characteristics and Identification
(1983) Cambridge University Press, London.
2. DEAK T., BEUCHAT L.R.
Comparison of the SIM, **API 20 C** and **ID 32 C** Systems for
Identification of Yeasts Isolated from Fruit Juice Concentrates
and Beverages.
(1993) Journal of Food Protection, 56, 585-592.
3. GUTIERREZ J., MARTIN E., LOZANO C., CORONILLA J.,
NOGALES C.
Evaluation of the ATB 32 C, Automicrobic system and
API 20 C using clinical yeast isolates.
(1994) Ann. Biol. Clin., 50, 443-446.
4. KREGER VAN RIJ N.J.W.
The Yeasts : A Taxonomic Study.
(1984) Elsevier, Amsterdam.
5. McGINNIS M.R. and al.
Taxonomic and Nomenclatural Evaluation of the genera
Candida and *Torulopsis*
(1984) J. Clin. Microbiol., 20, 813-814.
6. MONGET D., CANIAUX I., DESMONCEAUX M.,
GUICHERD M.
ATB 32 C : A New Automated Method for the Identification of
Yeasts.
(1987) Florence, Fifth International Symposium On Rapid
Methods and Automation in Microbiology and Immunology,
4-6 Nov. 1987.
7. WARREN N.G., SHADOMY H.J.
Yeasts of medical importance
in : BALOWS A., HAUSLER W.J., HERMANN K.L.,
ISENBERG H.D., SHADOMY H.J.
Fifth edition,
(1991) Manual of Clinical Microbiology, 617-629.
8. WICKERHAM L.J., BURTON K.A.
Carbon Assimilation Tests for the Classification of Yeasts.
(1948) J. Bact., 56, 363-371.

**TABLE DES SYMBOLES / INDEX OF SYMBOLS / SIMBOLE / CUADRO DE SIMBOLOS /
TABELLA DEI SIMBOLI / QUADRO DE SÍMBOLOS / ΠΙΝΑΚΑΣ ΣΥΜΒΟΛΩΝ /
SYMBOLER / SYMBOLFORTEGNELSE / TABELA SYMBOLI**

Symbole / Symbol Símbolo / Simbolo Σύμβολο	Signification / Meaning / Bedeutung Significado / Significato / Επεξήγηση Betydelse / Betydning / Znaczenie
	Référence du catalogue Catalogue number (GB) / Catalog number (US) Bestellnummer / Número de catálogo / Numero di catalogo Referência de catálogo / Αριθμός καταλόγου Katalognummer / Katalognummer / Numer katalogowy
	Dispositif médical de diagnostic in vitro In Vitro Diagnostic Medical Device In Vitro Diagnostikum Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro Dispositivo médico para diagnóstico in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Medicintekniska produkter för in vitro diagnostik Medicinsk udstyr til in vitro-diagnostik Wyrób do diagnostyki In Vitro
	Fabricant / Manufacturer / Hersteller / Fabricante Fabbicante / Κατασκευαστής / Tillverkare / Producent
	Limites de température / Temperature limitation Temperaturbegrenzung / Limite de temperatura Limiti di temperatura / Limites de temperatura Περιορισμοί θερμοκρασίας / Temperaturbegränsning Temperaturbegrænsning Przestrzeżać zakresu temperatury
	Utiliser jusque / Use by / Verwendbar bis Fecha de caducidad / Utilizzare entro / Prazo de validade Ημερομηνία λήξης / Använd före / Holdbar til / Użyć przed
	Code du lot / Batch code Chargenbezeichnung / Código de lote Codice del lotto / Código do lote Αριθμός Παρτίδας / Lot nummer / Lotnummer Kod partii
	Consulter les instructions d'utilisation Consult Instructions for Use Gebrauchsanweisung beachten Consulte las instrucciones de uso Consultare le istruzioni per l'uso Consulte as instruções de utilização Συμβουλευτείτε τις οδηγίες χρήσης Se handhavandebeskrivningen / Se brugsanvisning Sprawdź w instrukcji obsługi
	Contenu suffisant pour "n" tests Contains sufficient for <n> tests Inhalt ausreichend für <n> Prüfungen Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi Conteúdo suficiente para "n" ensaios Περιεχόμενο επαρκές για «ν» εξετάσεις Räcker till "n" antal tester Indeholder tilstrækkeligt til "n" test Wystarczy na wykonanie <n> testów

FICHE DE RESULTATS / RESULT SHEET / ERGEBNISBLATT / HOJA DE RESULTADOS / SCHEDA PER LA REGISTRAZIONE DEI RISULTATI /
FICHA DE RESULTADOS / ΦΥΛΛΟ ΑΠΟΤΕΛΕΣΜΑΤΩΝ / RAPPORTBLAD / RESULTATARK / KARTA WYNIKÓW

ID 32 C

REF 32 200

Origine / Source / Herkunft / Origen / Origem / Προέλευση / Ursprung / Oprindelse / Pochodzenie

ID32C	1	0	1	2	3	4	5	6	7	8	9	A	B
	0	0	0	0	0	0	0	0	0	0	0	0	0
	1	GAL	ACT	SAC	NAG	LAT	ARA	CEL	RAF	MAL	TRE	2KG	MDG
	0	SOR	XYL	RIB	GLY	RHA	PLE	ERY	MEL	GRT	MLZ	GNT	LVT
	1	2	4	1	2	4	1	2	4	1	2	4	

ID32C	1	0	1	2	3	4	5	6	7	8	9	A	B
	0	0	0	0	0	0	0	0	0	0	0	0	0
	1	GAL	ACT	SAC	NAG	LAT	ARA	CEL	RAF	MAL	TRE	2KG	MDG
	0	SOR	XYL	RIB	GLY	RHA	PLE	ERY	MEL	GRT	MLZ	GNT	LVT
	1	2	4	1	2	4	1	2	4	1	2	4	

0	0	0	0	0
1	C	D	E	
0	0	0	0	0
1	MAN	LAC	INO	
0	GLU	SBE	GLN	
	1	2	4	

0	0	0	0	0
1	C	D	E	F
0	0	0	0	0
1	MAN	LAC	INO	O
0	GLU	SBE	GLN	ESC
	1	2	4	



Incub. / Inkub. / Θερμοκρασία επώασης :
24:00 48:00

Autres tests / Other tests / Andere Tests /
Otras pruebas / Altri test / Outros testes /
Άλλες εξετάσεις / Andra tester /
Andre tests / Inne testy

Ident. / Ταυτοποίηση :



bioMérieux SA
RCS LYON 673 620 399
69280 Marcy-l'Etoile / France
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
www.biomerieux.com

bioMérieux, Inc
Box 15969,
Durham, NC 27704-0969 / USA
Tel. (1) 919 620 20 00
Fax (1) 919 620 22 11

