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Thioglycollate broth with resazurin (THIO-T)

IVD

Sterility testing. Enrichment of non-fastidious bacteria.

SUMMARY AND EXPLANATION

Thioglycollate broth with resazurin (THIO-T) is intended for the detection of anaerobic bacteria but is also used for the detection of aerobic bacteria.

It is recommended for the non-selective enrichment of non-fastidious bacteria in human specimens (4, 5).

This medium complies with the European, American, and Japanese pharmacopoeias (1, 2, 3) for sterility testing.

PRINCIPLE

The broth contains a mixture of peptones which encourage the growth of most bacteria.

The reducing agents (L-cystine and thioglycolic acid) and yeast extract included in the medium, favor the growth of anaerobic bacteria.

The redox indicator (resazurin) reveals the presence of oxygen (pink to mauve color).

CONTENT OF THE KIT**Ready-to-use medium**

REF 42 074 20 x 9 ml tubes

COMPOSITION**Theoretical formula.**

This medium can be adjusted and/or supplemented according to the performance criteria required.

Casein peptone (bovine).....	15 g
L-cystine	0.5 g
Anhydrous dextrose	5 g
Yeast extract	5 g
Sodium chloride	2.5 g
Sodium thioglycollate.....	0.5 g
or Thioglycolic acid	0.3 ml
Resazurine	0.001 g
Agar	0.75 g
Purified water.....	1 l

pH 7.1

MATERIAL REQUIRED BUT NOT PROVIDED

- Bacteriology incubator.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).

- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI M29-A, Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline – Current revision". For additional information on handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories – CDC/NIH – Latest edition", or the current regulations in the country of use.
- Culture media should not be used as manufacturing material or components.
- Do not use reagents after the expiry date.
- Do not use tubes which show signs of contamination.
- The broth may be slightly opalescent but this does not alter the performance.
- Do not use tubes which have a pink-mauve halo greater than 1 cm on the surface (observe on media which has been left to stand).
- Microscopic elements, possibly coming from dead micro-organisms, may be observed in the broth, but this does not alter the performance of the medium.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient's history, the source of the specimen, colonial and microscopic morphology and, if necessary, the results of any other tests performed.
- According to the European Pharmacopoeia, the broth must be handled under Class A conditions, which should be situated within a Class B environment, or in an isolator.
- The outside of the tubes is not sterile and the culture medium may be accidentally contaminated when the tube cap is unscrewed. Disinfect the tube before unscrewing the cap in a controlled environment.

STORAGE CONDITIONS

- **Store the tubes in their box at 2-25°C until the expiry date.**
- **Do not freeze.**

SPECIMENS**For use in medical bacteriology:**

The broth can be inoculated using all types of clinical specimens requiring non-selective enrichment.

For use in industrial bacteriology:

Follow the recommendations of the pharmacopoeias to perform specimen preparation.

INSTRUCTIONS FOR USE

For use in medical bacteriology:

1. Inoculate the specimen to be tested.
2. Incubate at 37°C. The user is responsible for choosing the appropriate temperature for the intended use, in accordance with current standards. The cultures are generally examined after 24 hours of incubation. Incubation time varies according to the type of specimen and the microorganisms being tested for. In certain cases (no growth at 24 hours), it may be necessary to prolong incubation.
3. In case of growth, subculture on an appropriate medium.

For use in industrial bacteriology:

Follow the methods described in the pharmacopoeias.

Note:

- For sterility testing by membrane filtration, it is recommended to use the ready-to-use medium in bottles (Ref. 44 021: see the corresponding technical sheet).

READING AND INTERPRETATION

After incubation, observe the bacterial growth, generally associated with turbidity of the broth.

QUALITY CONTROL

Protocol:

The nutrient capacity of the medium can be tested using the following strain(s):

- *Clostridium sporogenes* ATCC 11437
- *Staphylococcus aureus* ATCC 25923

Range of expected results:

Strain	Results	
<i>Clostridium sporogenes</i> ATCC 11437	Growth after 72 hours at 30-35°C	Appearance of turbidity
<i>Staphylococcus aureus</i> ATCC 25923	Growth after 24 hours at 33-37°C	

Note:

It is the responsibility of the user to perform Quality Control taking into consideration the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature, etc.).

The Thioglycollate broth with resazurin is designed and developed to meet the strictest quality requirements. The results of the strains tested in the batch by batch quality control are given on the quality control certificate available on request.

The quality control is performed according to the recommendations of the standard XP CEN ISO/TS 11133-2.

LIMITATIONS OF THE METHOD

- Growth depends on the requirements of each individual microorganism. It is therefore possible that certain strains which have specific requirements, (substrate, temperature, incubation conditions etc.) may not develop.
- Given the wide variety of specimens tested, it is the responsibility of the user to validate this medium for its specific intended use.

PERFORMANCE

Nutrient capacity:

234 human specimens were studied, including 64 deep wounds, 40 abscesses, 37 pus specimens, 23 wounds, 20 abdominal specimens, 13 drainage specimens, 6 smears, 4 puncture specimens, 3 pleural-type specimens, 24 varied specimens.

177 of the 234 specimens studied showed microbial growth following primary culture and/or enrichment in THIO-T broth. A total of 280 microorganisms were isolated (Gram+ and Gram- aerobes, facultative anaerobes or anaerobes; yeasts).

Number of microorganisms isolated using primary culture and enrichment in THIO-T broth after incubation at 37°C:

Primary culture alone	Primary culture + enrichment
259	280

The combination of primary culture and enrichment in THIO-T broth enabled the recovery of 21 microorganisms more than with primary culture alone, including 7 anaerobic or aero-anaerobic strains.

Organism recovery rate following enrichment in THIO-T at 37°C according to the broth incubation time:

Incubation time	Organism recovery rate
1 day	88%
2 days	96%
3 days	98%
4 days	100%
Total 5 days	100%

WASTE DISPOSAL

Dispose of used and unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazard and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

LITERATURE REFERENCES

1. European Pharmacopoeia 5.
2. USP 28.
3. Japanese Pharmacopoeia 15.
4. BREWER J.M. - Clear liquid mediums for the « aerobic » cultivation of anaerobes - *J.A.M.A.* , 1940, Vol. 115, p. 598-600.
5. VERA H.D. - A comparative study of materials suitable for the cultivation of clostridia - *J. Bacteriol.*, 1944, Vol. 47, p. 59.
6. RUOFF K.L. – Manual of clinical microbiology – 6e Ed., ASM Press, 1995, p. 299-307
7. XP CEN ISO/TS 11133-2 (January 2004) - Microbiology of food and animal feeding stuffs. Guidelines on preparation and production of culture media. Part 2: Practical guidelines on performance testing of culture media.

INDEX OF SYMBOLS

Symbol	Meaning
	Catalogue number
	In vitro diagnostic medical device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult instructions for use
	Contains sufficient for <n> tests

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

Any other name or trademark is the property of its respective owner.

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bioMérieux and the blue logo are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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