

Yeast identification system

**SUMMARY AND EXPLANATION**

API 20 C AUX is a system for the precise identification of the most frequently encountered yeasts. The complete list of those species that it is possible to identify with this system is given in the Identification Table at the end of this package insert.

**PRINCIPLE**

The API 20 C AUX strip consists of 20 cupules containing dehydrated substrates which enable the performance of 19 assimilation tests. The cupules are inoculated with a semi-solid minimal medium. The yeasts will only grow if they are capable of utilizing each substrate as the sole carbon source.

The reactions are read by comparing them to growth controls. Identification is obtained by referring to the Analytical Profile Index or using the identification software.

**CONTENT OF THE KIT (Kit for 25 tests)**

- 25 API 20 C AUX strips
- 25 incubation boxes
- 25 ampules of API C Medium
- 25 result sheets
- 1 package insert

**COMPOSITION****Strip**

The composition of the API 20 C AUX strip is given below in the list of tests:

| TESTS | ACTIVE INGREDIENTS                 | QTY<br>(mg/cup.) |
|-------|------------------------------------|------------------|
| 0     | None                               | -                |
| GLU   | D-GLUcose                          | 1.2              |
| GLY   | GLYcerol                           | 1.2              |
| 2KG   | calcium 2-Keto-Gluconate           | 1.2              |
| ARA   | L-ARAbinose                        | 1.2              |
| XYL   | D-XYLose                           | 1.2              |
| ADO   | ADOnitol                           | 1.2              |
| XLT   | XyLiTol                            | 1.2              |
| GAL   | D-GALactose                        | 1.9              |
| INO   | INOsitol                           | 2.36             |
| SOR   | D-SORbitol                         | 1.2              |
| MDG   | Methyl- $\alpha$ D-Glucopyranoside | 1.2              |
| NAG   | N-Acetyl-Glucosamine               | 1.2              |
| CEL   | D-CELlobiose                       | 1.2              |
| LAC   | D-LACtose<br>(bovine origin)       | 1.2              |
| MAL   | D-MALtose                          | 1.2              |
| SAC   | D-SACcharose (sucrose)             | 1.2              |
| TRE   | D-TREhalose                        | 1.2              |
| MLZ   | D-MeLeZitose                       | 1.2              |
| RAF   | D-RAFFinose                        | 1.9              |

**Medium**

| API C Medium<br>7 ml |                                 |                 |
|----------------------|---------------------------------|-----------------|
|                      | Ammonium sulfate                | 5 g             |
|                      | Monopotassium phosphate         | 0.31 g          |
|                      | Dipotassium phosphate           | 0.45 g          |
|                      | Disodium phosphate              | 0.92 g          |
|                      | Sodium chloride                 | 0.1 g           |
|                      | Calcium chloride                | 0.05 g          |
|                      | Magnesium sulfate               | 0.2 g           |
|                      | L-Histidine                     | 0.005 g         |
|                      | L-Tryptophan                    | 0.02 g          |
|                      | L-Methionine                    | 0.02 g          |
|                      | Gelling agent                   | 0.5 g           |
|                      | Vitamin solution                | 1 ml            |
|                      | Trace elements                  | 10 ml           |
|                      | DeminerIALIZED water            | to make 1000 ml |
|                      | final pH : 6.4-6.8 (at 20-25°C) |                 |

Although API C Medium contains gelling agent, it **requires no prior heating** and may be as easily pipetted as a liquid medium. It is preferable to warm it at room temperature a few hours before use. **Do not shake.**

**REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED****Reagents / Instrumentation**

- API Suspension Medium, 2 ml (Ref. 70 700) or API NaCl 0.85 % Medium, 2 ml (Ref. 20 070)
- Sabouraud Medium (Ref. 42 026 or 43 171 or equivalent)
- McFarland Standard (Ref. 70 900) No. 2
- API 20 C AUX Analytical Profile Index (Ref. 20 290), **apiweb™** identification software (Ref. 40 011), ATB™ instrument or **mini API®** (consult bioMérieux)
- RAT Medium [Rice Agar Tween]

**Material**

- Pipettes or PSipettes
- Ampule protector
- Ampule rack
- General microbiology laboratory equipment

**WARNINGS AND PRECAUTIONS**

- **For *in vitro* diagnostic use and microbiological control.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, yeast cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling yeasts should be observed throughout this procedure. Refer to "CLSI® M29-A, *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Current revision*". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories- CDC/NIH - Latest edition", or to the regulations currently in use in each country.



## QUALITY CONTROL

The strips and media are systematically quality controlled at various stages of their manufacture.

**Streamlined quality control** may be used to confirm acceptable performance of the API 20 C AUX system after shipping/storage. This methodology may be performed by following the instructions above for testing and meeting the criteria stated in CLSI® M50-A Quality Control for Commercial Microbial Identification Systems.

As there are no substrates that are consistently sensitive to degradation during shipping conditions, streamlined quality control may be conducted by testing two strains: ***Cryptococcus laurentii* ATCC® 18803** that is mostly positive and ***Candida glabrata* ATCC 15126**, which is mostly negative for reactions on the API 20 C AUX system.

For those users who are required to perform **comprehensive quality control** testing with the strip, the following three strains should be tested to demonstrate positive and negative reactivity for the most of the API 20 C AUX tests.

- |                                  |            |                                  |           |
|----------------------------------|------------|----------------------------------|-----------|
| 1. <i>Cryptococcus laurentii</i> | ATCC 18803 | 3. <i>Candida guilliermondii</i> | ATCC 6260 |
| 2. <i>Candida glabrata</i>       | ATCC 15126 |                                  |           |

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

|    | 0 | GLU | GLY | 2KG | ARA | XYL | ADO | XLT | GAL | INO | SOR | MDG | NAG | CEL | LAC | MAL | SAC | TRE | MLZ | RAF |
|----|---|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 1. | - | +   | -   | +   | +   | +   | +   | +   | +   | +   | +   | +   | +   | +   | +   | +   | +   | -   | +   | +   |
| 2. | - | +   | V   | -   | -   | -   | -   | -   | -   | -   | -   | -   | -   | -   | -   | -   | -   | +   | -   | -   |
| 3. | - | +   | +   | +   | +   | +   | +   | +   | +   | -   | +   | +   | +   | +   | -   | +   | +   | -   | +   | +   |

Profiles obtained after 48 hours of incubation after culture on Sabouraud agar.

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

### LIMITATIONS OF THE METHOD

- The API 20 C AUX system is intended uniquely for the identification of yeasts included in the database (see Identification Table at the end of this package insert). It cannot be used to identify any other microorganisms or to exclude their presence.
- Only pure cultures of a single organism should be used.

### RANGE OF EXPECTED RESULTS

Consult the Identification Table at the end of this package insert for the range of expected results for the various biochemical reactions.

### PERFORMANCES

5156 collection strains and strains of various origins belonging to species included in the database were tested :

- 89.7 % of the strains were correctly identified (with or without supplementary tests).
- 6.1 % of the strains were not identified.
- 4.2 % of the strains were misidentified.

### WASTE DISPOSAL

Unused ampules of API C Medium may be considered as non hazardous waste and disposed of accordingly. Dispose of all used or unused reagents (other than ampules of API C Medium) as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

### WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

|                       |        |
|-----------------------|--------|
| PROCEDURE             | p. I   |
| IDENTIFICATION TABLE  | p. II  |
| LITERATURE REFERENCES | p. III |
| INDEX OF SYMBOLS      | p. IV  |

BIOMERIEUX, the blue logo, API, ATB and **apiweb** are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

ATCC is a trademark belonging to American Type Culture Collection.

Any other name or trademark is the property of its respective owner.



**bioMérieux SA**  
 RCS LYON 673 620 399  
 69280 Marcy-l'Etoile / France  
 Tél. 33 (0)4 78 87 20 00  
 Fax 33 (0)4 78 87 20 90  
[www.biomerieux.com](http://www.biomerieux.com)

**bioMérieux, Inc**  
 Box 15969,  
 Durham, NC 27704-0969 / USA  
 Tél. (1) 919 620 20 00  
 Fax (1) 919 620 22 11  
 Printed in France

