

## Cryo-Babies® and CryoTags® Low Temperature Labels

For secure uncompromised labeling

- Temperature resistance -196°C to 150°C. Liquid and vapor phase LN<sub>2</sub>, water baths to 100°C and autoclaving and moderate ovens to 150°C
- Chemically inert resistant to most solvents and caustic agents. Adhere to most plastics (100% virgin PP tubes recommended), glass and metals without cracking, peeling, or degrading
- Label sheets will not jam printers. Flexible polyolefin material and acrylic adhesive. Unique combo expands and contracts together with temperature fluctuations



**A** autoclavable

**D** disposable

### Cryo-Babies® and CryoTags® on a Roll of 1000

code	color	label size	ideal use
HEACRY-0500	White	0.94" x 0.5"/24 x 13 mm	0.5 mL Tubes
HEACRY-1000	White	1.28" x 0.5"/33 x 13 mm	1.5 mL Tubes
HEACRY-1000(x)	<b>Specify</b>	1.28" x 0.5"/33 x 13 mm	1.5 mL Tubes
HEACRY-2000	White	1.50" x 0.75"/38 x 19 mm	General
HEACRY-2000(x)	<b>Specify</b>	1.50" x 0.75"/38 x 19 mm	General

Cryo-Babies®, CryoTags®, Tough-Spots®, Tough-Tags®, PetriSEAL™, ContainerSEAL™ and DuraSEAL™ are registered trademarks of Diversified Biotech.

### Laser Cryo-Babies® and CryoTags® on sheets (20 sheets per Pack)

code	color	label size	sheet layout	labels/ sheet	ideal use	no. labels
HEALCRY-2380(x)	<b>Specify</b>	0.94" x 0.5"/24 x 13 mm	7 Across/17 Down	119	0.5 mL Tubes	2,380
HEARNBW-2300	Rainbow*	0.94" x 0.5"/24 x 13 mm	7 Across/17 Down	119	0.5 mL Tubes	2,380
HEALCRY-1700(x)	<b>Specify</b>	1.28" x 0.5"/33 x 13 mm	5 Across/17 Down	85	1.5 mL Tubes	1,700
HEARNBW-2200	Rainbow*	1.28" x 0.5"/33 x 13 mm	5 Across/17 Down	85	1.5 mL Tubes	1,700
HEAMINI-1700	White	1.28" x 0.5"/33 x 13 mm	3 Across/17 Down	51	1.5 mL Tubes	1,020
HEALCRY-1200(x)	<b>Specify</b>	1.50" x 0.75"/38 x 19 mm	5 Across/12 Down	60	General	1,200
HEARNBW-2100	Rainbow*	1.50" x 0.75"/38 x 19 mm	5 Across/12 Down	60	General	1,200
HEALCRY-1100(x)	<b>Specify</b>	1.69" x 0.75"/43 x 19 mm	4 Across/13 Down	52	Cryogenic Vials	1,040
HEARNBW-2000	Rainbow*	1.69" x 0.75"/43 x 19 mm	4 Across/13 Down	52	Cryogenic Vials	1,040
HEASIDE-1000(x)	<b>Specify</b>	1.50" x 0.25"/38 x 6 mm	4 Across/39 Down	156	Microplates	3,120
HEARNBW-2400	Rainbow*	1.50" x 0.25"/38 x 6 mm	4 Across/39 Down	156	Microplates	3,120
HEALCRY-1258(x)	<b>Specify</b>	2.625" x 1.0"/67 x 25 mm	3 Across/10 Down	30	Racks & Boxes	600
HEARNBW-2500	Rainbow*	2.625" x 1.0"/67 x 25 mm	3 Across/10 Down	30	Racks & Boxes	600

To Specify color, substitute (x) with a letter - B=Blue, G=Green, O=Orange, R=Red, Y=Yellow

\*Rainbow packs contain 4 sheets of each Blue, Green, Orange, Red, Yellow