

Tough-Spots® Tube Labels

Pre-cut peel off labels for tube tops

- Resistant to most solvents and caustic agents. Adhere to most plastics, glass, and metals without cracking, peeling or degrading
- Label sheets will not jam printers. Flexible polyolefin material and acrylic adhesive

A autoclavable

D disposable

Tough-Spots® on a roll

part no.	color	label size	ideal use	labels/roll
HEATS-TNY	White	3/16" / 4.8 mm	0.2 mL Tubes	1500
HEAT-SPOTS	White	3/8" / 9.5 mm	0.5 mL Tubes	1000
HEAT-SPOTS-(x)	Specify	3/8" / 9.5 mm	0.5 mL Tubes	1000
HEATSPOTS-50	White	1/2" / 13 mm	1.5/2.0 mL Tubes	1000
HEATSPOTS-50(x)	Specify	1/2" / 13 mm	1.5/2.0 mL Tubes	1000
HEAT-SPOTS-44	White	7/16" / 11 mm	1.5/2.0 mL Tubes	1000
HEAT-SPOTS-75	White	3/4" / 19 mm	15-50 mL Tubes	500
HEAT-SPOTS-100	White	1.0" / 25 mm	50 mL Tubes	500

To **Specify** color, substitute (x) with a letter: B=Blue, C=Coffee, D=Gold, G=Green, L=Lavender, M=Maroon, O=Orange, P=Pink, R=Red, S=Silver, Y=Yellow

Resistance
-196°C to
+150°C



Laser Tough-Spots® On sheets. 20 sheets per pack

part no.	color	label size	sheet layout	labels/sheet	ideal use	no. labels
HEASPOT-1(x)00	Specify	3/8" / 9.5 mm	12 Across/16 Down	192	0.5 mL Tubes	3,840
HEARNBW-1000	Rainbow*	3/8" / 9.5 mm	12 Across/16 Down	192	0.5 mL Tubes	3,840
HEASPOT-2(x)00	Specify	1/2" / 13 mm	12 Across/16 Down	192	1.5-2.0 mL Tubes	3,840
HEARNBW-1100	Rainbow*	1/2" / 13 mm	12 Across/16 Down	192	1.5-2.0 mL Tubes	3,840
HEASPOT-5000	White	7/16" / 11 mm	13 Across/18 Down	234	1.5-2.0 mL Tubes	4,680
HEASPOT-3000	White	3/4" / 19 mm	9 Across/12 Down	108	15-50 mL Tubes	2,160
HEASPOT-4000	White	1.0" / 25 mm	7 Across/8 Down	63	50 mL Tubes	1,260

To **Specify** color, substitute (x) with a number: 0=White, 1=Blue, 2=Green, 3=Orange, 4=Red, 5=Yellow; *Rainbow packs contain 4 sheets of each.

Tough-Tags® On a Roll of 1000

part no.	color	label size	ideal use
HEATT-TNY	White	0.81" x 0.28" / 21 x 7 mm	0.2 mL Tubes
HEATTSW-1000	White	0.94" x 0.5" / 24 x 13 mm	0.5 mL Tubes
HEATTS(x)-1000	Specify	0.94" x 0.5" / 24 x 13 mm	0.5 mL Tubes
HEATTLW-1000	White	1.28" x 0.5" / 33 x 13 mm	1.5/2.0 mL Tubes
HEATTL(x)-1000	Specify	1.28" x 0.5" / 33 x 13 mm	1.5/2.0 mL Tubes
HEATTGP-1000	White	1.50" x 0.75" / 38 x 19 mm	General
HEATT-SWALL	White	1.50" x 0.25" / 38 x 6 mm	Microplates

To **Specify** color, substitute (x) with a letter: B=Blue, G=Green, O=Orange, R=Red, Y=Yellow



Laser Tough-Tags® On sheets. 25 sheets per pack

part no.	color	label size	sheet layout	labels/sht	ideal use	no. labels
HEATTSW-2240	White	0.94" x 0.5" / 24 x 13 mm	7 Across/17 Down	119	0.5 mL Tubes	2,975
HEATTSW-2240(x)	Specify	0.94" x 0.5" / 24 x 13 mm	7 Across/17 Down	119	0.5 mL Tubes	2,975
HEARNBW-3100	Rainbow*	0.94" x 0.5" / 24 x 13 mm	7 Across/17 Down	119	0.5 mL Tubes	2,975
HEATTLW-2016	White	1.28" x 0.5" / 33 x 13 mm	5 Across/17 Down	85	1.5/2.0 mL Tubes	2,125
HEATTLW-2016(x)	Specify	1.28" x 0.5" / 33 x 13 mm	5 Across/17 Down	85	1.5/2.0 mL Tubes	2,125
HEARNBW-3000	Rainbow*	1.28" x 0.5" / 33 x 13 mm	5 Across/17 Down	85	1.5/2.0 mL Tubes	2,125
HEATTGP-1050	White	1.50" x 0.75" / 38 x 19 mm	5 Across/12 Down	60	General	1,500
HEATTGP-1050(x)	Specify	1.50" x 0.75" / 38 x 19 mm	5 Across/12 Down	60	General	1,500
HEARNBW-3200	Rainbow*	1.50" x 0.75" / 38 x 19 mm	5 Across/12 Down	60	General	1,500
HEAMISL-1000	White	0.875" x 0.875" / 22 x 22 mm	8 Across/12 Down	96	Microscope Slides	2,400
HEAMISL-1000(x)	Specify	0.875" x 0.875" / 22 x 22 mm	8 Across/12 Down	96	Microscope Slides	2,400
HEARNBW-3300	Rainbow*	0.875" x 0.875" / 22 x 22 mm	8 Across/12 Down	96	Microscope Slides	2,400

To **Specify** color, substitute (x) with a letter: B=Blue, G=Green, O=Orange, R=Red, Y=Yellow; *Rainbow packs contain 5 sheets of each.