

smartBOATS™

A truly green solution to sample weighing

SmartBoats® are made from FSC® Certified (FSC® Mix) paper (Eco-IN) and when disposed of in the standard waste stream they will naturally biodegrade (Eco-OUT)

- Control contents like a polystyrene weigh boat
- An alternative to regular weighing papers as they are much easier to handle
- Made of grease proof paper that meets FDA 21 CFR 176.170 (components of paper in contact with aqueous and fatty foods) and with 21 CFR 176.180 (contact with dry foods)
- Compatible with heat applications with a temperature limitation of 210°C
- Available in a small and medium in both brown for light colored samples and white for dark samples

Certifications:

USA: ISO 9002 certified facility in accordance with BGV, Haccp and FDA regulations. FSC® Certification

EU: ISO 14001 Certification - Environmental policy, BRC/IOP Certification, ISO 9001 Certification, FSC® Certification



| Part No. | Size | L x W x H | Pkg. | Color |
|-----------------|--------|-----------------------|--------|---------|
| HS120312 | Small | 35 x 35 x 15.8 mm | 500/pk | Brown ■ |
| HS120313 | Medium | 47.6 x 47.6 x 14.2 mm | 500/pk | Brown ■ |
| HS120314 | Small | 35 x 35 x 15.8 mm | 500/pk | White □ |
| HS120315 | Medium | 47.6 x 47.6 x 14.2 mm | 500/pk | White □ |



Weighing Boats

Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss.


- Select white, black, or blue for best color contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Can be used with microwave moisture analyzers
- All are FDA/food-grade compliant
- Disposable



Standard Weighing Boats

- Boats are shallow with a wide, flat bottom to resist tipping
- Rounded corners simplify transfer

| Regular Material Part No. | Antistatic Material Part No. | Size | L xW x H | Pkg. | Color |
|---------------------------|------------------------------|--------|-------------------|--------|---|
| HS1420A | HS1420AA | Small | 46 x 46 x 8 mm | 500/pk | White <input type="checkbox"/> |
| HS1420B | HS1420BB | Medium | 85 x 85 x 24 mm | 500/pk | White <input type="checkbox"/> |
| HS1420C | HS1420CC | Large | 140 x 140 x 22 mm | 500/pk | White <input type="checkbox"/> |
| HS1421A | HS120222 | Small | 46 x 46 x 8 mm | 500/pk | Blue <input checked="" type="checkbox"/> |
| HS1421B | HS120223 | Medium | 85 x 85 x 24 mm | 500/pk | Blue <input checked="" type="checkbox"/> |
| HS1421C | HS120224 | Large | 140 x 140 x 22 mm | 500/pk | Blue <input checked="" type="checkbox"/> |
| HS1423A | HS1423AA | Small | 46 x 46 x 8 mm | 500/pk | Black <input checked="" type="checkbox"/> |
| HS1423B | HS1423BB | Medium | 85 x 85 x 24 mm | 500/pk | Black <input checked="" type="checkbox"/> |
| HS1423C | HS1423CC | Large | 140 x 140 x 22 mm | 500/pk | Black <input checked="" type="checkbox"/> |



Hexagonal Weighing Boats

- Bend flexible boats at points to form pour spouts
- Flat bottoms and sloping sides provide accurate transfer


| Regular Material Part No. | Antistatic Material Part No. | Size | L xW x H | Pkg. | Color |
|---------------------------|------------------------------|---------|------------------|--------|--------------------------------|
| — | HS14251A | Small | 38 x 25 x 9 mm | 500/pk | White <input type="checkbox"/> |
| — | HS14251B | Medium | 64 x 51 x 19 mm | 500/pk | White <input type="checkbox"/> |
| — | HS14251C | Large | 102 x 76 x 22 mm | 500/pk | White <input type="checkbox"/> |
| — | HS14251D | X-large | 127 x 89 x 31 mm | 500/pk | White <input type="checkbox"/> |



Pour-Boat Weighing Dishes

- Pour-spout design easily dispenses samples without flexing the dish
- Flat-bottom dishes are made from 12-mil thick polystyrene for better strength, stability, and handling

| Regular Material Part No. | Antistatic Material Part No. | Size | L xW x H | Pkg. | Color |
|---------------------------|------------------------------|--------|-------------------|--------|--------------------------------|
| HS1419A | HS120225 | Small | 43 x 58 x 13 mm | 500/pk | White <input type="checkbox"/> |
| HS1419B | HS120226 | Medium | 83 x 132 x 26 mm | 500/pk | White <input type="checkbox"/> |
| HS1419C | HS120227 | Large | 108 x 183 x 26 mm | 500/pk | White <input type="checkbox"/> |



Square Weighing Boats

- Shallow boat with a wide, flat bottom to resist tipping
- Ideal for solid product weighing

| Regular Material Part No. | Antistatic Material Part No. | Size | L xW x H | Pkg. | Color |
|---------------------------|------------------------------|--------|----------------|--------|--------------------------------|
| HS1422 | — | Medium | 90 x 90 x 7 mm | 500/pk | White <input type="checkbox"/> |

Diamond Weighing Boats

- Diamond shape that easily flexes for accurate pouring
- Wide, flat bottom to resist tipping



| Regular Material Part No. | Antistatic Material Part No. | Size | L xW x H | Pkg. | Color |
|---------------------------|------------------------------|--------|------------------|--------|---|
| HS1424AA | HS1425AA | 5 mL | 55 x 35 x 6 mm | 500/pk | White <input type="checkbox"/> |
| HS1424A | HS1425A | 25 mL | 71 x 46 x 14 mm | 500/pk | White <input type="checkbox"/> |
| HS1424B | HS1425B | 30 mL | 78 x 56 x 14 mm | 500/pk | White <input type="checkbox"/> |
| HS1424C | HS1425C | 100 mL | 119 x 90 x 19 mm | 500/pk | White <input type="checkbox"/> |
| HS1427AA | HS1426AA | 5 mL | 55 x 35 x 6 mm | 500/pk | Black <input checked="" type="checkbox"/> |
| HS1427A | HS1426A | 25 mL | 71 x 46 x 14 mm | 500/pk | Black <input checked="" type="checkbox"/> |
| HS1427B | HS1426B | 30 mL | 78 x 56 x 14 mm | 500/pk | Black <input checked="" type="checkbox"/> |
| HS1427C | HS1426C | 100 mL | 119 x 90 x 19 mm | 500/pk | Black <input checked="" type="checkbox"/> |

Flexible boats bend to form pour spout