

# A AUTOCLAVES Vertical, Economy

## UTFB/UTKBS, Full Stainless steel Hand Wheel Type of Quick-open Door Structure Pressure Safety Lock



UTKBS-50LV

Computer control auto recycle sterilization, LCD screen indicate working state, Temp. scope 50-126°C, Time scope 0-99h, Overpressure auto-discharging 0.145-0.165MPa, Automatic shut down with beep reminding after sterilization, Hand wheel type of quick-open door structure, Pressure safety lock.



UTKBS-50LV

**ECONOMY**



UTFB-30LV

### Vertical Pressure Autoclave

- Computer control auto recycle sterilization.
- LCD Screen indicate working state.
- Temp. scope 50-134°C.
- Time scope 0-99h.
- Overpressure auto-discharging 0.145-0.165Mpa.
- Automatic shut down with beep reminding after sterilization.

**Optional: Printer.**



UTFB-30LV

Model	UTFB-30LV	UTKBS-50LV	UTKBS-75LV	UTKBS-100LV	UTKBS-150LV
Capacity	30L	50L	75L	100L	150L
Power source	220V	220V	220/380V	380V	380V
Consumed power	3kw	3.5kw	3.5/4.5kw	4.5kw	6kw
Chamber Site	Ø350x330	Ø350x520	Ø388x650	Ø500x550	Ø500x760
Weight	60kg	80kg	90kg	100kg	120kg
Outside Size (mm)	580x580x1000	600x580x1250	600x600x1400	630x630x1300	630x630x1450