

gloveon VinBlue

Vinyl Exam Gloves Powder Free, Standard Cuff

VinBlue is a light-weight, powder-free vinyl glove, with a beaded cuff for comfortable fit. Made from 100% PVC, VinBlue is light and flexible, reducing hand and muscle fatigue. Food and chemical tested, VinBlue is also free of latex, eliminating the risk and discomfort associated with latex allergies. The blue colour makes VinBlue instantly recognizable and familiar to food process workers.



		GloveOn VinBlue	
Length (mm)		≥ 240	
Thickness Measurements (mm)			
Palm (centre of Palm)		0.08 ± 0.02	
Finger (13mm ± 3mm from tip)		0.10 ± 0.02	
Physical Properties		Before Ageing	After Ageing
Tensile Strength (MPa)		≥ 16	≥ 14
Elongation (%)		≥ 350	≥ 350
Inspection Levels & AQL		Inspection Level	AQL
Watertightness		G1	1.50
Physical Dimensions		S2	4.00
Tensile Strength		S2	4.00
Visual Inspection (Major)		S4	4.40
Visual Inspection (Minor)		S4	4.00
Particulate Residue		N = 5	≤ 2mg/glove

REORDER CODE

VNY21XS XSMALL
VNY21SS SMALL
VNY21MM MEDIUM
VNY21LL LARGE
VNY21XL X-LARGE

FEATURES

- Beaded Cuff
- Powder free
- Not made with natural rubber latex
- Food contact tested
- Lab chemical tested
- Ambidextrous
- Blue colour

PACKAGING

200 gloves per box for XS to XL
10 boxes per carton

REGULATORY COMPLIANCE

TGA - ARTG 100687,

STANDARDS

EN 13130-1:2004
EN 1186-1:2002
AS/NZs 4011-1997

MANUFACTURING ACCREDITATIONS

ISO 9001:2008
ISO 13485:2003
EN ISO 13485:2003

A brand by

mun
www.munglobal.com

Mun Australia
Toll free: 1800 456 837 (GLOVES)
info.au@munglobal.com

 **Birthing Kit Foundation**
Australia
Proud Supporter of
Birthing Kit Foundation (Australia)

 Protection Always On