

gloveon Matador



Nitrile Exam Gloves Powder Free, Long Cuff

With quality, consistency and reliability built into the GloveOn's 8 mil nitrile glove, even a bull fighter will have the courage, trust and confidence to go into the ring to face the challenges while donning the **Matador**.



GloveOn Matador		
Length (mm)	≥ 280	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.20 ± 0.02	
Finger (13mm ± 3mm from tip)	0.24 ± 0.02	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 400
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.50
Physical Dimensions	S2	4.00
Tensile Strength	S2	4.00
Visual Inspection (Major)	S4	2.50
Visual Inspection (Minor)	S4	4.00
Particulate Residue	N = 5	≤ 2mg/glove

Tested in accordance with EN 374-3:2003 with the following chemical permeation performance levels:

Chemical	Performance Level	Permeation Performance Level	Measured Breakthrough Time (Minutes)
4% Chlorhexidine gluconate	6		
40% Sodium hydroxide	6	1	> 10
35% Gluteraldehyde	6	2	> 30
1% Gluteraldehyde	6	3	> 60
4 % Gluteraldehyde	6	4	> 120
10% Formalin (Formaldehyde)	6	5	> 240
35% Ethanol	3	6	> 480
n-Hexane	2		
Isopropanol	3		
Acetone	Not recommended		
3% Hydrogen Peroxide	6		
n-Heptane (J)	5		
96% Sulphuric Acid (K)	2		
40% Sodium Hydroxide (L)	6		

REORDER CODE

NTR18XS	X-SMALL
NTR18SS	SMALL
NTR18MM	MEDIUM
NTR18LL	LARGE
NTR18XL	X-LARGE

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Long cuff
- Sky blue colour

PACKAGING

50 gloves per box for XS to XL
10 boxes per carton

REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k),
MDD 93/42/EEC, REACH,
ROHS Directive 2002/95/EC,
PPE 89/686/EEC

STANDARDS

ASTM D6319, EN 455 part 1, 2, 3 & 4,
EN 374 part 1, 2, 3, AS/NZS 4011:1997

MANUFACTURING ACCREDITATIONS

ISO 9001:2008
ISO 13485:2003
EN ISO 13485:2003

A brand by

mun
www.munglobal.com

Mun Australia
Toll free: 1800 456 837 (GLOVES)
info.au@munglobal.com

**Birthing Kit
Foundation**
Australia
Proud Supporter of
Birthing Kit Foundation (Australia)

gn Protection Always On