

# gloveon Hammer

## Nitrile Exam Gloves Powder Free, Standard Cuff



**Hammer** is the nitrile glove meant to deliver an impact for the tough jobs that demand high tensile strength, elongation and a definitive colour. This glove is suitable for security, tattooist, beauty salon, motor mechanic, cleaner and general maintenance.



GloveOn Hammer		
<b>Length (mm)</b>	≥ 240	
<b>Thickness Measurements (mm)</b>		
Palm (centre of Palm)	0.08 ± 0.02	
Finger (13mm ± 3mm from tip)	0.12 ± 0.02	
<b>Physical Properties</b>	<b>Before Ageing</b>	<b>After Ageing</b>
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 500
Elongation after Isopropanol immersion (%)	≥ 500	N/A
<b>Inspection Levels &amp; AQL</b>	<b>Inspection Level</b>	<b>AQL</b>
Watertightness	G1	1.50
Physical Dimensions	S2	4.00
Tensile Strength	S2	4.00
Visual Inspection (Major)	S4	2.50
Visual Inspection (Minor)	S4	4.00
Particulate Residue	N = 5	≤ 2mg/glove

### REORDER CODE

NTR23SS	SMALL
NTR23MM	MEDIUM
NTR23LL	LARGE
NTR23XL	X-LARGE

### FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- UL certified to meet NFPA 1999
- Ambidextrous
- Standard cuff
- Black colour

### PACKAGING

100 gloves per box for S to L  
90 gloves per box for XL  
10 boxes per carton

### REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k),  
MDD 93/42/EEC, REACH,  
ROHS Directive 2002/95/EC

### STANDARDS

ASTM D6319, ASTM D6978,  
EN455 part 1, 2, 3 & 4, NFPA 1999

### MANUFACTURING ACCREDITATIONS

ISO 9001:2008  
ISO 13485:2003  
EN ISO 13485:2012

A brand by

**mun**  
www.munglobal.com

**Mun Australia**  
Toll free: 1800 456 837 (GLOVES)  
info.au@munglobal.com

 **Birthing Kit  
Foundation  
Australia**  
Proud Supporter of  
Birthing Kit Foundation (Australia)

 Protection Always On