

# gloveon Maverick

## Nitrile Exam Gloves Powder Free, Standard Cuff



**Maverick** the All-Rounder nitrile glove has proven itself in medical, laboratory and industrial use. Tested for tensile strength and chemical resistance, Maverick's elastic formulation mimics the donning properties of latex gloves in softness, stretchability, tactility and comfort fitting.



GloveOn Maverick		
<b>Length (mm)</b>	≥ 230	
<b>Thickness Measurements (mm)</b>		
Palm (centre of Palm)	≥ 0.08	
Finger (13mm ± 3mm from tip)	≥ 0.11	
<b>Physical Properties</b>	<b>Before Ageing</b>	<b>After Ageing</b>
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 400
<b>Inspection Levels &amp; AQL</b>	<b>Inspection Level</b>	<b>AQL</b>
Watertightness	G1	1.50
Physical Dimensions	S2	4.00
Tensile Strength	S2	4.00
Visual Inspection (Major)	S4	2.50
Visual Inspection (Minor)	S4	4.00
Particulate Residue	N = 5	≤ 2mg/glove

### REORDER CODE

NTR28XS	X-SMALL
NTR28SS	SMALL
NTR28MM	MEDIUM
NTR28LL	LARGE
NTR28XL	X-LARGE

### FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Standard cuff
- Aqua blue colour

### PACKAGING

200 gloves per box for XS to L  
180 gloves per box for XL  
10 boxes per carton

### REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k), MDD 93/42/EEC, REACH, ROHS Directive 2002/95/EC, EC 10/2011, REGULATION (EC) No 1935/2004, EC 1935/2004, PPE 89/686/EEC

### STANDARDS

ASTM D6319, ASTM D6978, EN 455 part 1, 2, 3 & 4, EN 1186, EN 13130, CEN/TS 14234, EN 420, EN 374 part 1, 2 & 3, AS/NZS 4011:1997

### MANUFACTURING ACCREDITATIONS

ISO 9001:2008  
ISO 13485:2003  
EN ISO 13485:2003

### Chemotherapy Drugs and Concentration (Tested for Resistance to Permeation by Chemotherapy Drugs as per ASTM D6978-05-Test Report PN 104145B)

Minimum Breakthrough  
Detection Time (minutes)

Carmustine (BCNU), 3.3mg/ml (3,300 ppm)	15.3 minutes
Cisplatin, 1.0mg/ml (1,000 ppm)	>240 minutes
Cyclophosphamide (Cytoxan), 20.0mg/ml (20,000 ppm)	>240 minutes
Dacarbazine (DTIC), 10.0mg/ml (10,000 ppm)	>240 minutes
Doxorubicin Hydrochloride, 2.0mg/ml (2,000 ppm)	>240 minutes
Etoposide (Toposar), 20.0mg/ml (20,000 ppm)	>240 minutes
Fluorouracil, 50.0mg/ml (50,000 ppm)	>240 minutes
Methotrexate, 25.0mg/ml (25,000 ppm)	>240 minutes
Mitomycin C, 0.5mg/ml (500 ppm)	>240 minutes
Paclitaxel (Taxol), 6.0mg/ml (6,000 ppm)	>240 minutes
Thiotepa, 10.0mg/ml (10,000 ppm)	46.3 minutes
Vincristine Sulfate, 1.0mg/ml (1,000 ppm)	>240 minutes

**WARNING:** Carmustine and Thiotepa, at the tested concentration, degraded Maverick nitrile glove at 15.3 minutes and 46.3 minutes, respectively. The safe use of gloves in chemotherapy treatment is solely the decision of clinicians authorised to make such decision.

A brand by

**dmun**  
www.munglobal.com

**Mun Australia**  
Toll free: 1800 456 837 (GLOVES)  
info.au@munglobal.com

 **Birthing Kit  
Foundation  
Australia**  
Proud Supporter of  
Birthing Kit Foundation (Australia)

 Protection Always On