



gloveon

Aegis

Sterile Nitrile Exam Gloves Powder Free, Long Cuff

Aegis, like an armour, offers superb hand barrier protection and it dons like a second skin. The extended cuff length provides extra protection, high tensile strength and resistance to chemicals as expected from a high performance nitrile glove. This is our recommended sterile exam glove for procedures and jobs that demand tactile sensitivity.



GloveOn Aegis		
Length (mm)	≥ 280	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.08 ± 0.02	
Finger (13mm ± 3mm from tip)	0.12 ± 0.02	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 400
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.50
Physical Dimensions	S2	4.00
Tensile Strength	S2	4.00
Visual Inspection (Major)	S4	2.50
Visual Inspection (Minor)	S4	4.00
Particulate Residue	N = 5	≤ 2mg/glove

REORDER CODE

SNT56SS	SMALL	(6.0 - 6.5)
SNT56MM	MEDIUM	(7.0 - 7.5)
SNT56LL	LARGE	(8.0 - 8.5)
SNT56XL	X-LARGE	(9.0)

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Long cuff
- Violet blue colour
- Sterilised using gamma irradiation

PACKAGING

1 pair per pouch
50 pairs per box
4 boxes per carton

REGULATORY COMPLIANCE

TGA - ARTG 203002, FDA 510(k),
MDD 93/42/EEC, REACH,
ROHS Directive 2002/95/EC

STANDARDS

ASTM D6319, EN455 part 1,2,3 & 4,
ISO 11193-1, ISO 11137

MANUFACTURING ACCREDITATIONS

ISO 9001:2008
ISO 13485:2003
EN ISO 13485:2012

Chemotherapy Drugs and Concentration (Tested for Resistance to Permeation by Chemotherapy Drugs as per ASTM D6978-05-Test Report PN 104145A)

Minimum Breakthrough
Detection Time (minutes)

Carmustine (BCNU), 3.3mg/ml (3,300 ppm)	1.2 minutes
Cisplatin, 1.0mg/ml (1,000 ppm)	>240 minutes
Cyclophosphamide (Cytoxan), 20.0mg/ml (20,000 ppm)	>240 minutes
Dacarbazine (DTIC), 10.0mg/ml (10,000 ppm)	>240 minutes
Doxorubicin Hydrochloride, 2.0mg/ml (2,000 ppm)	>240 minutes
Etoposide (Toposar), 20.0mg/ml (20,000 ppm)	>240 minutes
Fluorouracil, 50.0mg/ml (50,000 ppm)	>240 minutes
Methotrexate, 25.0mg/ml (25,000 ppm)	>240 minutes
Mitomycin C, 0.5mg/ml (500 ppm)	>240 minutes
Paclitaxel (Taxol), 6.0mg/ml (6,000 ppm)	>240 minutes
Thiotepa, 10.0mg/ml (10,000 ppm)	30.9 minutes
Vincristine Sulfate, 1.0mg/ml (1,000 ppm)	>240 minutes

WARNING: Carmustine and Thiotepa, at the tested concentration, degraded Aegis sterile nitrile glove at 1.2 minutes and 30.9 minutes, respectively. The safe use of gloves in chemotherapy treatment is solely the decision of clinicians authorised to make such decision.

A brand by

D:mun
www.munglobal.com

Mun Australia
Toll free: 1800 456 837 (GLOVES)
info.au@munglobal.com

**Birthing Kit
Foundation**
Australia
Proud Supporter of
Birthing Kit Foundation (Australia)

gn Protection Always On