

PXi range—chemiluminescence/fluorescence gel and blot imaging



PXi versions	System pixels	Lens	Max field of view
PXi 4EZ	4.0 mp	Motor driven zoom f1.2	18 x 13.5 cm
PXi 4 (Touch)	4.0 mp	Fixed f0.95	15 x 11 cm
PXi 6 (Touch)	6.0 mp	Fixed f1.4	13 x 10 cm
PXi 9 (Touch)	9.0 mp	Fixed f1.4	13 x 10 cm
PXi 6 Access (Touch)	6.0 mp	Fixed f0.95	21 x 17 cm
PXi 9 Access (Touch)	9.0 mp	Fixed f0.95	21 x 17 cm