

NEW



pH Calibration Solutions

230 mL Value Series

 **HANNA**[®]
instruments



Designed for value.

- 230 mL size
- New packaging with clear labeling of pH value
- New twist cap and foil inner seal with pull tab
- Certificate of analysis included

Hanna has designed three new pH calibration solutions in redesigned 230 mL bottles. The new twist off cap and inner foil seal with pull tab allows for easier opening. Each bottle label is clearly marked with their respective value and includes a pH/temperature reference chart plus a designated spot to mark when the bottle was opened. Each solution is supplied with a certificate of analysis.



HANNA
instruments

Certificate of Analysis

Product name: Buffer Solution, pH 10.01 ± 0.02 @ 25 °C
Product code: HI7010-023
Lot number: 0285
Valid until: June 2018
Date of calibration: 2016-06-13
pH: 10.01 ± 0.02 @ 25 °C

DATE OPENED

°C	°F	pH
0	32	7.13
5	41	7.10
10	50	7.07
15	59	7.05
20	68	7.03
25	77	7.01
30	86	7.00
35	95	6.99
40	104	6.98
45	113	6.98
50	122	6.98
55	131	6.98
60	140	6.99
65	149	6.99
	158	7.00



Value Series

All new solutions and packaging in 230 mL bottles.



HI7004-023

pH 4.01 Calibration Solution

±0.02 pH @ 25°C (77°F)

- includes certificate of analysis



HI7007-023

pH 7.01 Calibration Solution

±0.02 pH @ 25°C (77°F)

- includes certificate of analysis



HI7010-023

pH 10.01 Calibration Solution

±0.02 pH @ 25°C (77°F)

- includes certificate of analysis

HANNA
instruments

hannainst.com

HANNA
instruments