


# RAL STAINER

## Technical features

- Bath staining machine
- LCD touch screen
- Utilizes RAL Staining Kits
- RFID tags
- Autonomy from 1 to 20 slides (2 racks of 10 slides)
- 2 loading/unloading/drying zones
- 4 staining stations, 1 rinsing station
- Rinsing liquid reservoir
- Automatic waste removal in 10L container

- Height: 58.7cm – Diameter: 40cm
- Weight: 28 kg
- Supply voltage range: 110-230V/50-60Hz/150 VA

Reference : 414 195



## Kit RAL Stainer Cold ZN

- Suitable for 500 slides
- Life Span: 14 days after opening

Reference : 414 194



Fixative 1 x 200 ml  
Fuchsin 1 x 200 ml  
Discolouring Solution 1 x 200 ml  
Methylene Blue 1 x 200 ml

## Kit RAL Stainer Fluo-RAL

- Suitable for 300 slides
- Life Span: 7 days after opening

Reference : 414 193



Fixative Fluo-RAL 1 x 200 ml  
Auramine Fluo-RAL 1 x 100 ml  
Solution for Auramine 1 x 100 ml  
Discolouring Solution Fluo-RAL 1 x 200 ml  
Thiazine Red Fluo-RAL 1 x 100 ml  
Solution for Thiazine Red 1 x 100 ml

Distributed by

bioMérieux S.A.  
69280 Marcy l'Etoile  
France  
Tel.: 33 (0)4 78 87 20 00  
Fax: 33 (0)4 78 87 20 90  
[www.biomerieux.com](http://www.biomerieux.com)

Manufactured by



Site Montesquieu - 33650 Martillac  
France



# RAL STAINER

Standardize and secure your staining protocols in tuberculosis

**SIMPLE | RELIABLE | FLEXIBLE | SAFE**

