



POWERFUL PARTNERS FOR MICROBIAL DETECTION & STERILITY TESTING

Culture media range

*Choose a product and
a supplier you can rely on!*



BIOMÉRIEUX
INDUSTRY

PIONEERING DIAGNOSTICS


Sterility testing culture media range

With more than 50 years of design and manufacturing experience **bioMérieux understands culture media.**

Detecting contaminants is a critical step for the release of sterile pharmaceuticals – for this purpose bioMérieux offers complete peace of mind through **the sensitivity and performance of our sterility testing culture media range.**


Designed with the pharmaceutical industry in mind all our products conform to the harmonized requirements of the USP, EP and JP in terms of both formulation and QC Growth Performance Testing methods.



 By controlling the entire manufacturing process, all the way from the raw materials to the final product, bioMérieux can ensure our products give consistent performance from batch to batch – everytime!

 ISO 9001 compliant manufacturing facilities that work in accordance with GMP requirements.

 Broad variety of presentation formats, with septa, screw cap design and a range of tubes.

 Thio-clear medium formulated using the unique bioMérieux know-how for ultra-clarity and ease of reading.



Clean Room range

Sterile inside and out – Ethylene Oxide vapour treated bottles and dust-minimal plastic packaging ensure these media are suitable for use in clean rooms and isolators.

| Media | Reference | Format |
|-------------------------------------|-----------|---|
| Trypcase Soy Broth | 44011 | 12 x 100 ml bottles for membrane filtration |
| | 42633 | 6 x 100 ml bottles for membrane filtration (CLEAN ROOM) |
| | 42100 | 20 x 9 ml tubes |
| Fluid Thioglycollate Medium (FTM) | 44001 | 12 x 100 ml bottles for membrane filtration |
| | 42074 | 20 x 9 ml tubes |
| Thioglycollate (Thio-clear) | 44021 | 12 x 100 ml bottles for membrane filtration |
| | 42634 | 6 x 100 ml bottles for membrane filtration (CLEAN ROOM) |
| | 42018 | 20 x 9 ml tubes |
| Fluid A Rinsing solution | 42044 | 12 x 100 ml bottles for membrane filtration |
| | 42616 | 4 x 300 ml bottles for membrane filtration |
| | 42640 | 6 x 100 ml bottles for membrane filtration (CLEAN ROOM) |
| Fluid D Rinsing solution with tween | 42624 | 4 x 300 ml bottles for membrane filtration |
| | 42644 | 6 x 100 ml bottles for membrane filtration (CLEAN ROOM) |

bioMérieux provides comprehensive documentation, training and validation support.

bioMérieux S.A.
69280 Marcy l'Etoile
France
Tel.: +33 (0)4 78 87 20 00
Fax: +33 (0)4 78 87 20 90

www.biomerieux.com
www.biomerieux-industry.com

