





COMBOURG RANGE

RANGE	TYPE DE CONTENANT	CAPS	SIZES (Dia X h) mm	PACKAGING	OPENING DIAMETER mm	MINIMUM VOLUME validated per container	MAXIMUM VOLUME validated per container	MINIMUM STANDARD BATCH SIZE (l)	MAXIMUM STANDARD BATCH SIZE (l)	COMMENTS	PHOTOS
TUBES	Tube 17X125	Blue, green, red or gold	17X125	20 & 100 TUBES	11	2,5	15	26	150	Double thread	STANDARD TUBES
	Tube 17,75X180	Blue, green, red or gold	17,75X180	100 TUBES	11	9	20	50	120	Double thread	
	Tube 19,75X200	Blue, green, red or gold	19,75X200	100 TUBES	12,25	18	20	50	80	Round-bottomed ; simple f	
PACKAGING : printed tubes & labelled boxes											
SIROP BOTTLES	125ml	Blue, green, red or gold	49X114,4	6 bottles	20	20	104	100	260	/	STANDARD SIROP BOTTLES
	150ml	Blue, green, red or gold	51,6X122,8	6 bottles	20	90	100	100	220	/	
	250ml	Blue, green, red or gold	61,7X142,9	6 bottles	20	90	200	200	500	/	
	300ml	Blue, green, red or gold	65,6X151,3	6 bottles	20	200	286	160	490	/	
	500ml	Blue, green, red or gold	76,8X176,8	6 bottles	20	300	460	275	320	/	
PACKAGING : printed bottles & labelled boxes											
SEPTUM PERFORABLE BOTTLES	120ml	Red or gold	48,5X106,1	6 bottles	33	50	100	55	260	Can be simple, double or three wrapped	 <p>1&2=120 and 250 3= Bibby bottle (no listed in the sheet) 4&5=500 and 1000</p>
	250ml	Red or gold	66,2X130,5	6 bottles	30	170	250	80	280		
	500ml	Red or gold	90,5X152	6 bottles	30	250	500	140	490		
	1000ml	Red or gold	90,5X247	6 bottles	30	500	1000	95	560		
PACKAGING : printed bottles & labelled boxes											
SI-MAX BOTTLES	1000ml	/	101X225	6 bottles	29,8mm	490	990	120	560		
PACKAGING : printed bottles & labelled boxes											
Twist off jars	156ml	Gold	56,4X85,5	6 jars	42,5	90	102	40	140	/	
	212ml	Gold	67,4X84,5	12 jars	48,5	20	120	40	140	/	
	370ml	Gold	71,4X119	6 jars	58,9	225	225	60	210	/	
	660ml	Gold	86,8X136,8	12 jars	66,5	100	450	60	210	/	
	1700ml	Gold		8 jars	66,5	900	1000	195	550	/	
PACKAGING : printed jars & labelled boxes, except for 1700ml jars : labelled											
Twist off bottles	1000ml	Gold	88,3X254,5	6 bottles	43	600	1000	150	480	/	
PACKAGING : printed jars & labelled boxes											
BAGS	3000 ml	/	/	4 BAGS	/	3000	3000	550	1095	/	STANDARD BAGS
	5000 ml	/	/	2 BAGS	/	3,375	5000	693	1385	/	
PACKAGING : bags and boxes are labelled											
PLATES	90mm	/	DIA 90MM	10, 20 or 120	/	16	30	55	440	can be double or three wrapped	STANDARD PETRI DISHES PLATES
	140mm	/	DIA 140MM	10	/	60	70	40	70	/	
	55mm	/	DIA 55MM	10 or 120	/	9	9	50	110	/	
	Count Tact	/	DIA 65MM	10, 60 or 120	/	/	/	45	220	Can be double or three wrapped	
	Square box	/	120X120MM	10	/	50	70	50	70	/	
PACKAGING : Petri dishes are printed only on the bottom and boxes are labelled											
MINIBAGS	MINIBAG 225	/	/	10 MBGS	12	225	225	504	504		
PACKAGING : Mini-Bags are pre-printed w REF, NAME of the product... Study need to be done for others media /Annual quantity											
Bijous & universals bottles	BIJOUS BOTTLES	Grey	20,25X47,25	50 bottles	11,8	3	5	51	261	/	
	UNIVERSAL BOTTLES	Grey	27,75X83	50 bottles	18	9	15	261	541	/	
PACKAGING : Bijous and universals bottles and boxes are labelled											
General comments	<p>_Minimum and maximum volume / tube or bottle are those currently used. Others can be studied.</p> <p>_Minimum and maximum batch size can be adjusted, depending of the volume per bottle/tube.</p> <p>_Glass balls can be put in tubes and bottles. "Durham tube" can be put in tube.</p>										