

Introductory Offer



More Capacity

Centrifuge 5920 R

The new Centrifuge 5920 R delivers extraordinary high capacity in a very compact and ergonomic product design. Its state-of-the-art refrigeration system provides excellent cooling performance and keeps your temperature sensitive samples safe.

- > Max. capacity: 4 × 1000 mL or 52 × 50 mL conical
- > Outstanding rotor versatility
- > Space-saving design
- > Advanced temperature management
- > NEW: Eppendorf Conical Tubes, 15 mL and 50 mL



Ordering Information

| Cat. No. | Centrifuge 5920R without rotor | Introductory Promotional Price |
|---------------------------------|--|--------------------------------|
| 5948 000.085 | Centrifuge 5920 R, 230 V, 50-60 Hz, without rotor, with AU-plug | \$16,537.00 |
| Centrifuge 5920R Bundles | | |
| 5948 000.280 | Centrifuge 5920 R, 230 V, 50-60 Hz, incl. rotor S-4x1000, high-capacity buckets and adapter for 15 mL/50 mL conical tubes | \$20,468.00 |
| 5948 000.387 | Centrifuge 5920 R, 230 V, 50-60 Hz, incl. rotor S-4x1000, round buckets and adapter for 15 mL/50 mL conical tubes | \$20,136.00 |
| 5948 000.484 | Centrifuge 5920 R, 230 V, 50-60 Hz, incl. rotor S-4x1000, plate/tube buckets and adapter for 15 mL/50 mL conical tubes | \$20,136.00 |
| 5948 000.581 | Centrifuge 5920 R, 230 V, 50-60 Hz, incl. rotor S-4x750, round buckets and adapter for 15 mL/50 mL conical tubes | \$20,468.00 |
| 5948 000.689 | Centrifuge 5920 R, 230 V, 50-60 Hz, incl. rotor S-4x1000, plate/tube buckets and adapter for 13mm/16mm round-bottom tubes | \$20,136.00 |
| 5948 000.786 | Centrifuge 5920 R, 230 V, 50-60 Hz, incl. rotor S-4x1000, high-capacity buckets and adapter for 13mm/16mm round-bottom tubes | \$20,468.00 |

For more ordering information please visit www.eppendorf.com/centrifugation

General Terms and Conditions

- > Introductory promotional prices are available to you from the Authorised Distributor where you purchased the goods.
- > All prices are shown in NZD and do not include packing and freight charges.
- > Offers are only valid for customers residing in New Zealand.
- > Promotions cannot be combined with any other offers and are not applicable to demo equipment.
- > Contact Eppendorf or your Authorised Distributor for other rotor/bucket and adaptor combinations.

Authorised Distributors in New Zealand



Eppendorf South Pacific Pty. Ltd.
 Phone: +61 2 9889 5000
 E-mail: info@eppendorf.com.au
 Website: www.eppendorf.com

Visit: www.eppendorf.com/dealers

www.eppendorf.com/centrifugation