



The new generation enhances your lab

Develop your lab with **VIDAS® 3, the new generation of VIDAS®.**

Increased performance and confidence

Automation of the analysis procedure and complete traceability create confidence and peace of mind for each result you sign.

Greater autonomy

With access to a wide range of reagents in single-test format, you handle each request on demand, for same-day diagnostic results.

Fast adoption and a focus on productivity

Full connectivity to internal and external networks, a short training period and minimized operator maintenance time allow your lab technicians to focus on high level tasks.

Strength in numbers

VIDAS® gives you access to a worldwide community and a select range of services to consolidate your lab's performance and results.

The new generation is more connected



You benefit from **bioMérieux's experience** with VIDAS® services, ensuring your lab's continued development.

Instrument services

- Remote support services provide solutions, immediately
- Extensive network of local subsidiaries and distributors: skilled systems engineers and application specialists available wherever you are

Certification support and integrated quality control

- VIDAS® validation package for your instrument and methodology
- Support in obtaining ISO 15189 certification
- Compare your internal quality with your peers from around the world on the myQC platform

Online training

- Online training sessions complement VIDAS® on-site service
- Continually develop your lab's skills and expertise

A worldwide VIDAS® community

- The myvidas.com platform strengthens your expertise and makes the most of your VIDAS® solution
- Share information and experiences with your VIDAS® peers



Scan this code to discover more on the new generation and get access to **myvidas.com**

0713 / 0305585/089/GB/A / Document and/or pictures not legally binding. Modifications by bioMérieux can be made without prior notice / BIO MÉRIEUX, the blue logo and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies. Any other name or trademark is the property of its respective owner / BIO MÉRIEUX S.A. RCS Lyon 073 620 399 / Printed in France / MC CANV / ACS Numéro 682 068 941

bioMérieux S.A.
69280 Marcy l'Etoile
France
Tel. : + 33 (0)4 78 87 20 00
Fax : +33 (0)4 78 87 20 90
www.biomerieux.com
www.biomerieux-diagnostics.com



The new generation empowers your lab



Meet any request at any time

Your first priority is to serve clinicians and their patients. The VIDAS® brand is trusted around the world as a leader in on-demand immunoassay testings.

VIDAS® 3, the new generation empowers your lab.

VIDAS® 3 is designed with your specific needs and workflow in mind to provide test results that are precise and accountable.

The system is equipped with automation, traceability and continual access, for greater productivity.

Vidas 3

> VIDAS® 3 focuses on productivity

- **Facilitated lab workflow** with a simple answer to complex needs
- **Reduced operator time** for maintenance, to better focus on high-level tasks
- **Remote access:** a bioMérieux specialist provides both prompt troubleshooting and data for on-site instrument maintenance via screen sharing



> VIDAS® 3 allows fast adoption

- **Intuitive automation:** test loading is streamlined, ensuring standardization within your lab
- **Fully automated pipetor action** with single-use, disposable tips, to avoid contamination
- Specific **management of emergency tests**
- Any **qualified lab technician** can run routine operations




> VIDAS® 3 is designed for continual access

- Three separate racks each containing three configurable segments can hold a total of **27 samples, diluents and calibrators** in the temperature-controlled loading bay
- **Four independent compartments** can each run 3 tests
- VIDAS® 3 can process up to **36 tests per hour**



> VIDAS® 3 handles test security

- **Less human manipulation** from the sample to the result
- **Sample and reagent barcodes** are scanned upon test loading for complete traceability of all test components
- **Automated calibration, on-board dilution**
- **Integrated quality control management**
- **Bi-directional LIS connectivity**



> VIDAS® 3 uses VIDAS® reagents

- **Over 100 diagnostic tests** in single-test format for cardiovascular and infectious diseases, cancer, fertility, thyroid and allergies
- **On-demand testing:** 1 patient, 1 test, 1 result
- **Ready-to-use reagents**
- **One order reference per test:** kits contain all components needed for one test



Height: 61 cm; 24 in / Width: 75 cm; 29.5 in / Depth: 64 cm; 25.5 in