

Centri-HUGE Deals with Eppendorf this Winter! See Inside!

Eppendorf Centrifuges

eppendorf



15%
OFF LIST

See page 2 for more details.

Eppendorf Mobile Centrifuge Tables

eppendorf



10%
OFF LIST

See page 3 for more details.

Eppendorf ThermoMixer & SmartBlocks



15%
OFF LIST

20%
OFF LIST

See page 3 for more details.

Hanna Meter & Calibration Solutions

Buy a new meter,
Get 50% Off any
calibration solutions!

Get 20% Off when
you trade in your old
testing equipment!



HANNA
instruments

See page 4 for more details.

Greiner Serological Pipettes

greiner bio-one

20%
OFF LIST



See page 4 for more details.

Centri-HUGE Deals with Epp

Centrifuges

Centrifuge 5425 is your new flagship microcentrifuge!

It boasts a wide range of rotors, providing excellent versatility, while running at low noise levels. With its 24-place capacity and speed up to 21,300 x g, it is perfectly equipped for all modern molecular biology applications.

Description	Code	Pack Size	List Price (excl. GST)	Promo Price (excl. GST)
Centrifuge 5425, Foil with Rotor FA-24 x 2	EP5405000085	Ea	\$4,054.66	\$3,446.46
Centrifuge 5425, Knob with Rotor FA-24 x 2	EP5405000484	Ea	\$4,054.66	\$3,446.46
Best of Eppendorf Bundle 1: Centrifuge 5420 with Research Plus Pipette 3pk (Option 2), 2 Boxes of SafeLock 1.5ml Tubes Ep Quality and Eppendorf Tube Racks for 1.5/2.0 mL Tubes	EP2232000125	Ea	\$5,551.81	\$4,720.00



Centrifuge 5420 - designed with your molecular biology routine applications in mind!

Centrifuge 5420 supports your bread-and-butter spin column-based nucleic acid purification and gel-extraction workflows. It is the 5425's little sibling and is a perfect solution for labs that don't require the extra versatility of the 5425.

Description	Code	Pack Size	List Price (excl. GST)	Promo Price (excl. GST)
Centrifuge 5420 with Rotor FA-24 x 2	EP5420000083	Ea	\$3,339.48	\$2,838.56
Best of Eppendorf Bundle 2: Centrifuge 5420 with Research Plus Pipette 3pk (Option 2), 2 Boxes of SafeLock 1.5ml Tubes Ep Quality and Eppendorf Tube Racks for 1.5/2.0 mL Tubes	EP2232000124	Ea	\$4,849.15	\$4,122.00

Centrifuge 5430R - user-friendly, durable and versatile!

Unique cross-over Centrifuge 5430 R combines the best features of a refrigerated microcentrifuge and multipurpose centrifuge in one instrument. This is a perfect model for labs looking at something more comprehensive than your average microcentrifuge, but less comprehensive than a large work-horse.

Description	Code	Pack Size	List Price (excl. GST)	Promo Price (excl. GST)
Centrifuge 5430R, Foil with Rotor FA-45-30-11	EP542800008	Ea	\$12,289.65	\$10,446.20
Centrifuge 5430R, Knob with Rotor FA-45-30-11	EP542800048	Ea	\$12,289.65	\$10,446.20



*While Stocks Last

Eppendorf this Winter!*

eppendorf

NEW Space Saver!

Mobile Table with every 58XX & 59XX Centrifuge!

The new Eppendorf Mobile Tables save precious benchtop space and provide the freedom and flexibility to use your Centrifuge 5804/5804 R, 5810/5810 R, 5910 R or 5920 R wherever you need it.

Description	Code	Pack Size	List Price (excl. GST)	Promo Price (excl. GST)
Low Mobile Table for Centrifuge	EP5921000010	Ea	\$2,990.93	\$2,691.84
High Mobile Table for Centrifuge	EP5921000012	Ea	\$3,738.67	\$3,364.81



ThermoMixer™ C

ThermoMixer™ C Bundles

Improve your assay results by dry incubation and mixing of your samples simultaneously - the Eppendorf ThermoMixer™ C is for heating, cooling, and mixing almost any of your lab vessels.

Description	Code	Pack Size	List Price (excl. GST)	Promo Price (excl. GST)
ThermoMixer Bundle with 1.5mL Block- NZ	EP2232000089	Ea	\$6,400.75	\$5,440.64
ThermoMixer Plate Bundle with Thermotop- NZ	EP2232000090	Ea	\$7,818.95	\$6,646.10



ThermoMixer™ C SmartBlocks

These individually sensor-controlled SmartBlocks are the best solution for maintaining accurate, uniform temperatures. A wide variety of these flexible thermoblocks for all common lab vessels support a broad range of applications - for use with Eppendorf ThermoMixer™ C (heating/ cooling and mixing) and ThermoStat C (heating/cooling).

Description	Code	Pack Size	List Price (excl. GST)	Promo Price (excl. GST)
SmartBlock PCR 96, includes Lid	EP5306000006	Ea	\$1,619.76	\$1,295.81
SmartBlock PCR 384, includes Lid	EP5307000000	Ea	\$1,619.76	\$1,295.81
ThermoTop with Condenser Protect-Technology	EP5308000003	Ea	\$1,101.44	\$ 881.15
SmartBlock 8 x 5mL Tubes	EP5309000007	Ea	\$1,303.00	\$1,042.40
SmartBlock DWP 1000 µL	EP5310000002	Ea	\$1,619.76	\$1,295.81
SmartBlock DWP 500 µL	EP5316000004	Ea	\$1,619.76	\$1,295.81
SmartExtender 1.5 mL	EP5322000008	Ea	\$1,713.34	\$1,370.67
SmartBlock 24 x 1.5mL Tubes	EP5360000038	Ea	\$1,303.00	\$1,042.40
SmartBlock 24 x 0.5mL Tubes	EP5361000031	Ea	\$1,303.00	\$1,042.40
SmartBlock 24 x 2.0mL Tubes	EP5362000035	Ea	\$1,303.00	\$1,042.40
SmartBlock Plates, includes Lid	EP5363000039	Ea	\$1,619.76	\$1,295.81
SmartBlock 24 x 12mm Tubes	EP5364000024	Ea	\$1,619.76	\$1,295.81
SmartBlock 4 x 50mL Conical Tubes	EP5365000028	Ea	\$1,303.00	\$1,042.40
SmartBlock 8 x 15mL Conical Tubes	EP5366000021	Ea	\$1,303.00	\$1,042.40
SmartBlock Cryo 24 x 1, 5mL-2mL	EP5367000025	Ea	\$1,303.00	\$1,042.40



*While Stocks Last

Buy a Hanna Meter and Get 50% Off any Associated Calibration Solutions!

Buy any meter and receive 50% off any calibration solutions.

To take advantage of this promotion, please email info@mediray.co.nz with your desired meter and calibration solution.



Terms and Conditions


- The 50% off discount only applies for up to three units of a calibration solution.
- If you would like a range of different standards, the 50% discount will be applied to each (for up to 3 units in total).
- The discount cannot be applied retrospectively.

Make the most of your old testing equipment by trading it in!

If you have a meter or titrator for measuring standard analytical parameters that's old or isn't working as it should be, then why not consider trading it in?

In return, you'll get 20% off Hanna instruments testers, meters and titrators.*

Contact info@mediray.co.nz with details to determine if your meter qualifies.



Terms and Conditions

- *The validity of the discounted item will be judged on a case-by-case basis due to the scope and complexity of products on offer, in other words, like must be traded for like.
- The trade-in discount will only apply to a single item, and not the entirety of the order if more than one item is to be purchased.
- Trade-in discounts are not cumulative. E.g. two separate instruments traded in simultaneously do not equate to a 40% discount on a single item.
- This discount does not apply to consumables or accessories (including electrodes/probes).

Greiner CELLSTAR Serological Pipettes

Greiner pipettes are precision extruded from the highest quality, optical-grade polystyrene. All volumes are identified by the standard international color code convention. Accuracy is maintained by high surface polish, both internally and externally, to reduce meniscus formation and ensure complete liquid expulsion with minimum carry-over.

Precision is aided by fine-line graduation printing providing accuracy of better than $\pm 2\%$.

- 1 mL, 2 mL and 5 mL sizes incorporate drawn tips while larger sizes feature molded tips to provide maximum accuracy
- Packed individually, then double packed for extra protection
- Designed to fit all standard pipetting devices
- Cotton plugged for safety
- Sterilised by γ -irradiation
- DNase, RNase and human DNA free

20%
OFF LIST


greiner bio-one



Description	Code	Pack Size	List Price (excl. GST)	Promo Price (excl. GST)
CELLSTAR Serological Pipette 5ml	GR606180	200	\$84.53	\$67.60
CELLSTAR Serological Pipette 10ml	GR607180	200	\$84.53	\$67.60
CELLSTAR Serological Pipette 25ml	GR760180	200	\$162.75	\$130.20
CELLSTAR Serological Pipette 50ml	GR768180	100	\$154.88	\$123.90